

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA AKZENT Plus PASTE FLUID

Revision date: 02.06.2022

Product code: 277

Page 1 of 8

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

##### 1.1. Product identifier

VITA AKZENT Plus PASTE FLUID

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

###### Use of the substance/mixture

Use as laboratory reagent

##### 1.3. Details of the supplier of the safety data sheet

|                         |                                       |                             |
|-------------------------|---------------------------------------|-----------------------------|
| Company name:           | VITA Zahnfabrik H.Rauter GmbH & Co.KG |                             |
| Street:                 | Spitalgasse 3                         |                             |
| Place:                  | D-79713 Bad Säckingen                 |                             |
| Post-office box:        | 1338                                  |                             |
|                         | D-79704 Bad Säckingen                 |                             |
| Telephone:              | +49(0)7761-562-0                      | Telefax: +49(0)7761-562-299 |
| e-mail:                 | info@vita-zahnfabrik.com              |                             |
| Contact person:         | regulatory affairs                    |                             |
| e-mail:                 | info@vita-zahnfabrik.com              |                             |
| Internet:               | www.vita-zahnfabrik.com               |                             |
| Responsible Department: | Regulatory Affairs                    |                             |

##### 1.4. Emergency telephone number:

+49-(0)761-19240

##### Further Information

medical device

#### SECTION 2: Hazards identification

##### 2.1. Classification of the substance or mixture

###### Regulation (EC) No 1272/2008

This mixture is not classified as hazardous in accordance with Regulation (EC) No 1272/2008.

##### 2.2. Label elements

###### Regulation (EC) No 1272/2008

###### Special labelling of certain mixtures

|        |   |
|--------|---|
| EUH208 | Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. |
|--------|---|

|        |   |
|--------|---|
| EUH210 | Safety data sheet available on request. |
|--------|---|

##### 2.3. Other hazards

No information available.

#### SECTION 3: Composition/information on ingredients

##### 3.2. Mixtures

###### Chemical characterization

Substance, organic

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA AKZENT Plus PASTE FLUID

Revision date: 02.06.2022

Product code: 277

Page 2 of 8

#### Hazardous components

| CAS No     | Chemical name  |              |          | Quantity |
|------------|--|--------------|----------|----------|
|            | EC No  | Index No     | REACH No |          |
|            | Classification (Regulation (EC) No 1272/2008)  |              |          |          |
| 55965-84-9 | reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1)  |              |          | < 0.1 %  |
|            | -  | 613-167-00-5 |          |          |
|            | Acute Tox. 2, Acute Tox. 2, Acute Tox. 3, Skin Corr. 1C, Eye Dam. 1, Skin Sens. 1A, Aquatic Acute 1, Aquatic Chronic 1; H330 H310 H301 H314 H318 H317 H400 H410 EUH071 |              |          |          |

Full text of H and EUH statements: see section 16.

#### Specific Conc. Limits, M-factors and ATE

| CAS No     | EC No   | Chemical name   | Quantity |
|------------|---|---|----------|
|            | Specific Conc. Limits, M-factors and ATE  |   |          |
| 55965-84-9 | -   | reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) | < 0.1 %  |
|            | inhalation: ATE = 0.5 mg/l (vapours); inhalation: ATE = 0.05 mg/l (dusts or mists); dermal: ATE = 50 mg/kg; oral: ATE = 100 mg/kg Skin Corr. 1C; H314: >= 0.6 - 100 Skin Irrit. 2; H315: >= 0.06 - < 0.6 Eye Dam. 1; H318: >= 0.6 - 100 Eye Irrit. 2; H319: >= 0.06 - < 0.6 Skin Sens. 1A; H317: >= 0.0015 - 100<br>M acute; H400: M=100<br>M chron.; H410: M=100 |   |          |

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

##### After inhalation

Provide fresh air. If breathing is irregular or stopped, administer artificial respiration. Medical treatment necessary.

##### After contact with skin

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention.

##### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

##### After ingestion

Rinse mouth immediately and drink plenty of water.

#### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

##### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

#### 5.2. Special hazards arising from the substance or mixture

Non-flammable.

#### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

#### Additional information

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA AKZENT Plus PASTE FLUID

Revision date: 02.06.2022

Product code: 277

Page 3 of 8

#### SECTION 6: Accidental release measures

##### 6.1. Personal precautions, protective equipment and emergency procedures

###### General advice

Use personal protection equipment.

##### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

##### 6.3. Methods and material for containment and cleaning up

###### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

##### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

#### SECTION 7: Handling and storage

##### 7.1. Precautions for safe handling

###### Advice on safe handling

No special measures are necessary.

###### Advice on protection against fire and explosion

No special fire protection measures are necessary.

###### Advice on general occupational hygiene

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat, drink, smoke, sniff.

##### 7.2. Conditions for safe storage, including any incompatibilities

###### Requirements for storage rooms and vessels

Keep container tightly closed.

###### Hints on joint storage

No special measures are necessary.

##### 7.3. Specific end use(s)

Use as laboratory reagent

#### SECTION 8: Exposure controls/personal protection

##### 8.1. Control parameters

##### 8.2. Exposure controls



###### Individual protection measures, such as personal protective equipment

###### Eye/face protection

Wear eye/face protection.

###### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA AKZENT Plus PASTE FLUID

Revision date: 02.06.2022

Product code: 277

Page 4 of 8

specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Suitable material: NBR (Nitrile rubber)

#### Skin protection

Use of protective clothing.

#### Respiratory protection

Technical ventilation of workplace Open windows to ensure natural ventilation.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

|                 |                |
|-----------------|----------------|
| Physical state: | Liquid         |
| Colour:         | colourless     |
| Odour:          | characteristic |

#### Changes in the physical state

|   |                |
|---|----------------|
| Melting point/freezing point:                             | not determined |
| Boiling point or initial boiling point and boiling range: | 207 °C         |
| Flash point:  | 109 °C         |

#### Flammability

|               |                |
|---------------|----------------|
| Solid/liquid: | not applicable |
| Gas:          | not applicable |

#### Explosive properties

The product is not: Explosive.

|                         |                |
|-------------------------|----------------|
| Lower explosion limits: | not determined |
| Upper explosion limits: | not determined |

#### Self-ignition temperature

|                            |                |
|----------------------------|----------------|
| Solid:                     | not applicable |
| Gas:                       | not applicable |
| Decomposition temperature: | not determined |
| pH-Value:                  | not determined |

#### Solubility in other solvents

not determined

|  |                           |
|--|---------------------------|
| Partition coefficient n-octanol/water: | not determined            |
| Vapour pressure:                       | not determined            |
| Density:                               | 1,06200 g/cm <sup>3</sup> |
| Relative vapour density:               | not determined            |

### 9.2. Other information

#### Information with regard to physical hazard classes

Oxidizing properties  
Not oxidising.

#### Other safety characteristics

|                   |                |
|-------------------|----------------|
| Solid content:    | 0,0 %          |
| Evaporation rate: | not determined |

#### Further Information

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA AKZENT Plus PASTE FLUID

Revision date: 02.06.2022

Product code: 277

Page 5 of 8

#### SECTION 10: Stability and reactivity

##### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

##### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

##### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

##### 10.4. Conditions to avoid

none

##### 10.5. Incompatible materials

No information available.

##### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

#### SECTION 11: Toxicological information

##### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

###### Acute toxicity

| CAS No     | Chemical name   |               |         |        |        |
|------------|---|---------------|---------|--------|--------|
|            | Exposure route  | Dose          | Species | Source | Method |
| 55965-84-9 | reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) |               |         |        |        |
|            | oral  | ATE 100 mg/kg |         |        |        |
|            | dermal  | ATE 50 mg/kg  |         |        |        |
|            | inhalation vapour   | ATE 0.5 mg/l  |         |        |        |
|            | inhalation dust/mist  | ATE 0.05 mg/l |         |        |        |

###### Additional information on tests

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

#### SECTION 12: Ecological information

##### 12.1. Toxicity

The product is not: Ecotoxic.

##### 12.2. Persistence and degradability

The product has not been tested.

##### 12.3. Bioaccumulative potential

The product has not been tested.

##### 12.4. Mobility in soil

The product has not been tested.

##### 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

The product has not been tested.

##### 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA AKZENT Plus PASTE FLUID

Revision date: 02.06.2022

Product code: 277

Page 6 of 8

#### 12.7. Other adverse effects

No information available.

#### **Further information**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### **Disposal recommendations**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### **Contaminated packaging**

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

### SECTION 14: Transport information

#### **Land transport (ADR/RID)**

|  |  |
|--|--|
| <u>14.1. UN number or ID number:</u>     | No dangerous good in sense of this transport regulation. |
| <u>14.2. UN proper shipping name:</u>    | No dangerous good in sense of this transport regulation. |
| <u>14.3. Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| <u>14.4. Packing group:</u>              | No dangerous good in sense of this transport regulation. |

#### **Inland waterways transport (ADN)**

|  |  |
|--|--|
| <u>14.1. UN number or ID number:</u>     | No dangerous good in sense of this transport regulation. |
| <u>14.2. UN proper shipping name:</u>    | No dangerous good in sense of this transport regulation. |
| <u>14.3. Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| <u>14.4. Packing group:</u>              | No dangerous good in sense of this transport regulation. |

#### **Marine transport (IMDG)**

|  |  |
|--|--|
| <u>14.1. UN number or ID number:</u>     | No dangerous good in sense of this transport regulation. |
| <u>14.2. UN proper shipping name:</u>    | No dangerous good in sense of this transport regulation. |
| <u>14.3. Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| <u>14.4. Packing group:</u>              | No dangerous good in sense of this transport regulation. |

#### **Air transport (ICAO-TI/IATA-DGR)**

|  |  |
|--|--|
| <u>14.1. UN number or ID number:</u>     | No dangerous good in sense of this transport regulation. |
| <u>14.2. UN proper shipping name:</u>    | No dangerous good in sense of this transport regulation. |
| <u>14.3. Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| <u>14.4. Packing group:</u>              | No dangerous good in sense of this transport regulation. |

#### 14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

#### 14.6. Special precautions for user

No information available.

#### 14.7. Maritime transport in bulk according to IMO instruments

not applicable

### SECTION 15: Regulatory information

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

##### **EU regulatory information**

Restrictions on use (REACH, annex XVII):

Entry 75

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA AKZENT Plus PASTE FLUID

Revision date: 02.06.2022

Product code: 277

Page 7 of 8

2004/42/EC (VOC): 0,014 % (0,143 g/l)  
Information according to 2012/18/EU (SEVESO III): Not subject to 2012/18/EU (SEVESO III):

#### National regulatory information

Water hazard class (D): 2 - obviously hazardous to water  
Skin resorption/Sensitization: Causes allergic hypersensitivity reactions.

#### 15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

### SECTION 16: Other information

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%  
CLP: Classification, labelling and Packaging  
REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
EmS: Emergency Schedules  
MFAG: Medical First Aid Guide  
ICAO: International Civil Aviation Organization  
MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
IBC: Intermediate Bulk Container  
VOC: Volatile Organic Compounds  
SVHC: Substance of Very High Concern  
For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA AKZENT Plus PASTE FLUID

Revision date: 02.06.2022

Product code: 277

Page 8 of 8

#### Relevant H and EUH statements (number and full text)

|        |   |
|--------|---|
| H301   | Toxic if swallowed.   |
| H310   | Fatal in contact with skin.   |
| H314   | Causes severe skin burns and eye damage.  |
| H317   | May cause an allergic skin reaction.  |
| H318   | Causes serious eye damage.  |
| H330   | Fatal if inhaled.   |
| H400   | Very toxic to aquatic life.   |
| H410   | Very toxic to aquatic life with long lasting effects.   |
| EUH071 | Corrosive to the respiratory tract.   |
| EUH208 | Contains reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction. |
| EUH210 | Safety data sheet available on request.   |

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*