

## Safety Data Sheet

according to WHMIS

### VITA AMBRIA INVEST F

Revision date: 30.03.2023

Product code: 329

Page 1 of 6

## 1. Identification

### Product identifier

VITA AMBRIA INVEST F

### Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Auxiliary for manufacture of dental products  
Laboratory chemicals

### Details of the supplier of the safety data sheet

|                         |                                       |                             |
|-------------------------|---------------------------------------|-----------------------------|
| Company name:           | VITA Zahnfabrik H.Rauter GmbH & Co.KG |                             |
| Street:                 | Spitalgasse 3                         |                             |
| Place:                  | D-79713 Bad Säckingen                 |                             |
| Post-office box:        | 1338                                  |                             |
|                         | D-79704 Bad Säckingen                 |                             |
| Telephone:              | +49(0)7761-562-0                      | Telefax: +49(0)7761-562-299 |
| e-mail:                 | info@vita-zahnfabrik.com              |                             |
| Contact person:         | regulatory affairs                    |                             |
| e-mail:                 | info@vita-zahnfabrik.com              |                             |
| Internet:               | www.vita-zahnfabrik.com               |                             |
| Responsible Department: | Regulatory Affairs                    |                             |

**Emergency telephone number:** +49-(0)761-19240

## 2. Hazard identification

### Classification of the substance or mixture

#### WHMIS 2015

This mixture is not classified as hazardous in accordance with WHMIS 2015.

### Label elements

#### Other hazards

No information available.

## 3. Composition/information on ingredients

### Mixtures

#### Chemical characterization

Mixtures Product/Substance is inorganic.

#### Hazardous components

none (according to WHMIS)

## 4. First-aid measures

### Description of first aid measures

#### After inhalation

Provide fresh air.

#### After contact with skin

Wash with plenty of water. Take off contaminated clothing and wash it before reuse.

#### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

## Safety Data Sheet

according to WHMIS

### VITA AMBRIA INVEST F

Revision date: 30.03.2023

Product code: 329

Page 2 of 6

#### After ingestion

Rinse mouth immediately and drink plenty of water.

#### Most important symptoms and effects, whether acute or delayed

No information available.

#### Indication of immediate medical attention and special treatment needed

Treat symptomatically.

## 5. Fire-fighting measures

#### Extinguishing media

##### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

#### Specific hazards arising from the hazardous product

Non-flammable.

#### Special protective equipment and precautions for fire-fighters

In case of fire: Wear self-contained breathing apparatus.

## 6. Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

##### General advice

Use personal protection equipment.

#### Environmental precautions

No special environmental measures are necessary. Clean contaminated articles and floor according to the environmental legislation.

#### Methods and material for containment and cleaning up

##### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## 7. Handling and storage

#### Precautions for safe handling

##### Advice on safe handling

No special measures are necessary.

##### Advice on protection against fire and explosion

No special fire protection measures are necessary.

##### Advice on general occupational hygiene

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat, drink, smoke, sniff.

#### Conditions for safe storage, including any incompatibilities

##### Requirements for storage rooms and vessels

Keep container tightly closed.

##### Hints on joint storage

No special measures are necessary.

## Safety Data Sheet

according to WHMIS

### VITA AMBRIA INVEST F

Revision date: 30.03.2023

Product code: 329

Page 3 of 6

#### Further information on storage conditions

storage temperature: 5 - 35° C

### 8. Exposure controls/Personal protection

#### Control parameters

#### Exposure controls

#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

##### Skin protection

Use of protective clothing.

##### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

### 9. Physical and chemical properties

#### Information on basic physical and chemical properties

Physical state: Liquid  
 Colour: white

#### Changes in the physical state

Melting point/freezing point: 0 °C  
 Boiling point or initial boiling point and boiling range: 100 °C  
 Flash point: not applicable

#### Flammability

Solid/liquid: not applicable  
 Gas: not applicable

#### Explosive properties

The product is not: Explosive.

Lower explosive limits: not determined  
 Upper explosive limits: not determined

#### Self-ignition temperature

Solid: not applicable  
 Gas: not applicable

Decomposition temperature: not determined

pH-Value: 9,5

Viscosity / dynamic: < 20 mPa·s  
 (at 20 °C)

Water solubility: easily soluble

#### Solubility in other solvents

not determined

## Safety Data Sheet

according to WHMIS

### VITA AMBRIA INVEST F

Revision date: 30.03.2023

Product code: 329

Page 4 of 6

|  |                           |
|--|---------------------------|
| Partition coefficient n-octanol/water: | not determined            |
| Vapour pressure:<br>(at 20 °C)         | 23 hPa                    |
| Density:                               | 1,1-1,4 g/cm <sup>3</sup> |
| Relative vapour density:               | not determined            |

#### Other information

##### Information with regard to physical hazard classes

Oxidizing properties  
Not oxidising.

##### Other safety characteristics

|                   |                |
|-------------------|----------------|
| Solid content:    | not determined |
| Evaporation rate: | not determined |

##### Further Information

## 10. Stability and reactivity

#### Reactivity

No hazardous reaction when handled and stored according to provisions.

#### Chemical stability

The product is stable under storage at normal ambient temperatures.

#### Possibility of hazardous reactions

No known hazardous reactions.

#### Conditions to avoid

Steel, Aluminium, Copper

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

#### Information on toxicological effects

##### Additional information on tests

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

The product is not: Ecotoxic.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

#### Mobility in soil

The product has not been tested.

#### Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

## Safety Data Sheet

according to WHMIS

### VITA AMBRIA INVEST F

Revision date: 30.03.2023

Product code: 329

Page 5 of 6

#### Other adverse effects

No information available.

#### **Further information**

Avoid release to the environment.

### 13. Disposal considerations

#### Waste treatment methods

##### **Disposal recommendations**

Dispose of waste according to applicable legislation.

##### **Contaminated packaging**

Wash with plenty of water. Completely emptied packages can be recycled.

### 14. Transport information

#### Hazard classes:

##### **Marine transport (IMDG)**

##### UN number or ID number:

No dangerous good in sense of this transport regulation.

##### United Nations proper shipping name:

No dangerous good in sense of this transport regulation.

##### Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### Packing group:

No dangerous good in sense of this transport regulation.

##### **Air transport (ICAO-TI/IATA-DGR)**

##### UN number or ID number:

No dangerous good in sense of this transport regulation.

##### United Nations proper shipping name:

No dangerous good in sense of this transport regulation.

##### name:

##### Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### Packing group:

No dangerous good in sense of this transport regulation.

#### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

### 15. Regulatory information

#### Canadian regulations

### 16. Other information

#### **Abbreviations and acronyms**

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

CAS: Chemical Abstracts Service

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

LL50: Lethal loading, 50%

EL50: Effect loading, 50%

EC50: Effective Concentration 50%

## Safety Data Sheet

according to WHMIS

### VITA AMBRIA INVEST F

Revision date: 30.03.2023

Product code: 329

Page 6 of 6

ErC50: Effective Concentration 50%, growth rate

NOEC: No Observed Effect Concentration

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

RID: Regulations concerning the international carriage of dangerous goods by rail

ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

IMDG: International Maritime Code for Dangerous Goods

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

IATA: International Air Transport Association

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*