


VITA ENAMIC® Polishing Set clinical/technical

For simple pre- and high-gloss polishing of VITA ENAMIC hybrid ceramic



Date of issue: 2022-09

VITA – perfect match. **VITA**

Product description

The VITA ENAMIC polishing sets were developed for material-specific surface treatment of restorations made of hybrid ceramics in dental practices and laboratories. The sets include various polishing instru-

ments for pre- and high-gloss polishing. With these instruments, occlusal surfaces, cusps, fissures and contact points of the restoration can be polished to a high gloss in accordance with the material.

Product benefits

Effective final results

These instruments enable the user to achieve a high gloss finish of ENAMIC restorations in a very short time. Accurate concentricity, matched grit sizes and the individual geometries of the instruments ensure good results.

Simple and safe handling

The instruments provide superior material removal performance. Efficient handling without the mandatory use of polishing paste allows for easy and fast polishing. The ability to sterilize the clinical instruments ensures safety.

Material-specific treatment

These instruments, which have been developed especially for VITA ENAMIC, can be used for material-specific reworking.

Product data

VITA ENAMIC Polishing Set clinical: Includes a total of six instruments for the contra-angle. It includes three instruments each for pre- and high-gloss polishing.

VITA ENAMIC Polishing Set technical: Includes a total of six instruments for the handpiece. It includes three instruments each for pre- and high-gloss polishing.

VITA ENAMIC® Polishing Sets

| Designation | Details | Content/pieces | REF |
|-------------------------------------|--|----------------|------------|
| VITA ENAMIC Polishing Set clinical | Two-stage polishing system with pre- and high-gloss polishers for the contra-angle | 6 | EENPSETCV1 |
| VITA ENAMIC Polishing Set technical | Two-stage polishing system with pre- and high-gloss polishers for the handpiece | 6 | EENPSETTV1 |

VITA ENAMIC® Polishing Set clinical - refill packs

| Designation | Details | Content/pieces | REF |
|------------------------------------|-----------------------------------|----------------|-----------|
| Pre-polishing (pink) | SiC instruments (silicon carbide) | | |
| Brush, VI-EB14m | | 3 | ERWEB14M3 |
| lens, VI-EL10m | | 6 | ERWEL10M6 |
| cup, VI-EC7m | | 6 | ERWEC7M6 |
| High-gloss polishing (grey) | Diamond instruments | | |
| Brush, VI-EB14f | | 3 | ERWEB14F3 |
| lens, VI-EL10f | | 6 | ERWEL10F6 |
| cup, VI-EC6f | | 6 | ERWEC6F6 |

VITA ENAMIC® Polishing Set technical - refill packs

| Designation | Details | Content/pieces | REF |
|------------------------------------|-----------------------------------|----------------|-----------|
| Pre-polishing (pink) | SiC instruments (silicon carbide) | | |
| Brush, VI-EB17m | | 3 | ERHEB17M3 |
| tip, VI-ES5m | | 6 | ERHES5M6 |
| wheel, VI-ER15m | | 6 | ERHER15M6 |
| High-gloss polishing (grey) | Diamond instruments | | |
| Brush, VI-EB17f | | 3 | ERHEB17F3 |
| tip, VI-ES5f | | 6 | ERHES5F6 |
| wheel, VI-ER15f | | 6 | ERHER15F6 |

■ EVE Ernst Vetter GmbH · Neureutstr. 6 · 75210 Keltern · Germany

EVE Ernst Vetter GmbH has been certified and the following products bear the mark **CE 0483**:

VITA ENAMIC Polishing Set clinical and the corresponding refill packs.



VITA

Vertrieb: VITA Zahnfabrik H. Rauter GmbH & Co. KG · Spitalgasse 3 · 79713 Bad Säckingen · Germany · Tel. +49/7761/562-0 · Fax +49/7761/562-299 · www.vita-zahnfabrik.com