

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA FIRING PASTE

Revision date: 21.01.2020

Product code: 276

Page 1 of 8

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

##### 1.1. Product identifier

VITA FIRING PASTE

CAS No: -  
 Index No: 650-017-00-8  
 EC No: -  
 UFI: VV00-70P2-1006-FDGP

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

###### Use of the substance/mixture

Auxiliary for manufacture of dental products

##### 1.3. Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H.Rauter GmbH & Co.KG  
 Street: Spitalgasse 3  
 Place: D-79713 Bad Säckingen  
 Post-office box: 1338  
 D-79704 Bad Säckingen  
 Telephone: +49(0)7761-562-0  
 e-mail: info@vita-zahnfabrik.com  
 Contact person: regulatory affairs  
 e-mail: info@vita-zahnfabrik.com  
 Internet: www.vita-zahnfabrik.com  
 Responsible Department: Regulatory Affairs

Telefax: +49(0)7761-562-299

##### 1.4. Emergency telephone number:

+49-(0)761-19240

#### SECTION 2: Hazards identification

##### 2.1. Classification of the substance or mixture

###### Regulation (EC) No 1272/2008

Skin Irrit. 2; H315  
 Carc. 1B; H350i

Full text of hazard statements: see SECTION 16.

##### 2.2. Label elements


###### Regulation (EC) No 1272/2008

###### Hazard components for labelling

Aluminosilicate fibre, dioxosilane

Signal word: Danger

###### Pictograms:



###### Hazard statements

H315 Causes skin irritation.  
 H350i May cause cancer by inhalation.

###### Precautionary statements

P202 Do not handle until all safety precautions have been read and understood.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA FIRING PASTE

Revision date: 21.01.2020

Product code: 276

Page 2 of 8

P280 Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.  
 P362+P364 Take off contaminated clothing and wash it before reuse.  
 P405 Store locked up.

#### 2.3. Other hazards

No information available.

### SECTION 3: Composition/information on ingredients

#### 3.2. Mixtures

##### Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (Regulation (EC) No 1272/2008)			
142844-00-6	Aluminosilicate fibre, dioxosilane			20 - < 25 %
	604-314-4		01-2119458050-50	
	Carc. 1B, Skin Irrit. 2; H350i H315			

Full text of H and EUH statements: see section 16.

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

##### General information

When in doubt or if symptoms are observed, get medical advice.

##### After inhalation

Provide fresh air.

##### After contact with skin

Wash with plenty of water. Immediately remove any contaminated clothing, shoes or stockings.

##### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

##### After ingestion

Rinse mouth immediately and drink plenty of water.

#### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

##### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

#### 5.2. Special hazards arising from the substance or mixture

Non-flammable.

#### 5.3. Advice for firefighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

#### Additional information

Suppress gases/vapours/mists with water spray jet.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA FIRING PASTE

Revision date: 21.01.2020

Product code: 276

Page 3 of 8

#### SECTION 6: Accidental release measures

##### 6.1. Personal precautions, protective equipment and emergency procedures

###### General advice

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

##### 6.2. Environmental precautions

No special environmental measures are necessary.

##### 6.3. Methods and material for containment and cleaning up

###### Other information

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

##### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

#### SECTION 7: Handling and storage

##### 7.1. Precautions for safe handling

###### Advice on safe handling

If handled uncovered, arrangements with local exhaust ventilation have to be used. Avoid dust formation. Do not breathe dust.

###### Advice on protection against fire and explosion

No special fire protection measures are necessary.

###### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

##### 7.2. Conditions for safe storage, including any incompatibilities

###### Requirements for storage rooms and vessels

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Provide adequate ventilation as well as local exhaust at critical locations.

###### Hints on joint storage

No special measures are necessary.

##### 7.3. Specific end use(s)

Auxiliary for manufacture of dental products

#### SECTION 8: Exposure controls/personal protection

##### 8.1. Control parameters

##### 8.2. Exposure controls



###### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA FIRING PASTE

Revision date: 21.01.2020

Product code: 276

Page 4 of 8

#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatril P NBR (Nitrile rubber) Breakthrough time: 480 min

##### Skin protection

Use of protective clothing.

##### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Avoid dust formation. Dust must be exhausted directly at the point of origin.

### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state:

Colour: white

Odour: none

##### Changes in the physical state

Melting point/freezing point: 1752 °C

Boiling point or initial boiling point and boiling range: ?

Flash point: ?

##### Flammability

Solid/liquid: not determined

Gas: not applicable

##### Explosive properties

The product is not: Explosive.

Lower explosion limits: not determined

Upper explosion limits: not determined

##### Self-ignition temperature

Solid: not determined

Gas: not applicable

Decomposition temperature: not determined

pH-Value: not determined

Water solubility: No

##### Solubility in other solvents

not determined

Partition coefficient n-octanol/water: not determined

Vapour pressure: <=1100 hPa  
(at 50 °C)

Density: 2,5-2,72 g/cm³

Relative vapour density: not determined

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA FIRING PASTE

Revision date: 21.01.2020

Product code: 276

Page 5 of 8

#### 9.2. Other information

##### Information with regard to physical hazard classes

Oxidizing properties

Not oxidising.

##### Other safety characteristics

Solid content:

50 %

Evaporation rate:

not determined

##### Further Information

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

#### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

#### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

#### 10.4. Conditions to avoid

none

#### 10.5. Incompatible materials

No information available.

#### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

### SECTION 11: Toxicological information

#### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

##### Acute toxicity

Based on available data, the classification criteria are not met.

##### Irritation and corrosivity

Causes skin irritation.

Serious eye damage/eye irritation: Based on available data, the classification criteria are not met.

##### Sensitising effects

Based on available data, the classification criteria are not met.

##### Carcinogenic/mutagenic/toxic effects for reproduction

May cause cancer by inhalation. (Aluminosilicate fibre, dioxosilane)

Germ cell mutagenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

##### STOT-single exposure

Based on available data, the classification criteria are not met.

##### STOT-repeated exposure

Based on available data, the classification criteria are not met.

##### Aspiration hazard

Based on available data, the classification criteria are not met.

##### Additional information on tests

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP]. Special hazards arising from the substance or mixture!

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA FIRING PASTE

Revision date: 21.01.2020

Product code: 276

Page 6 of 8

#### SECTION 12: Ecological information

##### 12.1. Toxicity

The product is not: Ecotoxic.

##### 12.2. Persistence and degradability

The product has not been tested.

##### 12.3. Bioaccumulative potential

The product has not been tested.

##### 12.4. Mobility in soil

The product has not been tested.

##### 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

The product has not been tested.

##### 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

##### 12.7. Other adverse effects

No information available.

##### **Further information**

Avoid release to the environment.

#### SECTION 13: Disposal considerations

##### 13.1. Waste treatment methods

###### **Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

###### **Contaminated packaging**

Hazardous waste according to Directive 2008/98/EC (waste framework directive). Handle contaminated packages in the same way as the substance itself. Dampen dust and place it in a properly closed receptacle and dispose of it safely.

#### SECTION 14: Transport information

##### **Land transport (ADR/RID)**

###### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

###### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

###### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

###### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

##### **Inland waterways transport (ADN)**

###### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

###### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

###### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

###### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

##### **Marine transport (IMDG)**

###### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

###### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

###### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

###### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

##### **Air transport (ICAO-TI/IATA-DGR)**

###### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA FIRING PASTE

Revision date: 21.01.2020

Product code: 276

Page 7 of 8

**14.2. UN proper shipping name:** No dangerous good in sense of this transport regulation.

**14.3. Transport hazard class(es):** No dangerous good in sense of this transport regulation.

**14.4. Packing group:** No dangerous good in sense of this transport regulation.

#### 14.5. Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

#### 14.6. Special precautions for user

No information available.

#### 14.7. Maritime transport in bulk according to IMO instruments

not applicable

### SECTION 15: Regulatory information

#### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

##### EU regulatory information

Authorisations (REACH, annex XIV):

Substances of very high concern, SVHC (REACH, article 59):

Aluminosilicate fibre, dioxosilane

Information according to 2012/18/EU (SEVESO III): Not subject to 2012/18/EU (SEVESO III)

##### National regulatory information

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC). Observe employment restrictions under the Maternity Protection Directive (92/85/EEC) for expectant or nursing mothers.

Water hazard class (D): 3 - highly hazardous to water

#### 15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

### SECTION 16: Other information

#### Changes

This data sheet contains changes from the previous version in section(s): 1,2,4,6,7,9,11.

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA FIRING PASTE

Revision date: 21.01.2020

Product code: 276

Page 8 of 8

LL50: Lethal loading, 50%  
 EL50: Effect loading, 50%  
 EC50: Effective Concentration 50%  
 ErC50: Effective Concentration 50%, growth rate  
 NOEC: No Observed Effect Concentration  
 BCF: Bio-concentration factor  
 PBT: persistent, bioaccumulative, toxic  
 vPvB: very persistent, very bioaccumulative  
 RID: Regulations concerning the international carriage of dangerous goods by rail  
 ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
 EmS: Emergency Schedules  
 MFAG: Medical First Aid Guide  
 ICAO: International Civil Aviation Organization  
 MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
 IBC: Intermediate Bulk Container  
 SVHC: Substance of Very High Concern  
 For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Classification for mixtures and used evaluation method according to regulation (EC) No 1272/2008 [CLP]

Classification	Classification procedure
Skin Irrit. 2; H315	Calculation method
Carc. 1B; H350i	Calculation method

#### Relevant H and EUH statements (number and full text)

H315 Causes skin irritation.  
 H350i May cause cancer by inhalation.

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*