

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 1 of 8

## 1. Identification

### Product identifier

VITA VM LC flow

### Recommended use of the chemical and restrictions on use

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H.Rauter GmbH &amp; Co.KG

Street: Spitalgasse 3

Place: D-79713 Bad Säckingen

Post-office box: 1338

D-79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Contact person: regulatory affairs

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

Responsible Department: Regulatory Affairs

### Emergency phone number:

+49-(0)761-19240

### Further Information

medical device

## 2. Hazard(s) identification

### Classification of the chemical

#### 29 CFR Part 1910.1200

Respiratory or skin sensitization: Skin Sens. 1

### Label elements

#### 29 CFR Part 1910.1200

Signal word: Warning

Pictograms:



### Hazard statements

May cause an allergic skin reaction

### Precautionary statements

Avoid breathing dust/fume/gas/mist/vapors/spray.

Contaminated work clothing must not be allowed out of the workplace.

Wear protective gloves/protective clothing/eye protection/face protection.

If on skin: Wash with plenty of water.

### Hazards not otherwise classified

No information available.

## 3. Composition/information on ingredients

### Mixtures

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 2 of 8

#### Hazardous components

CAS No	Components	Quantity
72869-86-4	7,7,9-Trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)	17.55 %
109-16-0	2,2'-ethylenedioxydiethyl dimethacrylate	16.3 %
94108-97-1	Ditrimethylolpropane Tetraacrylate	1.222 %
2867-47-2	2-dimethylaminoethyl methacrylate	0.223 %

#### 4. First-aid measures

##### Description of first aid measures

###### **After inhalation**

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

###### **After contact with skin**

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

###### **After contact with eyes**

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

###### **After ingestion**

Rinse mouth immediately and drink plenty of water.

##### Most important symptoms and effects, both acute and delayed

No information available.

##### Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

#### 5. Fire-fighting measures

##### Extinguishing media

###### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

##### Specific hazards arising from the chemical

Non-flammable.

##### Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

##### Additional information

Suppress gases/vapors/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

#### 6. Accidental release measures

##### Personal precautions, protective equipment and emergency procedures

###### **General advice**

Provide adequate ventilation. Do not breathe gas/fume/vapor/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

##### Environmental precautions

Do not allow to enter into surface water or drains.

##### Methods and material for containment and cleaning up

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 3 of 8

#### Other information

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment (PPE): see section 8

Disposal: see section 13

## 7. Handling and storage

### Precautions for safe handling

#### Advice on safe handling

No special measures are necessary.

#### Advice on protection against fire and explosion

No special fire protection measures are necessary.

#### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

### Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed.

#### Hints on joint storage

No special measures are necessary.

## 8. Exposure controls/personal protection

### Control parameters

#### Exposure limits

CAS No	Substance	ppm	mg/m <sup>3</sup>	f/cc	Category	Origin
7631-86-9	Silica, amorphous	-	6		TWA (8 h)	REL

### Exposure controls



#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatril P NBR (Nitrile rubber) Breakthrough time: 10 min

##### Skin protection

Use of protective clothing.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 4 of 8

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Fresh air (open windows and doors) is necessary.

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state: viscous  
 Color:  
 Odor: characteristic

#### Changes in the physical state

Melting point/freezing point: not determined  
 Boiling point or initial boiling point and boiling range: > 283 °C  
 Flash point: 151 °C

#### Flammability

Solid/liquid: not determined  
 Gas: not applicable

#### Explosive properties

The product is not: Explosive.

#### Self-ignition temperature

Solid: not determined  
 Gas: not applicable  
 Decomposition temperature: not determined  
 pH-Value: not determined  
 Water solubility: No

#### Solubility in other solvents

not determined

Partition coefficient n-octanol/water: not determined  
 Vapor pressure: <=1100 hPa  
 (at 50 °C)  
 Density: not determined  
 Relative vapour density: not determined

### Other information

#### Information with regard to physical hazard classes

Oxidizing properties  
 Not oxidising.

#### Other safety characteristics

Solid content: 62,84 %  
 Evaporation rate: not determined

#### Further Information

## 10. Stability and reactivity

### Reactivity

No hazardous reaction when handled and stored according to provisions.

### Chemical stability

The product is stable under storage at normal ambient temperatures.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 5 of 8

#### Possibility of hazardous reactions

No known hazardous reactions.

#### Conditions to avoid

none

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

### Information on toxicological effects

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

CAS No	Components				
	Exposure route	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE 500 mg/kg			
	dermal	ATE 1100 mg/kg			

#### **Irritation and corrosivity**

Based on available data, the classification criteria are not met.

#### **Sensitizing effects**

May cause an allergic skin reaction

(7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers);

2,2'-ethylenedioxydiethyl dimethacrylate; 2-dimethylaminoethyl methacrylate)

#### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

#### **Specific target organ toxicity (STOT) - single exposure**

Based on available data, the classification criteria are not met.

#### **Specific target organ toxicity (STOT) - repeated exposure**

Based on available data, the classification criteria are not met.

Carcinogenicity (IARC): Silica, amorphous (CAS 7631-86-9) is listed in group 3.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### **Additional information on tests**

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

Harmful to aquatic life with long lasting effects.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 6 of 8

#### Mobility in soil

The product has not been tested.

#### Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### Other adverse effects

No information available.

#### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

### 13. Disposal considerations

#### Waste treatment methods

##### **Disposal recommendations**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### **Contaminated packaging**

Handle contaminated packages in the same way as the substance itself.

### 14. Transport information

#### **Marine transport (IMDG)**

<b><u>UN number or ID number:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>UN proper shipping name:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>Transport hazard class(es):</u></b>	No dangerous good in sense of this transport regulation.
<b><u>Packing group:</u></b>	No dangerous good in sense of this transport regulation.

#### **Air transport (ICAO-TI/IATA-DGR)**

<b><u>UN number or ID number:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>UN proper shipping name:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>Transport hazard class(es):</u></b>	No dangerous good in sense of this transport regulation.
<b><u>Packing group:</u></b>	No dangerous good in sense of this transport regulation.

#### **Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

#### **Special precautions for user**

No information available.

#### **Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code**

not applicable

### 15. Regulatory information

#### **U.S. Regulations**

##### **National regulatory information**

SARA Section 311/312 Hazards:

- 7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers) (72869-86-4): Immediate (acute) health hazard
- 2,2'-ethylenedioxydiethyl dimethacrylate (109-16-0): Immediate (acute) health hazard
- Ditrimethylolpropane Tetraacrylate (94108-97-1): Immediate (acute) health hazard
- 2-dimethylaminoethyl methacrylate (2867-47-2): Immediate (acute) health hazard

#### **State Regulations**

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 7 of 8

#### Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)

This product can not expose you to chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

### 16. Other information

#### Changes

Revision date:

Revision No: 2

This data sheet contains changes from the previous version in section(s): 1,2,4,6,7,8,9,11,12,13,14,16.

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%  
CLP: Classification, labelling and Packaging  
REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
EmS: Emergency Schedules  
MFAG: Medical First Aid Guide  
ICAO: International Civil Aviation Organization  
MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
IBC: Intermediate Bulk Container  
VOC: Volatile Organic Compounds  
SVHC: Substance of Very High Concern  
For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Other data

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 8 of 8

adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*