

## Safety Data Sheet

according to WHMIS

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 1 of 7

## 1. Identification

### Product identifier

VITA VM LC flow

### Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H.Rauter GmbH & Co.KG  
Street: Spitalgasse 3  
Place: D-79713 Bad Säckingen  
Post-office box: 1338  
D-79704 Bad Säckingen  
Telephone: +49(0)7761-562-0  
e-mail: info@vita-zahnfabrik.com  
Contact person: regulatory affairs  
e-mail: info@vita-zahnfabrik.com  
Internet: www.vita-zahnfabrik.com  
Responsible Department: Regulatory Affairs

Telefax: +49(0)7761-562-299

**Emergency telephone number:** +49-(0)761-19240

### Further Information

medical device

## 2. Hazard identification

### Classification of the substance or mixture

#### WHMIS 2015

Respiratory or skin sensitization: Skin Sens. 1

### Label elements

#### WHMIS 2015

**Signal word:** Warning**Pictograms:**

#### **Hazard statements**

May cause an allergic skin reaction.

#### **Precautionary statements**

Avoid breathing dust/fume/gas/mist/vapours/spray.  
Contaminated work clothing should not be allowed out of the workplace.  
Wear protective gloves/protective clothing/eye protection/face protection.  
IF ON SKIN: Wash with plenty of water.

### Other hazards

No information available.

## 3. Composition/information on ingredients

### Mixtures

## Safety Data Sheet

according to WHMIS

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 2 of 7

#### Hazardous components

CAS No	Chemical name	Quantity
72869-86-4	7,7,9-Trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)	10 - < 30% (*)
109-16-0	2,2'-ethylenedioxydiethyl dimethacrylate	10 - < 30% (*)
94108-97-1	Ditrimethylolpropane Tetraacrylate	0.5 - < 1.5% (*)
2867-47-2	2-dimethylaminoethyl methacrylate	0.1 - < 1% (*)

(\*) The actual concentration is withheld as a trade secret.

## 4. First-aid measures

#### Description of first aid measures

##### After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

##### After contact with skin

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

##### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

##### After ingestion

Rinse mouth immediately and drink plenty of water.

#### Most important symptoms and effects, whether acute or delayed

No information available.

#### Indication of immediate medical attention and special treatment needed

Treat symptomatically.

## 5. Fire-fighting measures

#### Extinguishing media

##### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

#### Specific hazards arising from the hazardous product

Non-flammable.

#### Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

#### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

## 6. Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

##### General advice

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

#### Environmental precautions

Do not allow to enter into surface water or drains.

## Safety Data Sheet

according to WHMIS

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 3 of 7

#### Methods and material for containment and cleaning up

##### Other information

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## 7. Handling and storage

#### Precautions for safe handling

##### Advice on safe handling

No special measures are necessary.

##### Advice on protection against fire and explosion

No special fire protection measures are necessary.

##### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

#### Conditions for safe storage, including any incompatibilities

##### Requirements for storage rooms and vessels

Keep container tightly closed.

##### Hints on joint storage

No special measures are necessary.

## 8. Exposure controls/Personal protection

#### Control parameters

#### Exposure controls



#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatril P NBR (Nitrile rubber) Breakthrough time: 10 min

##### Skin protection

Use of protective clothing.

##### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Fresh air (open windows and doors) is necessary.

## Safety Data Sheet

according to WHMIS

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 4 of 7

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state: viscous  
 Colour:  
 Odour: characteristic

#### Changes in the physical state

Melting point/freezing point: not determined  
 Boiling point or initial boiling point and boiling range: > 283 °C  
 Flash point: 151 °C

#### Flammability

Solid/liquid: not determined  
 Gas: not applicable

#### Explosive properties

The product is not: Explosive.

#### Self-ignition temperature

Solid: not determined  
 Gas: not applicable

Decomposition temperature: not determined

pH-Value: not determined

Water solubility: No

#### Solubility in other solvents

not determined

Partition coefficient n-octanol/water: not determined

Vapour pressure: <=1100 hPa  
 (at 50 °C)

Density: not determined

Relative vapour density: not determined

### Other information

#### Information with regard to physical hazard classes

Oxidizing properties  
 Not oxidising.

#### Other safety characteristics

Solid content: 62,84 %  
 Evaporation rate: not determined

#### Further Information

## 10. Stability and reactivity

### Reactivity

No hazardous reaction when handled and stored according to provisions.

### Chemical stability

The product is stable under storage at normal ambient temperatures.

### Possibility of hazardous reactions

No known hazardous reactions.

## Safety Data Sheet

according to WHMIS

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 5 of 7

#### Conditions to avoid

none

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

### Information on toxicological effects

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE 500 mg/kg			
	dermal	ATE 1100 mg/kg			

#### **Irritation and corrosivity**

Based on available data, the classification criteria are not met.

#### **Sensitizing effects**

May cause an allergic skin reaction.

(7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers); 2,2'-ethylenedioxydiethyl dimethacrylate; 2-dimethylaminoethyl methacrylate)

#### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

#### **STOT-single exposure**

Based on available data, the classification criteria are not met.

#### **STOT-repeated exposure**

Based on available data, the classification criteria are not met.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### **Additional information on tests**

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### **Ecotoxicity**

Harmful to aquatic life with long lasting effects.

#### **Persistence and degradability**

The product has not been tested.

#### **Bioaccumulative potential**

The product has not been tested.

#### **Mobility in soil**

The product has not been tested.

#### **Endocrine disrupting properties**

## Safety Data Sheet

according to WHMIS

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 6 of 7

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### Other adverse effects

No information available.

#### **Further information**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

### 13. Disposal considerations

#### Waste treatment methods

##### **Disposal recommendations**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### **Contaminated packaging**

Handle contaminated packages in the same way as the substance itself.

### 14. Transport information

#### Hazard classes:

##### **Marine transport (IMDG)**

##### UN number or ID number:

No dangerous good in sense of this transport regulation.

##### United Nations proper shipping

No dangerous good in sense of this transport regulation.

##### name:

##### Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### Packing group:

No dangerous good in sense of this transport regulation.

##### **Air transport (ICAO-TI/IATA-DGR)**

##### UN number or ID number:

No dangerous good in sense of this transport regulation.

##### United Nations proper shipping

No dangerous good in sense of this transport regulation.

##### name:

##### Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### Packing group:

No dangerous good in sense of this transport regulation.

#### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

### 15. Regulatory information

#### Canadian regulations

### 16. Other information

#### **Changes**

This data sheet contains changes from the previous version in section(s): 1,2,4,6,7,8,9,11,12,13,14,16.

#### **Abbreviations and acronyms**

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

## Safety Data Sheet

according to WHMIS

### VITA VM LC flow

Revision date: 13.08.2019

Product code: 284

Page 7 of 7

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LL50: Lethal loading, 50%

EL50: Effect loading, 50%

EC50: Effective Concentration 50%

ErC50: Effective Concentration 50%, growth rate

NOEC: No Observed Effect Concentration

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative

RID: Regulations concerning the international carriage of dangerous goods by rail

ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

VOC: Volatile Organic Compounds

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*