


# The comprehensive polishing solution for all ceramics



## VITA CERAMICS Polishing Set

- Simple and quick polishing
- For all types of ceramics
- Better performance, longer usage



**Find out more**

# VITA CERAMICS Polishing Set

## Polishing instruments for all types of ceramics

Especially adapted to two-stage polishing (pre-polishing and high-gloss polishing) of all dental ceramic materials. Particularly efficient and durable, due to a high diamond content in the pre-polishers, and more usable volume. For the use in the dental practice and laboratory.

### VITA CERAMICS Polishing Sets are suitable for:

- **Feldspar ceramic**, such as VITABLOCKS, VITA LUMEX AC
- **Glass ceramic**, such as VITA SUPRINITY PC, VITA AMBRIA
- **Oxide ceramic**, such as VITA YZ

### Technical Set

for the handpiece (extraoral)



Pre-polishing

High-gloss polishing

### Clinical Set

for the contra-angle (intraoral)

sterilizable



Pre-polishing

High-gloss polishing

### Product lines:

**VITA CERAMICS Polishing Set technical:** Prod. No. ECEPSETT

**VITA CERAMICS Polishing Set clinical:** Prod. No. ECEPSETC

+ Refill packs for each instrument

### Related product:

**VITA Polish Cera** diamond polishing paste for extraoral high-gloss polishing



Sales:

**VITA Zahnfabrik H. Rauter GmbH & Co. KG**

Spitalgasse 3

79713 Bad Säckingen

Germany

Phone: +49 7761 562-0

Hotline: +49 7761 562-222

E-mail: info@vita-zahnfabrik.com



**EVE Ernst Vetter GmbH**

Neureutstr. 6

75210 Kelttern

Germany

EVE Ernst Vetter GmbH has been certified and the following products bear the mark:

**CE 0483 MD**

VITA CERAMICS Polishing Set clinical and the corresponding refill packs.