

## Safety Data Sheet

according to WHMIS

### VITA OPAQUE FLUID

Revision date: 10.07.2023

Product code: 001

Page 1 of 8

## 1. Identification

### Product identifier

VITA OPAQUE FLUID

### Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H.Rauter GmbH & Co.KG  
Street: Spitalgasse 3  
Place: D-79713 Bad Säckingen  
Post-office box: 1338  
D-79704 Bad Säckingen  
Telephone: +49(0)7761-562-0  
E-mail: info@vita-zahnfabrik.com  
Contact person: regulatory affairs  
E-mail: info@vita-zahnfabrik.com  
Internet: www.vita-zahnfabrik.com  
Responsible Department: Regulatory Affairs

Telefax: +49(0)7761-562-299

### Further Information

medical device

## 2. Hazard identification

### Classification of the substance or mixture

#### WHMIS 2015

Skin corrosion/irritation: Skin Corr. 1

Serious eye damage/eye irritation: Eye Dam. 1

### Label elements

#### WHMIS 2015

#### Signal word:

Danger

#### Pictograms:



### Hazard statements

Causes severe skin burns and eye damage.

### Precautionary statements

Do not breathe dust/fume/gas/mist/vapours/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Immediately call a POISON CENTER/doctor.

### Other hazards

No information available.

## 3. Composition/information on ingredients

## Safety Data Sheet

according to WHMIS

### VITA OPAQUE FLUID

Revision date: 10.07.2023

Product code: 001

Page 2 of 8

#### Mixtures

##### Hazardous components

CAS No	Chemical name	Quantity
64-17-5	ethyl alcohol	1 - < 5% (*)

(\*) The actual concentration is withheld as a trade secret.

#### 4. First-aid measures

##### Description of first aid measures

###### General information

First aider: Pay attention to self-protection! Remove affected person from the danger area and lay down.

###### After inhalation

Provide fresh air. Medical treatment necessary.

###### After contact with skin

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention.

###### After contact with eyes

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.

###### After ingestion

Rinse mouth immediately and drink plenty of water. Do NOT induce vomiting. Adverse human health effects and symptoms: Gastric perforation. Call a physician immediately. Do not allow a neutralisation agent to be drunk.

##### Most important symptoms and effects, whether acute or delayed

No information available.

##### Indication of immediate medical attention and special treatment needed

Treat symptomatically.

#### 5. Fire-fighting measures

##### Extinguishing media

###### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

##### Specific hazards arising from the hazardous product

Non-flammable.

##### Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

##### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

#### 6. Accidental release measures

##### Personal precautions, protective equipment and emergency procedures

###### General advice

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

##### Environmental precautions

Do not allow to enter into surface water or drains.

##### Methods and material for containment and cleaning up

## Safety Data Sheet

according to WHMIS

### VITA OPAQUE FLUID

Revision date: 10.07.2023

Product code: 001

Page 3 of 8

#### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## 7. Handling and storage

### Precautions for safe handling

#### Advice on safe handling

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Advice on protection against fire and explosion

No special fire protection measures are necessary.

#### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

### Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Provide adequate ventilation as well as local exhaustion at critical locations.

#### Hints on joint storage

No special measures are necessary.

## 8. Exposure controls/Personal protection

### Control parameters

### Exposure controls



#### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Suitable eye protection: goggles.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL, Dermatrill P NBR (Nitrile rubber) Breakthrough time: 480 min

## Safety Data Sheet

according to WHMIS

### VITA OPAQUE FLUID

Revision date: 10.07.2023

Product code: 001

Page 4 of 8

#### Skin protection

Use of protective clothing.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Technical ventilation of workplace Open windows to ensure natural ventilation.

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state:	Liquid	
Colour:	colourless	
Odour:	characteristic	
Melting point/freezing point:		not determined
Boiling point or initial boiling point and boiling range:		100 °C
Flammability:		not applicable
		not applicable
Lower explosive limits:		not determined
Upper explosive limits:		not determined
Flash point:		not determined
Decomposition temperature:		not determined
pH-Value:		12,7
Solubility in other solvents		
not determined		
Partition coefficient n-octanol/water:		not determined
Vapour pressure:		<=1100 hPa
(at 50 °C)		
Density:		1,0000 g/cm <sup>3</sup>
Relative vapour density:		not determined

### Other information

#### Information with regard to physical hazard classes

##### Explosive properties

The product is not: Explosive.

##### Self-ignition temperature

Solid:

not applicable

Gas:

not applicable

##### Oxidizing properties

Not oxidising.

#### Other safety characteristics

Evaporation rate:

not determined

Solid content:

0,4 %

## 10. Stability and reactivity

### Reactivity

Possibility of hazardous reactions.

### Chemical stability

The product is stable under storage at normal ambient temperatures.

### Possibility of hazardous reactions

Exothermic reaction with: Acid, Peroxides, Oxidizing agent.

## Safety Data Sheet

according to WHMIS

### VITA OPAQUE FLUID

Revision date: 10.07.2023

Product code: 001

Page 5 of 8

#### Conditions to avoid

none

#### Incompatible materials

Keep away from: Acid, Oxidizing agent, Peroxides.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

### Information on toxicological effects

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

#### **ATEmix calculated**

ATE (oral) > 2000 mg/kg; ATE (dermal) > 2000 mg/kg; ATE (inhalation vapour) > 20 mg/l; ATE (inhalation dust/mist) > 5 mg/l

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
64-17-5	ethyl alcohol				
	dermal	LD50 7060 mg/kg			

#### **Irritation and corrosivity**

Causes severe skin burns and eye damage. (On basis of test data)

Causes serious eye damage. (On basis of test data)

#### **Sensitizing effects**

Based on available data, the classification criteria are not met.

#### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

#### **STOT-single exposure**

Based on available data, the classification criteria are not met.

#### **STOT-repeated exposure**

Based on available data, the classification criteria are not met.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### **Additional information on tests**

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

The product is not: Ecotoxic.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

#### Mobility in soil

The product has not been tested.

#### Endocrine disrupting properties

## Safety Data Sheet

according to WHMIS

### VITA OPAQUE FLUID

Revision date: 10.07.2023

Product code: 001

Page 6 of 8

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### Other adverse effects

No information available.

#### **Further information**

Avoid release to the environment.

### 13. Disposal considerations

#### Waste treatment methods

##### **Disposal recommendations**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation. Handle contaminated packages in the same way as the substance itself.

##### **Contaminated packaging**

Wash with plenty of water. Completely emptied packages can be recycled.

### 14. Transport information

#### **Canadian TDG**

<b><u>UN number:</u></b>	UN 1824
<b><u>Proper shipping name:</u></b>	SODIUM HYDROXIDE SOLUTION
<b><u>Hazard classes:</u></b>	8
<b><u>Packing group:</u></b>	III
Hazard label:	8
Limited quantity:	5 L



#### **Marine transport (IMDG)**

<b><u>UN number or ID number:</u></b>	UN 1824
<b><u>United Nations proper shipping name:</u></b>	SODIUM HYDROXIDE SOLUTION
<b><u>Transport hazard class(es):</u></b>	8
<b><u>Packing group:</u></b>	III
Hazard label:	8



Special Provisions:	223
Limited quantity:	5 L
Excepted quantity:	E1
EmS:	F-A, S-B
Segregation group:	18 - alkalis

#### **Air transport (ICAO-TI/IATA-DGR)**

<b><u>UN number or ID number:</u></b>	UN 1824
<b><u>United Nations proper shipping name:</u></b>	SODIUM HYDROXIDE SOLUTION
<b><u>Transport hazard class(es):</u></b>	8
<b><u>Packing group:</u></b>	III
Hazard label:	8

## Safety Data Sheet

according to WHMIS

### VITA OPAQUE FLUID

Revision date: 10.07.2023

Product code: 001

Page 7 of 8



Special Provisions:	A3 A803	
Limited quantity Passenger:	1 L	
Passenger LQ:	Y841	
Excepted quantity:	E1	
IATA-packing instructions - Passenger:		852
IATA-max. quantity - Passenger:		5 L
IATA-packing instructions - Cargo:		856
IATA-max. quantity - Cargo:		60 L

#### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

#### 15. Regulatory information

##### Canadian regulations

#### 16. Other information

##### Changes

This data sheet contains changes from the previous version in section(s): 1.

## Safety Data Sheet

according to WHMIS

### VITA OPAQUE FLUID

Revision date: 10.07.2023

Product code: 001

Page 8 of 8

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LL50: Lethal loading, 50%

EL50: Effect loading, 50%

EC50: Effective Concentration 50%

ErC50: Effective Concentration 50%, growth rate

NOEC: No Observed Effect Concentration

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative

RID: Regulations concerning the international carriage of dangerous goods by rail

ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*