

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA Die Release

Revision date: 10.07.2023

Product code: 016

Page 1 of 8

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

VITA Die Release

#### Further trade names

VITA Stumpfisolierung

CAS No: 8042-47-5  
EC No: 232-455-8  
UFI: X300-P0FG-800R-GA83

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### 1.3. Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H.Rauter GmbH & Co.KG  
Street: Spitalgasse 3  
Place: D-79713 Bad Säckingen  
Post-office box: 1338  
D-79704 Bad Säckingen  
Telephone: +49(0)7761-562-0  
E-mail: info@vita-zahnfabrik.com  
Contact person: regulatory affairs  
E-mail: info@vita-zahnfabrik.com  
Internet: www.vita-zahnfabrik.com  
Responsible Department: Regulatory Affairs  
Telefax: +49(0)7761-562-299

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### Regulation (EC) No 1272/2008

Asp. Tox. 1; H304

Full text of hazard statements: see SECTION 16.

### 2.2. Label elements

#### Regulation (EC) No 1272/2008

Signal word: Danger

Pictograms:



#### Hazard statements

H304 May be fatal if swallowed and enters airways.

#### Precautionary statements

P301+P310 IF SWALLOWED: Immediately call a POISON CENTER/doctor.  
P331 Do NOT induce vomiting.

### 2.3. Other hazards

No information available.

## SECTION 3: Composition/information on ingredients

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA Die Release

Revision date: 10.07.2023

Product code: 016

Page 2 of 8

#### 3.1. Substances

##### Chemical characterization

Substance, organic

##### Hazardous components

| CAS No    | Chemical name                                 |          |                       | Quantity |
|-----------|---|----------|-----------------------|----------|
|           | EC No   | Index No | REACH No              |          |
|           | Classification (Regulation (EC) No 1272/2008) |          |                       |          |
| 8042-47-5 | Paraffin oils                                 |          |                       | 100 %    |
|           | 232-455-8                                     |          | 01-2119487078-27-XXXX |          |
|           | Asp. Tox. 1; H304                             |          |                       |          |

Full text of H and EUH statements: see section 16.

##### Specific Conc. Limits, M-factors and ATE

| CAS No    | EC No                                    | Chemical name | Quantity |
|-----------|--|---------------|----------|
|           | Specific Conc. Limits, M-factors and ATE |               |          |
| 8042-47-5 | 232-455-8                                | Paraffin oils | 100 %    |
|           | oral: LD50 = > 5000 mg/kg                |               |          |

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

##### After inhalation

Provide fresh air. Medical treatment necessary.

##### After contact with skin

Wash with plenty of water. Immediately remove any contaminated clothing, shoes or stockings. Medical treatment necessary.

##### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

##### After ingestion

Observe risk of aspiration if vomiting occurs. Rinse mouth immediately and drink plenty of water. Induce vomiting when the affected person is not unconscious. Medical treatment necessary.

#### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

#### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

##### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

#### 5.2. Special hazards arising from the substance or mixture

Non-flammable.

#### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus.

#### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA Die Release

Revision date: 10.07.2023

Product code: 016

Page 3 of 8

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### General advice

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

### 6.3. Methods and material for containment and cleaning up

#### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advice on safe handling

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Advice on protection against fire and explosion

No special fire protection measures are necessary.

#### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

### 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed. Keep locked up. Store in a place accessible by authorized persons only. Provide adequate ventilation as well as local exhaust at critical locations.

#### Hints on joint storage

No special measures are necessary.

### 7.3. Specific end use(s)

Use as laboratory reagent

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

### 8.2. Exposure controls



#### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA Die Release

Revision date: 10.07.2023

Product code: 016

Page 4 of 8

gas/fumes/vapour/spray.

#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles NBR (Nitrile rubber) KCL Dermatril P Breakthrough time: 480 min

##### Skin protection

Use of protective clothing.

##### Respiratory protection

Technical ventilation of workplace Open windows to ensure natural ventilation.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

|   |            |                                  |
|---|------------|----------------------------------|
| Physical state:   | Liquid     |                                  |
| Colour:   | colourless |                                  |
| Odour:  | odourless  |                                  |
| Melting point/freezing point:                             |            | not determined                   |
| Boiling point or initial boiling point and boiling range: |            | > 300 °C                         |
| Flammability:   |            | not applicable<br>not applicable |
| Flash point:  |            | 195 °C                           |
| Decomposition temperature:                                |            | not determined                   |
| pH-Value:   |            | not determined                   |
| Viscosity / kinematic:<br>(at 40 °C)                      |            | 16 mm <sup>2</sup> /s            |
| Water solubility:   |            | No                               |
| Solubility in other solvents<br>not determined            |            |                                  |
| Partition coefficient n-octanol/water:                    |            | not determined                   |
| Vapour pressure:  |            | not determined                   |
| Density:  |            | 0,85000 g/cm <sup>3</sup>        |
| Relative vapour density:                                  |            | not determined                   |

### 9.2. Other information

#### Information with regard to physical hazard classes

##### Explosive properties

The product is not: Explosive.

##### Self-ignition temperature

Solid:

not applicable

Gas:

not applicable

##### Oxidizing properties

Not oxidising.

#### Other safety characteristics

Evaporation rate:

not determined

Solid content:

0,0 %

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA Die Release

Revision date: 10.07.2023

Product code: 016

Page 5 of 8

Viscosity / dynamic:

35 mPa·s

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

### 10.4. Conditions to avoid

none

### 10.5. Incompatible materials

No information available.

### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

## SECTION 11: Toxicological information

### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

#### Acute toxicity

Based on available data, the classification criteria are not met.

| CAS No    | Chemical name  |                   |         |           |        |
|-----------|----------------|-------------------|---------|-----------|--------|
|           | Exposure route | Dose              | Species | Source    | Method |
| 8042-47-5 | Paraffin oils  |                   |         |           |        |
|           | oral           | LD50 > 5000 mg/kg | Ratte   | SDB Merck |        |

#### Irritation and corrosivity

Based on available data, the classification criteria are not met.

#### Sensitising effects

Based on available data, the classification criteria are not met.

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### STOT-single exposure

Based on available data, the classification criteria are not met.

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

May be fatal if swallowed and enters airways.

#### Additional information on tests

Special hazards arising from the substance or mixture. Classification according to Regulation (EC) No 1272/2008

[CLP]: Health hazard properties.

## SECTION 12: Ecological information

### 12.1. Toxicity

The product is not: Ecotoxic.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA Die Release

Revision date: 10.07.2023

Product code: 016

Page 6 of 8

| CAS No    | Chemical name       |              |           |         |                |           |
|-----------|---------------------|--------------|-----------|---------|----------------|-----------|
|           | Aquatic toxicity    | Dose         | [h]   [d] | Species | Source         | Method    |
| 8042-47-5 | Paraffin oils       |              |           |         |                |           |
|           | Acute fish toxicity | LC50<br>mg/l | >1000     | 96 h    | Leuciscus idus | SDB Merck |

#### 12.2. Persistence and degradability

The product has not been tested.

#### 12.3. Bioaccumulative potential

The product has not been tested.

#### 12.4. Mobility in soil

The product has not been tested.

#### 12.5. Results of PBT and vPvB assessment

This substance does not meet the PBT/vPvB criteria of REACH, annex XIII.

The product has not been tested.

#### 12.6. Endocrine disrupting properties

This substance does not have endocrine disrupting properties with respect to non-target organisms.

#### 12.7. Other adverse effects

No information available.

#### Further information

Avoid release to the environment.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Hazardous waste according to Directive 2008/98/EC (waste framework directive). Handle contaminated packages in the same way as the substance itself.

### SECTION 14: Transport information

#### Land transport (ADR/RID)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Inland waterways transport (ADN)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Marine transport (IMDG)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Air transport (ICAO-TI/IATA-DGR)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA Die Release

Revision date: 10.07.2023

Product code: 016

Page 7 of 8

**14.2. UN proper shipping name:** No dangerous good in sense of this transport regulation.

**14.3. Transport hazard class(es):** No dangerous good in sense of this transport regulation.

**14.4. Packing group:** No dangerous good in sense of this transport regulation.

#### **14.5. Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

#### **14.6. Special precautions for user**

No information available.

#### **14.7. Maritime transport in bulk according to IMO instruments**

not applicable

### SECTION 15: Regulatory information

#### **15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

##### **EU regulatory information**

Restrictions on use (REACH, annex XVII):

Entry 3

Information according to 2012/18/EU (SEVESO III): Not subject to 2012/18/EU (SEVESO III)

##### **National regulatory information**

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC).

Water hazard class (D): 1 - slightly hazardous to water

#### **15.2. Chemical safety assessment**

For this substance a chemical safety assessment has not been carried out.

### SECTION 16: Other information

#### **Changes**

This data sheet contains changes from the previous version in section(s): 1,2,8,9.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### VITA Die Release

Revision date: 10.07.2023

Product code: 016

Page 8 of 8

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LL50: Lethal loading, 50%

EL50: Effect loading, 50%

EC50: Effective Concentration 50%

ErC50: Effective Concentration 50%, growth rate

NOEC: No Observed Effect Concentration

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative

RID: Regulations concerning the international carriage of dangerous goods by rail

ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

Asp. Tox: Aspiration hazard

#### Relevant H and EUH statements (number and full text)

H304                      May be fatal if swallowed and enters airways.

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.