

## Safety Data Sheet

according to WHMIS

### VITAVM LC BASE DENTINE\_ENAMEL\_EFFECT-ENAMEL\_NEUTRAL\_GINGIVA

Revision date: 10.07.2023

Product code: 161

Page 1 of 7

## 1. Identification

### Product identifier

VITAVM LC BASE DENTINE\_ENAMEL\_EFFECT-ENAMEL\_NEUTRAL\_GINGIVA

### Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name:	VITA Zahnfabrik H.Rauter GmbH & Co.KG	
Street:	Spitalgasse 3	
Place:	D-79713 Bad Säckingen	
Post-office box:	1338	
	D-79704 Bad Säckingen	
Telephone:	+49(0)7761-562-0	Telefax: +49(0)7761-562-299
E-mail:	info@vita-zahnfabrik.com	
Contact person:	regulatory affairs	
E-mail:	info@vita-zahnfabrik.com	
Internet:	www.vita-zahnfabrik.com	
Responsible Department:	Regulatory Affairs	

### Further Information

medical device

## 2. Hazard identification

### Classification of the substance or mixture

#### WHMIS 2015

Respiratory or skin sensitization: Skin Sens. 1

### Label elements

#### WHMIS 2015

Signal word: Warning

Pictograms:



#### Hazard statements

May cause an allergic skin reaction.

#### Precautionary statements

Avoid breathing dust/fume/gas/mist/vapours/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

### Other hazards

No information available.

## 3. Composition/information on ingredients

### Mixtures

## Safety Data Sheet

according to WHMIS

### VITAVM LC BASE DENTINE\_ENAMEL\_EFFECT-ENAMEL\_NEUTRAL\_GINGIVA

Revision date: 10.07.2023

Product code: 161

Page 2 of 7

#### Hazardous components

CAS No	Chemical name	Quantity
72869-86-4	7,7,9-Trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)	10 - < 30% (*)
109-16-0	2,2'-ethylenedioxydiethyl dimethacrylate	10 - < 30% (*)

(\*) The actual concentration is withheld as a trade secret.

#### 4. First-aid measures

##### Description of first aid measures

###### **After inhalation**

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

###### **After contact with skin**

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary. IF ON SKIN: Wash with plenty of soap and water.

###### **After contact with eyes**

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

###### **After ingestion**

Rinse mouth immediately and drink plenty of water.

##### Most important symptoms and effects, whether acute or delayed

No information available.

##### Indication of immediate medical attention and special treatment needed

Treat symptomatically.

#### 5. Fire-fighting measures

##### Extinguishing media

###### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

##### Specific hazards arising from the hazardous product

Non-flammable.

##### Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

##### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

#### 6. Accidental release measures

##### Personal precautions, protective equipment and emergency procedures

###### **General advice**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

##### Environmental precautions

Do not allow to enter into surface water or drains.

##### Methods and material for containment and cleaning up

## Safety Data Sheet

according to WHMIS

### VITAVM LC BASE DENTINE\_ENAMEL\_EFFECT-ENAMEL\_NEUTRAL\_GINGIVA

Revision date: 10.07.2023

Product code: 161

Page 3 of 7

#### Other information

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## 7. Handling and storage

#### Precautions for safe handling

##### **Advice on safe handling**

No special measures are necessary.

##### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

##### **Advice on general occupational hygiene**

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

#### Conditions for safe storage, including any incompatibilities

##### **Requirements for storage rooms and vessels**

Keep container tightly closed.

##### **Hints on joint storage**

No special measures are necessary.

## 8. Exposure controls/Personal protection

#### Control parameters

#### Exposure controls



#### Individual protection measures, such as personal protective equipment

##### **Eye/face protection**

Suitable eye protection: goggles.

##### **Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatril P Breakthrough time: 10 min NBR (Nitrile rubber)

##### **Skin protection**

Use of protective clothing.

##### **Respiratory protection**

Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

## Safety Data Sheet

according to WHMIS

### VITAVM LC BASE DENTINE\_ENAMEL\_EFFECT-ENAMEL\_NEUTRAL\_GINGIVA

Revision date: 10.07.2023

Product code: 161

Page 4 of 7

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state:	
Colour:	
Odour:	characteristic
Melting point/freezing point:	not determined
Boiling point or initial boiling point and boiling range:	139 °C
Flammability:	not determined not applicable
Lower explosive limits:	not determined
Upper explosive limits:	not determined
Flash point:	151 °C
Decomposition temperature:	not determined
pH-Value:	not determined
Water solubility:	No
Solubility in other solvents	not determined
Partition coefficient n-octanol/water:	not determined
Vapour pressure: (at 50 °C)	<=1100 hPa
Density:	not determined
Relative vapour density:	not determined

### Other information

#### **Information with regard to physical hazard classes**

##### Explosive properties

The product is not: Explosive.

##### Self-ignition temperature

Solid:

not determined

Gas:

not applicable

##### Oxidizing properties

Not oxidising.

#### **Other safety characteristics**

Evaporation rate:

not determined

Solid content:

59,1 %

## 10. Stability and reactivity

### Reactivity

No hazardous reaction when handled and stored according to provisions.

### Chemical stability

The product is stable under storage at normal ambient temperatures.

### Possibility of hazardous reactions

No known hazardous reactions.

### Conditions to avoid

none

### Incompatible materials

No information available.

## Safety Data Sheet

according to WHMIS

### VITAVM LC BASE DENTINE\_ENAMEL\_EFFECT-ENAMEL\_NEUTRAL\_GINGIVA

Revision date: 10.07.2023

Product code: 161

Page 5 of 7

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

### Information on toxicological effects

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

#### **ATEmix calculated**

ATE (oral) &gt; 2000 mg/kg; ATE (dermal) &gt; 2000 mg/kg; ATE (inhalation vapour) &gt; 20 mg/l; ATE (inhalation dust/mist) &gt; 5 mg/l

#### **Irritation and corrosivity**

Based on available data, the classification criteria are not met.

#### **Sensitizing effects**

May cause an allergic skin reaction.

(7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers); 2,2'-ethylenedioxydiethyl dimethacrylate)

#### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

#### **STOT-single exposure**

Based on available data, the classification criteria are not met.

#### **STOT-repeated exposure**

Based on available data, the classification criteria are not met.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### **Additional information on tests**

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

### Ecotoxicity

The product is not: Ecotoxic.

### Persistence and degradability

The product has not been tested.

### Bioaccumulative potential

The product has not been tested.

### Mobility in soil

The product has not been tested.

### Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

### Other adverse effects

No information available.

### **Further information**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

## 13. Disposal considerations

### Waste treatment methods

## Safety Data Sheet

according to WHMIS

### VITAVM LC BASE DENTINE\_ENAMEL\_EFFECT-ENAMEL\_NEUTRAL\_GINGIVA

Revision date: 10.07.2023

Product code: 161

Page 6 of 7

#### Disposal recommendations

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

#### Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

### 14. Transport information

#### Marine transport (IMDG)

**UN number or ID number:** No dangerous good in sense of this transport regulation.

**United Nations proper shipping name:** No dangerous good in sense of this transport regulation.

**Transport hazard class(es):** No dangerous good in sense of this transport regulation.

**Packing group:** No dangerous good in sense of this transport regulation.

#### Air transport (ICAO-TI/IATA-DGR)

**UN number or ID number:** No dangerous good in sense of this transport regulation.

**United Nations proper shipping name:** No dangerous good in sense of this transport regulation.

**Transport hazard class(es):** No dangerous good in sense of this transport regulation.

**Packing group:** No dangerous good in sense of this transport regulation.

#### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

### 15. Regulatory information

#### Canadian regulations

### 16. Other information

#### Changes

This data sheet contains changes from the previous version in section(s): 1.

## Safety Data Sheet

according to WHMIS

### VITAVM LC BASE DENTINE\_ENAMEL\_EFFECT-ENAMEL\_NEUTRAL\_GINGIVA

Revision date: 10.07.2023

Product code: 161

Page 7 of 7

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LL50: Lethal loading, 50%

EL50: Effect loading, 50%

EC50: Effective Concentration 50%

ErC50: Effective Concentration 50%, growth rate

NOEC: No Observed Effect Concentration

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative

RID: Regulations concerning the international carriage of dangerous goods by rail

ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

VOC: Volatile Organic Compounds

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*