

VITAVM.11

Product information



VITA shade determination

VITA shade communication

VITA shade reproduction

VITA shade control

VITA – perfect match.

VITA

Fine-structure ceramic for individualizing restorations made of zirconia-reinforced lithium silicate glass ceramic (ZLS), especially for VITA SUPRINITY PC.



Fig. 1 VITA VM 11 effect materials

Fine-structure ceramic for individualizing VITA SUPRINITY PC

VITA VM 11 is a low-fusing, fine-structure feldspar ceramic that has been specially developed for individualizing restorations made of zirconia-reinforced lithium silicate glass ceramic.*

The perfectly matched CTE values of substructure and veneering materials help minimize stress to ensure superior bonding and veneering reliability that is free of warping.

Product benefits

Highly esthetic restorations:

The high translucency and warm shades of VITA VM 11, in combination with the opalescent effect of VITA SUPRINITY PC, create highly esthetic restorations with a vivid play of colors.

Reliable bonding:

Stress-free and reliable bonding is ensured through a perfect match of both CTE ranges.



Fig. 2

Good grinding and polishing properties:

Thanks to the proven fine structure of VITA VM 11, its smooth and densely sealed surface can be easily and quickly polished.



Fig. 3

Simple processing:

Excellent stability, minimal shrinkage and high edge stability are distinctive features of VITA VM 11. Thanks to the very good surface wettability of the new glass ceramic VITA SUPRINITY PC, multiple layering without liner firing or washbake is possible.



Fig. 4

Convincing firing stability:

The impressive firing properties of VITA VM 11 result in very high dimensional stability, even after several firings.

Figures 2 - 4 were provided with the kind support of Mrs. Heike Assmann (MDT).

* can also be used for CELTRA® Duo and CELTRA® CAD FC.
CELTRA is a registered trademark of DeguDent GmbH.



Assortments

VITAVM®11 DENTINE KIT 3D-MASTER, BV11DK

Quantity	Contents	Material
6	12 g	TRANSPA DENTINE 0M1, 1M1, 1M2, 2M2, 3M2, 4M2
2	12 g	ENAMEL ENL, END
1	50 ml	LOW FUSING MODELLING LIQUID
		Instructions for use

VITAVM®11 DENTINE KIT classical, BV11DKC

Quantity	Contents	Material
8	12 g	TRANSPA DENTINE A1, A2, A3, A3.5, B2, C2, D2
2	12 g	ENAMEL ENL, END
1	50 ml	LOW FUSING MODELLING LIQUID
		Instructions for use


Individual materials

VITAVM®11

Shade	Material	Content	Prod. No.
0M1	TRANSPA DENTINE	12 g	B4906112
1M1	TRANSPA DENTINE	12 g	B4906412
1M2	TRANSPA DENTINE	12 g	B4906512
2M2	TRANSPA DENTINE	12 g	B4906912
3M2	TRANSPA DENTINE	12 g	B4907612
4M2	TRANSPA DENTINE	12 g	B4908312
A1	TRANSPA DENTINE	12 g	B4945112
A2	TRANSPA DENTINE	12 g	B4945212
A3	TRANSPA DENTINE	12 g	B4945312
A3.5	TRANSPA DENTINE	12 g	B4945412
B2	TRANSPA DENTINE	12 g	B4945712
C2	TRANSPA DENTINE	12 g	B4946112
D2	TRANSPA DENTINE	12 g	B4946412
WIN	WINDOW	12 g	B4918112
NT	NEUTRAL	12 g	B4920112
ENL	ENAMEL light	12 g	B4919112
END	ENAMEL dark	12 g	B4919212
MM1	MAMELON	12 g	B4916112
MM3	MAMELON	12 g	B4916312

Shade	Material	Content	Prod. No.
EP1	EFFECT PEARL	12 g	B4915112
EO1	EFFECT OPAL	12 g	B4917112
EO2	EFFECT OPAL	12 g	B4917212
EO3	EFFECT OPAL	12 g	B4917312
EO5	EFFECT OPAL	12 g	B4917512
EC1	EFFECT CHROMA	12 g	B4912112
EC3	EFFECT CHROMA	12 g	B4912312
EC5	EFFECT CHROMA	12 g	B4912512
EE1	EFFECT ENAMEL	12 g	B4909112
EE3	EFFECT ENAMEL	12 g	B4909312
EE5	EFFECT ENAMEL	12 g	B4909512
EE7	EFFECT ENAMEL	12 g	B4909712
EE8	EFFECT ENAMEL	12 g	B4909812
EE9	EFFECT ENAMEL	12 g	B4909912
EE11	EFFECT ENAMEL	12 g	B4910112
SD1	SUN DENTINE	12 g	B4936112
SD2	SUN DENTINE	12 g	B4936212
SD3	SUN DENTINE	12 g	B4936312

With the unique VITA SYSTEM 3D-MASTER, all natural tooth shades can be systematically determined and perfectly reproduced.



Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information: 2024-04

After the publication of this information for use any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified and the following products bear the mark

CE0124

VITA^{VM}11 · VITA SUPRINITY[®] PC



VITA

 VITA Zahnfabrik H. Rauter GmbH & Co.KG
Spitalgasse 3 · 79713 Bad Säckingen · Germany
Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299
Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446
www.vita-zahnfabrik.com · info@vita-zahnfabrik.com
 facebook.com/vita.zahnfabrik