

True Strength is Multi-Layered.



VITABLOCS® TriLuxe forte

Mill. Polish. Done.

- four layers for a natural shade gradient and unique esthetics
- simple polishing for optimum gloss
- high durability, thanks to enamel-like abrasion



CAD/CAM

VITABLOCS® TriLuxe forte

The polychromatic CAD/CAM block: Mill. Polish. Done.

VITABLOCS TriLuxe forte with integrated four-layer shade gradient makes optimized reproduction of the natural play of colors possible in the esthetic zone. Chairside treatments can be performed easily and efficiently in just one session. The four-layer shade gradient also ensures highly esthetic results.

[Click here for Studies](#)

From inlays and veneers to full crowns, VITABLOCS TriLuxe forte offers you a wide range of design options to ensure the best possible care for your patients.



Shades. Geometries. System compatibility

Shades

- 5 VITA SYSTEM 3D-MASTER® shades: 0M1C, 1M1C, 1M2C, 2M2C, 3M2C
- 9 VITA classical A1–D4® shades: A1C, A2C, A3C, A3.5C, B1C, B2C, B3C, C2C, C3C

Geometries



VITABLOCS TriLuxe forte – SYSTEM SOLUTION

VITA offers VITABLOCS TriLuxe forte with the holder system for CEREC/inLab (Dentsply Sirona).

VITABLOCS TriLuxe forte – UNIVERSAL SOLUTIONS*

- CORiTEC line (imes-icore GmbH)
- DGSHAPE DWX line (DGSHAPE Corporation), CS 3000 /CS 3100 (Carestream Dental, Inc.)
- N/R/S/Z/E lines (vhf camfacture AG), Zfx Inhouse5x (Zfx GmbH)
- Organical Desktop line (R+K CAD/CAM Technologie GmbH & Co. KG)
- DMG ULTRASONIC series (DMG Mori AG)
- Ceramill DRS (Amann Girrbach)

*Note: The range of variants/geometries/shades of VITA CAD/CAM materials may vary for individual CAD/CAM system partners or systems.



As an expert in the field, I have tried every block out there. My hands down favorite is VITABLOCS TriLuxe forte. This block has impressive esthetics and perfectly imitates natural tooth enamel."

Dr. Andres Powditch, Dentist, Chile

[Click here for Webinar](#)



The integrated shade gradient in four layers of intensity facilitates a natural play of shade and light in the esthetic zone."

**Dr. Mon Li and Sally Hsieh
(Dental Technician)
CEREC Asia, Taipei, Taiwan**





VITABLOCS® TriLuxe forte

Please note:

Our products must be used according to the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of application. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Publication of this product brochure: 2024-06. On publication of this information for use, previous versions will no longer be valid.

The current version in each case is available at www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified, and the following products bear the CE mark

CE 0124: VITABLOCS® TriLuxe forte

The products/systems of other manufacturers mentioned in this document are registered trademarks of the respective manufacturers.

MD

VITA Zahnfabrik H. Rauter GmbH & Co. KG

Spitalgasse 3
79713 Bad Säckingen
Germany

Phone: +49 7761 562-0
Hotline: +49 7761 562-222

info@vita-zahnfabrik.com
www.vita-zahnfabrik.com

Follow us on
social media!

