

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.10)

### VITA VIONIC BOND I; VITA VIONIC BOND multiDose I

Revision date: 05.09.2024

Product code: 288

Page 1 of 7

## 1. Identification

### Product identifier

VITA VIONIC BOND I; VITA VIONIC BOND multiDose I

### Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name:	VITA Zahnfabrik H.Rauter GmbH & Co.KG	
Street:	Spitalgasse 3	
Place:	D-79713 Bad Säckingen	
Post-office box:	1338	
	D-79704 Bad Säckingen	
Telephone:	+49(0)7761-562-0	Telefax: +49(0)7761-562-299
E-mail:	info@vita-zahnfabrik.com	
Contact person:	regulatory affairs	
E-mail:	info@vita-zahnfabrik.com	
Internet:	www.vita-zahnfabrik.com	
Responsible Department:	Regulatory Affairs	

### Further Information

medical device

## 2. Hazard identification

### Classification of the substance or mixture

#### UN GHS (ST/SG/AC.10/11/Rev.10)

Flammable liquids: Flam. Liq. 2

Skin corrosion/irritation: Skin Irrit. 2

Respiratory or skin sensitization: Skin Sens. 1

Specific target organ toxicity - single exposure: STOT SE 3 (respiratory tract irritation)

### Label elements

#### UN GHS (ST/SG/AC.10/11/Rev.10)

#### Hazard components for labelling

2,2'-ethylenedioxydiethyl dimethacrylate

dibenzoyl peroxide; benzoyl peroxide

methyl methacrylate

Signal word: Danger

#### Pictograms:



#### Hazard statements

H225	Highly flammable liquid and vapour.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H335	May cause respiratory irritation.

#### Precautionary statements

P403+P233	Store in a well-ventilated place. Keep container tightly closed.
-----------	--

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.10)

### VITA VIONIC BOND I; VITA VIONIC BOND multiDose I

Revision date: 05.09.2024

Product code: 288

Page 2 of 7

P235	Keep cool.
P233	Keep container tightly closed.
P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P280	Wear protective gloves/protective clothing/eye protection/face protection.

#### Other hazards

No information available.

### 3. Composition/information on ingredients

#### Mixtures

##### Relevant ingredients

CAS No	Chemical name	Quantity
	Classification (UN GHS (ST/SG/AC.10/11/Rev.10))	
80-62-6	methyl methacrylate	75 - < 80 %
	Flam. Liq. 2, Skin Irrit. 2, Skin Sens. 1, STOT SE 3; H225 H315 H317 H335	
109-16-0	triethylene glycol dimethacrylate	5 - < 10 %
	Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1, STOT SE 3; H315 H319 H317 H335	
94-36-0	dibenzoyl peroxide; benzoyl peroxide	1 - < 5 %
	Org. Perox. B, Eye Irrit. 2A, Skin Sens. 1; H241 H319 H317	

### 4. First-aid measures

#### Description of first aid measures

##### After inhalation

Provide fresh air. Medical treatment necessary.

##### After contact with skin

After contact with skin, wash immediately with plenty of water and soap.

##### After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

##### After ingestion

Rinse mouth immediately and drink plenty of water.

#### Most important symptoms and effects, whether acute or delayed

No information available.

#### Indication of immediate medical attention and special treatment needed

Treat symptomatically.

### 5. Fire-fighting measures

#### Extinguishing media

##### Suitable extinguishing media

Carbon dioxide (CO<sub>2</sub>), Foam, Extinguishing powder.

##### Unsuitable extinguishing media

Water.

#### Specific hazards arising from the hazardous product

Highly flammable. Vapours can form explosive mixtures with air.

#### Special protective equipment and precautions for fire-fighters

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.10)

### VITA VIONIC BOND I; VITA VIONIC BOND multiDose I

Revision date: 05.09.2024

Product code: 288

Page 3 of 7

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

#### Additional information

Use water spray jet to protect personnel and to cool endangered containers. Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

## 6. Accidental release measures

### Personal precautions, protective equipment and emergency procedures

#### General advice

Remove all sources of ignition. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

#### Environmental precautions

Do not allow uncontrolled discharge of product into the environment. Danger of explosion

#### Methods and material for containment and cleaning up

#### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## 7. Handling and storage

### Precautions for safe handling

#### Advice on safe handling

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Advice on protection against fire and explosion

Keep away from sources of ignition - No smoking. Take precautionary measures against static discharges. Vapours can form explosive mixtures with air.

#### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

### Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed. Keep in a cool, well-ventilated place. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

#### Hints on joint storage

Do not store together with: Oxidising agent. Pyrophoric or self-heating substances.

## 8. Exposure controls/Personal protection

### Control parameters

### Exposure controls



## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.10)

### VITA VIONIC BOND I; VITA VIONIC BOND multiDose I

Revision date: 05.09.2024

Product code: 288

Page 4 of 7

#### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Wear eye/face protection.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Butoject Butyl caoutchouc (butyl rubber)  
Breakthrough time: 60 min

##### Skin protection

Wear suitable protective clothing.

##### Respiratory protection

Technical ventilation of workplace Provide adequate ventilation as well as local exhaustion at critical locations.

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state:	Liquid	
Colour:	colourless	
Odour:	stinging	
Melting point/freezing point:		not determined
Boiling point or initial boiling point and boiling range:		101 °C
Flammability:		not applicable
		not applicable
Lower explosive limits:		2,1 vol. %
Upper explosive limits:		12,5 vol. %
Flash point:		10 °C
Auto-ignition temperature:		430 °C
Decomposition temperature:		not determined
pH-Value:		not determined
Water solubility:		No
Solubility in other solvents		
not determined		
Partition coefficient n-octanol/water:		not determined
Vapour pressure:		<=1100 hPa
(at 50 °C)		
Density:		not determined
Relative vapour density:		not determined

### Other information

#### Information with regard to physical hazard classes

Self-ignition temperature

Solid: not applicable

Gas: not applicable

Oxidizing properties

Not oxidising.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.10)

### VITA VIONIC BOND I; VITA VIONIC BOND multiDose I

Revision date: 05.09.2024

Product code: 288

Page 5 of 7

#### Other safety characteristics

Evaporation rate: not determined  
Solid content: 1,0 %

### 10. Stability and reactivity

#### Reactivity

Highly flammable.

#### Chemical stability

The product is stable under storage at normal ambient temperatures.

#### Possibility of hazardous reactions

No known hazardous reactions.

#### Conditions to avoid

Keep away from sources of heat (e.g. hot surfaces), sparks and open flames. Vapours can form explosive mixtures with air.

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

### 11. Toxicological information

#### Information on toxicological effects

##### Acute toxicity

Based on available data, the classification criteria are not met.

##### ATEmix calculated

ATE (oral) > 5000 mg/kg; ATE (dermal) > 5000 mg/kg; ATE (inhalation vapour) > 50 mg/l; ATE (inhalation dust/mist) > 12.5 mg/l

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
80-62-6	methyl methacrylate				
	dermal	LD50 > 5000 mg/kg			

#### Irritation and corrosivity

Skin corrosion/irritation: Causes skin irritation.

Serious eye damage/eye irritation: Based on available data, the classification criteria are not met.

#### Sensitizing effects

May cause an allergic skin reaction. (methyl methacrylate; triethylene glycol dimethacrylate; dibenzoyl peroxide; benzoyl peroxide)

#### Carcinogenic/mutagenic/toxic effects for reproduction

Germ cell mutagenicity: Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

#### STOT-single exposure

May cause respiratory irritation. (methyl methacrylate)

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.10)

### VITA VIONIC BOND I; VITA VIONIC BOND multiDose I

Revision date: 05.09.2024

Product code: 288

Page 6 of 7

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### Additional information on tests

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

The product is not: Ecotoxic.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

#### Mobility in soil

The product has not been tested.

#### Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### Other adverse effects

No information available.

#### Further information

Avoid release to the environment.

## 13. Disposal considerations

#### Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

## 14. Transport information

#### Marine transport (IMDG)

<u>UN number or ID number:</u>	UN 1993
<u>United Nations proper shipping name:</u>	FLAMMABLE LIQUID, N.O.S. (methyl methacrylate)
<u>Transport hazard class(es):</u>	3
<u>Packing group:</u>	II
Hazard label:	3



Special Provisions:	274
Limited quantity:	1 L
Excepted quantity:	E2
EmS:	F-E, S-E
Segregation group:	ammonium compounds

#### Air transport (ICAO-TI/IATA-DGR)

<u>UN number or ID number:</u>	UN 1993
--------------------------------	---------

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.10)

### VITA VIONIC BOND I; VITA VIONIC BOND multiDose I

Revision date: 05.09.2024

Product code: 288

Page 7 of 7

**United Nations proper shipping name:**

FLAMMABLE LIQUID, N.O.S. (methyl methacrylate)

**Transport hazard class(es):**

3

**Packing group:**

II

Hazard label:

3



Special Provisions:

A3

Limited quantity Passenger:

1 L

Passenger LQ:

Y341

Excepted quantity:

E2

IATA-packing instructions - Passenger:

353

IATA-max. quantity - Passenger:

5 L

IATA-packing instructions - Cargo:

364

IATA-max. quantity - Cargo:

60 L

**Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

**Special precautions for user**

Warning: Combustible liquid.

**Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code**

not applicable

## 15. Regulatory information

**National regulatory information**

Employment restrictions:

Observe employment restrictions for young people.

Skin resorption/Sensitization:

Causes allergic hypersensitivity reactions.

## 16. Other information

**Abbreviations and acronyms**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

**Further Information**

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

*(The data for the relevant ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*