



VITA VIONIC® TRY-IN RESIN

Technical Datasheet

Shade

- Tooth-colored resin (≈A2)

Content

- 1 kg bottle

Chemical composition

Polyfunctional Acrylates	46 - 65 %
Urethane bismethacrylate polymer	30 - 35 %
Methacrylate monomers	2 - 10 %
Initiators, stabilizers, pigments	2 - 6 %
MMA, TPO and phthalat free	0 %

Physical properties

Property	Test Standard	Value
Flexural strength	DIN EN ISO 178*	> 100 MPa
Flexural modulus	DIN EN ISO 178*	> 2200 MPa
Water sorption	DIN EN ISO 10477*	< 40 µg/mm ³
Water solubility	DIN EN ISO 10477*	< 7.5 µg/mm ³
Hardness		≥ 75 ShoreD
Viscosity**		850 - 1750 mPas

*Applies to cured objects | according to the norm at room temperature

**Applies to liquid resin

VITA Zahnfabrik H. Rauter GmbH & Co. KG

Spitalgasse 3
79713 Bad Säckingen
Germany

Phone: +49 7761 562-0
Hotline: +49 7761 562-222
E-Mail: info@vita-zahnfabrik.com

More information about VITA VIONIC TRY-IN RESIN can be found in the instructions for use and under: www.vita-zahnfabrik.com/VionicTryInResin
detax GmbH, Carl-Zeiss-Str. 4, 76275 Ettlingen
The product has been certified and bears the CE mark

CE0483 **MD**