

**Residual risks of glass ceramic CAD/CAM materials (such as VITA SUPRINITY PC)**

*Possible complications and residual risks of glass ceramic CAD/CAM materials in general, and not specific to VITA materials, are:*

- Peeling or loosening of the crown, veneer, onlay or inlay resulting in possible choking
- Plaque deposits
- Post-operative hypersensitivity
- Irritation of the mucosa
- Gingivitis
- Caries in marginal gaps
- Intolerance and/or potential allergic reaction
- Breakage

*These risks must be communicated to the patient by trained personnel (e.g., dental professionals).*

003/09.2025