

## Which vPad fits your requirements?

I want to ...	VITA vPad® LITE	VITA vPad® PRO
... create my own profile and individual firing programs.	✓ (VITA + 1 User)	✓ (VITA + 11 User)
... set up multiple digital workstations/ desktops.	✓ (up to 10 desktops)	✓ (unlimited flexibility)
... use SD slot for photo transfer.	✓ (1 x SD-Slot)	✓ (1 x SD-Slot)
... use USB interfaces for data transfer.	✓ (2 x USB)	✓ (2 x USB)
... transfer profiles, firing programs, and device data to other VITA vPads or PCs.	✓	✓
... send logged data online to VITA Service.	✓ (with Wi-Fi connection)	✓ (with Wi-Fi connection)
... control of up to 4 furnaces simultaneously.	✗	✓
... use a digital assistant function with step-by-step guidance.	✗	✓
... receive notifications on my smartphone with the Furnace app.	✗	✓
... set a timer for program start.	✗	✓
... use the integrated web browser.	✗	✓
... establish a Wi-Fi connection between furnace and vPad.	✗	✓