

Produktinformation

VITA Machinable Materials

Von Sirona empfohlene Mengen Dentatec Schleifzusatz in ml pro Tankfüllung:

Werkstoff	inLab	inLab MC XL	CEREC	CEREC MC XL
VITA CAD-Temp	50	75	50	75
VITA CAD-Waxx	10	75	10	75
VITABLOCS Mark II, TriLuxe, TriLuxe forte, RealLife	50	75	50	75
VITA In-Ceram YZ, AL, ZIRCONIA, ALUMINA, SPINELL	50	75	50	75



VITA