

VITA LUMIN® VACUUM

VITA Teeth



VITA shade taking

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 08.14

VITA shade, VITA made.

VITA


Formenkarte · Keramikzähne

Mould Chart · Ceramic denture teeth


Carte de formes · Dents en céramique

Carta de formas · Dientes de cerámica


Tabella delle forme · Denti di ceramica



schaufelförmig
ovoid
en forme de pelle
en forma de pala
a forma di pala





dreieckig
tapering
triangulaire
triangular
triangolare








Lower anteriors




rechteckig
rectangular
rectangulaire
rectangular
rettangolare





quadratisch
square
carré
cuadrado
quadrato








CUSPIFORM



SYNOFORM



VITA LUMIN® VACUUM Combination Chart

LUMIN® VACUUM Kombinationstabelle
LUMIN® VACUUM Combination Chart
Tableau de Combinaison LUMIN® VACUUM
Tabla de Combinación LUMIN® VACUUM
Tabella di Combinazione LUMIN® VACUUM

VITA LUMIN VACUUM Anteriores und Posteriores sind in den VITA classical Farben A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3 erhältlich.

Für weitere Farben fallen Lieferzeiten an.

VITA LUMIN VACUUM Anteriores and Posteriores are available in the VITA classical shades A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3.

Further colours are subject to delivery times.

VITA LUMIN VACUUM Anteriores et Posteriores sont disponibles dans les teintes VITA classical A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3.

Pour toutes les couleurs supplémentaires s'applique un délai de livraison.

I denti VITA LUMIN VACUUM Anteriores e Posteriores sono disponibili in colorazione VITA classical A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3.


Colori supplementari hanno un tempo di consegna.

La gama VITA LUMIN VACUUM Anteriores y Posteriores está disponible en los colores VITA classical A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3.

Otros colores tendrán plazo de entrega.

No.	 mm	No.  mm	No.  mm	
180	44.7 50.7	639 33.6	41C 29.5 31.8	11S 29.3 31.2
182	43.0 48.6	539 32.3	40C 28.3 29.5	10S 28.2 29.7
185	46.1 52.2	730 35.6	44C 32.3 34.3	13S 30.3 32.7
187	44.6 50.2	760 33.6	41C 29.5 31.8	11S 29.3 31.2
190	42.6 48.1	539 32.3	40C 28.3 29.5	10S 28.2 29.7
513	40.1 45.8	539 32.3	40C 28.3 29.5	10S 28.2 29.7
565	45.3 51.4	730 35.6	44C 32.3 34.3	13S 30.3 32.7
612	43.9 49.9	539 32.3	41C 29.5 31.8	11S 29.3 31.2
622	44.1 50.1	639 33.6	41C 29.5 31.8	11S 29.3 31.2
627	42.3 47.8	760 33.6	40C 28.3 29.5	10S 28.2 29.7
632	42.5 48.3	730 35.6	41C 29.5 31.8	11S 29.3 31.2
722	43.7 49.8	639 33.6	41C 29.5 31.8	11S 29.3 31.2
732	44.8 50.5	760 33.6	41C 29.5 31.8	11S 29.3 31.2
734	43.3 49.4	760 33.6	41C 29.5 31.8	11S 29.3 31.2
737	45.7 51.6	730 35.6	44C 32.3 34.3	13S 30.3 32.7
742	46.2 52.2	740 35.1	44C 32.3 34.3	13S 30.3 32.7
743	46.7 52.9	740 35.1	44C 32.3 34.3	13S 30.3 32.7
773	47.1 53.2	639 33.6	44C 32.3 34.3	13S 30.3 32.7
843	47.6 53.6	639 33.6	44C 32.3 34.3	13S 30.3 32.7
854	48.0 54.6	55 36.1	44C 32.3 34.3	13S 30.3 32.7

With the unique VITA SYSTEM 3D-MASTER, all natural tooth shades can be systematically determined and perfectly reproduced.



Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of application. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product. Furthermore, our liability for the accuracy of this information is independent of the legal basis and, in as far as legally permissible, shall always be limited to the value as invoiced of the goods supplied, excluding value-added tax. In particular, as far as legally permissible, we do not assume any liability for loss of earnings, indirect damages, ensuing damages or for third-party claims against the purchaser. Claims for damages based on fault liability (culpa in contrahendo, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence. The VITA Modulbox is not necessarily a component of the product.

Date of issue of this information: 08.14

After the publication of these working instructions any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA LUMIN VACUUM Teeth fulfil the requirements of ISO 4824 : 1993.

CE 0124

VITA

VITA Zahnfabrik H. Rauter GmbH & Co.KG
Postfach 1338 · D-79704 Bad Säckingen · Germany
Tel. +49(0)7761-562-0 · Fax +49(0)7761-562-299
Hotline: Tel. +49(0)7761-562-222 · Fax +49(0)7761-562-446
www.vita-zahnfabrik.com · info@vita-zahnfabrik.com
[facebook.com/vita.zahnfabrik](https://www.facebook.com/vita.zahnfabrik)

