
VITA LINGOFORM®

Product Information



Date of issue 03.09





Available in the VITA SYSTEM 3D-MASTER and VITAPAN classical A1-D4 shades



Productivity:

 Thanks to the multifunctional occlusal surface design, VITA LINGOFORM posteriors are quick and simple to set up - enabling high productivity and quality that satisfies the patient's needs.




Function

• The functionally designed occlusal areas facilitate the ideal static alignment of the prosthesis in situ, help reduce horizontal shear forces to the prosthesis bed and offer a high degree of wearing comfort.

Setup concepts

• The VITA LINGOFORM functional tooth is suitable for all customary setup concepts including lingualized occlusion: One for all.



Material:

• VITA MRP (microfiller reinforced polyacrylic). A material which is homogeneous throughout, and with outstanding material properties combining high mechanical strength with excellent biocompatibility.







Esthetics:

• The multilayer buildup and the nature identical surface morphology ensure lifelike, harmonious light reflection behavior

Indications:

- Complete denture prosthetics
- Implant prosthetics
- Telescopic constructions
- Conus constructions
- Attachments
- Model casting prostheses
- Interim restorations

VITA LINGOFORM® 4 molds in all shades - ideal for combined restorations



For further information see www.vita-zahnfabrik.com

VITA LINGOFORM acrylic teeth are available in the VITA SYSTEM 3D-MASTER and VITAPAN classical A1 - D4 shades. Shade compatibility with all VITA 3D-MASTER and VITAPAN classical materials is guaranteed.

With the unique VITA SYSTEM 3D-MASTER all natural tooth shades can be systematically determined and completely reproduced.



Please note: Our products should be used according to the working instructions. We cannot be held liable for damages resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers which are not compatible or not authorized for use with our product. Furthermore, our liability for the correctness of this information is independent of the legal ground and, in as far as legally permissible, is limited to the invoiced value of the goods supplied excluding turnover tax. In particular, as far as legally permissible, we do not assume any liability for profit loss, for indirect damages, for consequential damages or for claims of third parties against the purchaser. Claims for damages based on fault liability (fault in making the contract, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence. Date of issue of these instructions for use: 03.09

VITA Zahnfabrik has been certified according to the Guideline for Medical Devices and the following products bear the CE mark $\zeta \in 0124$:

VITA LINGOFORM®

US 5498157 A · AU 659964 B2 · EP 0591958 B1

