


VITAPAN®

Acrylic Teeth




VITA

The advantages of the classic VITAPAN® Anteriors

We paid close attention in geometry class!


VITAPAN acrylic teeth get an A+ as far as the variety of shapes is concerned. More than 30 different upper anterior sets – divided into four mould groups – guarantee ergonomic, simple and natural restorations for each case. Suitable lower anteriors are available for each indication and all requirements. Moreover, the clearly-structured presentation of the moulds facilitates communication between dentist and dental technician.



Ovoid Tapering Rectangular Square

An absolute "hard case".

A harsh law of nature says: "Only the strongest will survive." You will be on the safe side with VITAPAN acrylic teeth. Thanks to VITA's special repressing technique and superior material quality, VITAPAN teeth have held the record in abrasion resistance for more than 20 years – hardness offered to you and hard to exceed for competitors.



A comparative in-vitro study,
Quintessenz Zahntech 29,4,510-521
(2003)

Utmost brilliance for your teeth!

Outstanding physical properties and excellent aesthetics are combined in VITAPAN teeth. The intricate manual layering technique as well as the individual characterization add a natural appearance to the teeth.

Shade stability even after grinding, minimum water absorption and extremely low affinity to plaque ensure perfect results for demanding patients.



A unique formation

VITAPAN® Posteriors

CUSPIFORM



Excellence in design for mastication. The right solution for each indication.

VITAPAN teeth fulfill all requirements of modern, demanding dentures by offering two different types of occlusal concepts for a high degree and a low degree of intercuspation. The fully anatomical occlusal surfaces have a cuspal inclination of 23° to 28° and provide an optimum degree of intercuspation in the centric position.

SYNOFORM



VITAPAN SYNOFORM teeth with abraded, narrowly designed occlusal surfaces are suitable for a low degree of intercuspation, especially for dentures in geriatric dentistry. They ensure perfect integration in abraded residual dentition and allow easy determination of the occlusion. Any type of occlusion can be set up without the need of previous grinding. Both occlusal surface designs are suitable to prepare natural restorations safely and quickly.

Versatile teeth for all demands

VITAPAN® Acrylic teeth

Adding customization to prefabricated teeth

You don't have to book an all-inclusive vacation to have everything included. VITA provides their customers with a series of accessories suitable to customize VITAPAN acrylic teeth at any time. These accessories include the microparticle composite VITAVM LC. Naturally, VITAPAN anteriors and VITAPAN CUSPIFORM are also available in the well proven VITA SYSTEM 3D-MASTER shades.




VITAPAN® – better than Nasdaq, DAX and NEMAX

Superior profitability with maximum patient satisfaction – no, this is not an unrealistic promise of some investment bank but daily practice in thousands of dental laboratories for more than 20 years. Their natural morphology, superior material quality, simple processing and good price-performance ratio render VITAPAN acrylic teeth a true “asset” for everybody.

VITAPAN® tooth assortment. Closes any type of gap in dentition.

Anteriors: 31 upper anterior sets
13 lower posterior sets

Posteriors: 5 upper and lower sets of SYNOFORM posterior teeth
6 upper and lower sets of CUSPIFORM posterior teeth




„I'll certainly stand up to any challenge.“



„You seem to be quite tough.“



The unique VITA SYSTEM 3D-MASTER allows the systematic determination and complete reproduction of all natural tooth shades.



Please note: Our products have been developed for the use in the dental sector and have to be processed in accordance with the directions for use. Illustrations and written information do not include any undertaking as to qualities and are not binding.

VITA Zahnfabrik is certified according to the Medical Device Directive and the following products mentioned in this brochure bear the mark  0124 and comply with the requirements of the standard ISO 3336:1993.

VITAPAN® acrylic teeth

VITAVM®LC

„When shape matters, just take it!“



VITA

VITA Zahnfabrik H. Rauter GmbH & Co.KG
Postfach 1338 · D-79704 Bad Säckingen · Germany
Tel. + 49/7761/562-222 · Fax + 49/7761/562-446
www.vita-zahnfabrik.com · info@vita-zahnfabrik.com