

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 11.04.2013

Version number 2

Revision: 11.04.2013

Hazardous according to criteria of Australian Safety and Compensation Council

1 Identification of the substance/mixture and of the company/undertaking

- **Product identifier**
- **Trade name:** VITA DUO CEMENT Catalyst
- **Article number:** FCDCK
- **Relevant identified uses of the substance or mixture and uses advised against** Not determined
- **Application of the substance / product** Auxiliary for manufacture of dental prothesis
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
VITA Zahnfabrik
H. Rauter GmbH & Co. KG
- Postfach 1338
D 79704 Bad Säckingen
e-mail: info@vita-zahnfabrik.com
- **Emergency telephone number:** Tel.: ++49-(0)761-19240

Tel.: 07761/562-0
Fax: 07761/562299

2 Hazards identification

- **Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**



GHS07

Skin Sens. 1 H317 May cause an allergic skin reaction.

- **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**
not applicable



Xi; Irritant

R36/37/38: Irritating to eyes, respiratory system and skin.



Xi; Sensitising

R43: May cause sensitisation by skin contact.

- **Classification system:**
The classification is according to the latest editions of the EU-lists, and extended by company and literature data.

- **Label elements**

- **Labelling according to EU guidelines:**

The product has been classified and marked in accordance with EU Directives / Ordinance on Hazardous Materials

Observe the general safety regulations when handling chemicals

- **Code letter and hazard designation of product:**



Xi Irritant

- **Hazard-determining components of labelling:**
dibenzoyl peroxide

(Contd. on page 2)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 11.04.2013

Version number 2

Revision: 11.04.2013

Trade name: VITA DUO CEMENT Catalyst













(Contd. of page 1)

- **Risk phrases:**
 - 36/37/38 Irritating to eyes, respiratory system and skin.
 - 43 May cause sensitisation by skin contact.
- **Information concerning to particular hazards to man and environment**
- **Safety phrases:**
 - 23 Do not breathe gas/fumes/vapour/spray.
 - 24 Avoid contact with skin.
 - 29 Do not empty into drains.
 - 37 Wear suitable gloves.
 - 60 This material and its container must be disposed of as hazardous waste.
- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of substances listed below with nonhazardous additions.

- **Dangerous components:**

CAS: 109-16-0	triethylene glycol dimethacrylate  Xi R36/37/38;  Xi R43  Skin Irrit. 2, H315;  Eye Irrit. 2, H319;  Skin Sens. 1, H317; STOT SE 3, H335	3-<10%
CAS: 94-36-0 EINECS: 202-327-6 Reg.nr.: 01-02119511472-50	dibenzoyl peroxide  Xi R36;  Xi R43;  E R2  Org. Perox. B, H241;  Aquatic Acute 1, H400;  Eye Irrit. 2, H319;  Skin Sens. 1, H317; Expl. 1.5, H205	< 2.5%

- **Additional information** For the wording of the listed risk phrases refer to section 16.

4 First aid measures

- **Description of first aid measures**
- **General information** No special measures required.
- **After inhalation** Supply fresh air.
- **After skin contact** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact** Rinse opened eye for several minutes under running water.
- **After swallowing** If symptoms persist consult doctor.
- **Information for doctor**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Firefighting measures

- **Extinguishing media**
- **Suitable extinguishing agents** Use fire extinguishing methods suitable to surrounding conditions.
- **Special hazards arising from the substance or mixture** No further relevant information available.
- **Advice for firefighters**
- **Protective equipment:** No special measures required.

AUS

(Contd. on page 3)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 11.04.2013

Version number 2

Revision: 11.04.2013

Trade name: VITA DUO CEMENT Catalyst

(Contd. of page 2)

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Not required.
- **Environmental precautions:** Do not allow product to reach sewage system or any water course.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections** No dangerous substances are released.

7 Handling and storage

- **Handling**
- **Precautions for safe handling** No special measures required.
- **Information about fire - and explosion protection:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage**
- **Requirements to be met by storerooms and receptacles:** Store in cool, dry place in tightly closed receptacles.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** None.
- **Storage class** 10
- **Specific end use(s)** No further relevant information available.

8 Exposure controls/personal protection

- **Additional information about design of technical facilities:** No further data; see item 7.
- **Control parameters**
- **Ingredients with limit values that require monitoring at the workplace:**
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists valid during the making were used as basis.
- **Exposure controls**
- **Personal protective equipment**
- **General protective and hygienic measures**
The usual precautionary measures are to be adhered to when handling chemicals.
- **Respiratory protection:** Not required.
- **Protection of hands:**
PVC gloves



Protective gloves.

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

- **Material of gloves**
Nitrile rubber, NBR

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

- **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 4)

AUS

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 11.04.2013

Version number 2

Revision: 11.04.2013

Trade name: VITA DUO CEMENT Catalyst

(Contd. of page 3)

· **Eye protection:**

Tightly sealed goggles.

9 Physical and chemical properties

· **Information on basic physical and chemical properties**· **General Information**· **Appearance:**

Form:	Pasty
Colour:	Whitish
Odour:	Characteristic

· **Change in condition**

Melting point/Melting range:	undetermined
Boiling point/Boiling range:	undetermined

Flash point:	Not applicable
--------------	----------------

Ignition temperature:	Not applicable
-----------------------	----------------

Decomposition temperature:	Not applicable
----------------------------	----------------

Self-igniting:	Product is not selfigniting.
----------------	------------------------------

Danger of explosion:	Product does not present an explosion hazard.
----------------------	---

· **Explosion limits:**

Oxidizing properties	Not determined
----------------------	----------------

Vapour pressure:	Not determined
------------------	----------------

Density at 20 °C:	2.1 g/cm ³
-------------------	-----------------------

· **Solubility in / Miscibility with**

Water:	Not miscible or difficult to mix
--------	----------------------------------

Partition coefficient (n-octanol/water):	Not determined
--	----------------

· **Viscosity:**

dynamic:	Not determined
----------	----------------

· **Solvent content:**

Organic solvents:	0.0 %
-------------------	-------

· **Other information**

No further relevant information available.

10 Stability and reactivity

· **Reactivity** No decomposition if used and stored according to specifications.· **Chemical stability** No decomposition if used and stored according to specifications.· **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.· **Possibility of hazardous reactions** No dangerous reactions known· **Conditions to avoid** No further relevant information available.· **Incompatible materials:** No further relevant information available.· **Hazardous decomposition products:** No dangerous decomposition products known

AUS

(Contd. on page 5)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 11.04.2013

Version number 2

Revision: 11.04.2013

Trade name: VITA DUO CEMENT Catalyst

(Contd. of page 4)

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **LD/LC50 values relevant for classification:** Not determined
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** Irritating effect.
- **Sensitization:** Sensitization possible through skin contact.

12 Ecological information

- **Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **Persistence and degradability** No further relevant information available.
- **Behaviour in environmental systems:**
- **Bioaccumulative potential** No further relevant information available.
- **Mobility in soil** No further relevant information available.
- **Additional ecological information:**
- **General notes:**
At present there are no ecotoxicological assessments.
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **Other adverse effects** No further relevant information available.

13 Disposal considerations

- **Waste treatment methods**
- **Recommendation** Disposal must be made according to official regulations.
- **Uncleaned packaging:**
- **Recommendation:** Disposal must be made according to official regulations.

14 Transport information

- | | |
|---------------------------------------|-----------------|
| · UN-Number | |
| · ADN, IMDG, IATA | Void |
| · UN proper shipping name | |
| · ADG, ADN, IMDG, IATA | Void |
| · Transport hazard class(es) | |
| · ADN/R Class: | Void |
| · Packing group | |
| · IMDG, IATA | Void |
| · Environmental hazards: | |
| · Marine pollutant: | No |
| · Special precautions for user | Not applicable. |

(Contd. on page 6)

AUS

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 11.04.2013

Version number 2

Revision: 11.04.2013

Trade name: VITA DUO CEMENT Catalyst

(Contd. of page 5)

- | | |
|--|-----------------|
| · Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code | Not applicable. |
| · UN "Model Regulation": | - |

15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture**

- **Australian Inventory of Chemical Substances**

All ingredients are listed.

- **Standard for the Uniform Scheduling of Drugs and Poisons**

94-36-0 dibenzoyl peroxide	S2, S4, S5
------------------------------	------------

- **National regulations**

- **Waterhazard class:** Water hazard class 1 (Self-assessment): slightly hazardous for water.
- **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**

H205 May mass explode in fire.

H241 Heating may cause a fire or explosion.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

H335 May cause respiratory irritation.

H400 Very toxic to aquatic life.

R2 Risk of explosion by shock, friction, fire or other sources of ignition.

R36 Irritating to eyes.

R36/37/38 Irritating to eyes, respiratory system and skin.

R43 May cause sensitisation by skin contact.

- **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent