

VITA SUPRINITY® – Polishing Set clinical/technical

For simple pre- and high-gloss polishing with VITA SUPRINITY.



VITA shade taking

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 06.13

VITA shade, VITA made.

VITA

Product description

The VITA SUPRINITY Polishing Sets were developed for reliable, efficient and material-specific surface treatment of zirconia reinforced lithium silicate ceramic (ZLS) restorations in dental practices and laboratories. The sets include various polishing instruments for pre- and high-gloss polishing.

These instruments are suitable for careful and gentle polishing of occlusal surfaces, cusps, fissures and restoration contact points. They ensure a brilliant shine on the finished restoration.

Product benefits

Excellent final results:

The instruments produce highly esthetic and plaque-resistant surfaces. Precise concentricity, matched grit sizes and the individual geometries of the instruments guarantee superior results.

Simple and safe handling:

The instruments provide removal capacity that can be easily controlled with minimal wear. Good handling and the ability to use them without polishing paste enables simple and fast processing. Safe use is also guaranteed since the instruments can be sterilized.

Gentle and careful processing:

These instruments, which were developed especially for VITA SUPRINITY, ensure gentle and careful processing. As a result, the risk of possible formation of microcracks is reduced.

Product data

VITA SUPRINITY Polishing Set clinical:

Includes a total of six instruments for the contra-angle. Includes three instruments each for pre- and high-gloss polishing.

VITA SUPRINITY Polishing Set technical:

Includes a total of eight instruments for the handpiece. Includes four instruments each for pre- and high-gloss polishing.

VITA SUPRINITY® assortments

Designation	Details	Content/piece	Prod. No.
VITA SUPRINITY Polishing Set clinical	2-stage polishing system with pre- and high-gloss polishers for the contra-angle	6	ELSPSETC
VITA SUPRINITY Polishing Set technical	2-stage polishing system with pre- and high-gloss polishers for the handpiece	8	ELSPSETT

VITA SUPRINITY® Polishing Set clinical – refill packs

Designation	Details	Content/piece	Prod. no.
Pre-polishing (pink)	Diamond-coated		
Polishing instrument, tip, S4m		6	ERWS4M6
Polishing instrument, lens, L10m		6	ERWL10M6
Polishing instrument, cup, C6m		6	ERWC6M6
High-gloss polishing (grey)	Diamond-coated		
Polishing instrument, tip, S4f		6	ERWS4F6
Polishing instrument, lens, L10f		6	ERWL10F6
Polishing instrument, cup, C6f		6	ERWC6F6

VITA SUPRINITY® Polishing Set technical – refill packs

Designation	Details	Content/piece	Prod. No.
Pre-polishing (pink)	Diamond-coated		
Polishing instrument, wheel, R15m		6	ERHR15M6
Polishing instrument, tip, S5m		6	ERHS5M6
Polishing instrument, lens, L16m		6	ERHL16M6
Polishing instrument, long pin, LPm		6	ERHLPm6
High-gloss polishing (grey)	Diamond-coated		
Polishing instrument, wheel, R15f		6	ERHR15F6
Polishing instrument, tip, S5f		6	ERHS5F6
Polishing instrument, lens, L16f		6	ERHL16F6
Polishing instrument, long pin, LPf		6	ERHLPf6

More information about VITA SUPRINITY is available at: www.vita-suprinity.de / www.vita-suprinity.com

VITA