


VITA MFT®

Multi Functional Teeth

Mould chart acrylic teeth



VITA shade determination

VITA shade communication

VITA shade reproduction


VITA shade control

Date of issue: 01.15


VITA shade, VITA made.

VITA

Available in VITA SYSTEM 3D-MASTER
and VITA classical A–D shades



Ovoid



Triangular

040



T41



044




T43



049



T46



T49



VITA MFT – The new dimension in aesthetics

The natural-looking anatomy of VITA MFT anteriors, with approximal contacts in the area of the incisal third being slightly placed towards the palate, enables a trouble-free and appearingly organic alignment of the teeth.

The anatomically correct labial curve leads to a physiologically perfect set-up of the teeth where discretely broadened tooth necks support an age-appropriate papilla design.




Lower anteriors



Rectangular

R42



L33



R45




L34



R47



L37



Square

L37LN



S47



L40



S50*



The mould **L37LN** (**Long Neck**) was specifically developed for simulating particular periodontal characteristics which, for example, occur in partial prosthetics.

PU29 – PL29



PU31 – PL31





PU33 – PL33



easy-centric® – the new, effortless way for setting up dentures

The occlusal surface design of VITA MFT is based on the innovative principle called easy-centric which supports a professional and economical set-up procedure. Moreover, it leads to increased productivity in a quality well adapted to the patient's needs.

No.		No.		
	mm	mm	No.	mm
O40	40,0	L33 32,6	PU29 28,8	PL29 29,6
O44	43,9	L34 34,0	PU29 28,8	PL29 29,6
O49	49,2	L40 39,5	PU33 33,0	PL33 33,6
T41	40,5	L33 32,6	PU29 28,8	PL29 29,6
T43	43,3	L34 34,0	PU31 30,6	PL31 31,1
T46	46,1	L37 36,8	PU31 30,6	PL31 31,1
T49	48,8	L40 39,5	PU33 33,0	PL33 33,6
R42	42,1	L34 34,0	PU29 28,8	PL29 29,6
R45	44,5	L37 36,8	PU31 30,6	PL31 31,1
R47	46,6	L37 36,8	PU33 33,0	PL33 33,6
S47	46,8	L40 39,5	PU33 33,0	PL33 33,6
S50	49,7	L40 39,5	PU33 33,0	PL33 33,6



VITA MFT – Anteriors and posteriors

The whole range of indications in the pleasantly small VITA MFT product line. A reduced and harmonious range of shades/moulds ensures for increased economic efficiency.

The multifunctional VITA MFT tooth is suited for:

- Full dentures
- Implant prosthetics
- Telescope and cone constructions
- Attachment Dentures
- Model cast prostheses
- Interim solutions

The shade range of VITA MFT covers


- 13 VITA SYSTEM 3D-MASTER shades

VITA SYSTEM 3D-MASTER												
1M1	2L1.5	2M1	2M2	3L1.5	3L2.5	3M2	3R2.5	4L1.5	4M2	5M1	0M1	0M3
~ B1	~ B2		~ A2	~ C2 ~ D4	~ B3		~ A3.5		~ A4	~ C4		

- 13 VITA classical A–D shades

VITA classical A–D												
A1	~A2	A3	~A3.5	~A4	~B1	~B2	~B3	~C2	C3	~C4	D3	~D4
	2M2		3R2.5	4M2	1M1	2L1.5	3L2.5	3L1.5		5M1		3L1.5

With the unique VITA SYSTEM 3D-MASTER, all natural tooth shades can be systematically determined and perfectly reproduced.



Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of application. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product. Furthermore, our liability for the accuracy of this information is independent of the legal basis and, in as far as legally permissible, shall always be limited to the value as invoiced of the goods supplied, excluding value-added tax. In particular, as far as legally permissible, we do not assume any liability for loss of earnings, indirect damages, ensuing damages or for third-party claims against the purchaser. Claims for damages based on fault liability (culpa in contrahendo, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence. The VITA Modulbox is not necessarily a component of the product.

Date of issue of this information: 01.15

After the publication of these information for use any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified in accordance to the Medical Device Directive and the following products bear the CE mark **CE** 0124 :

VITA MFT®

VITA

VITA Zahnfabrik H. Rauter GmbH & Co.KG
 Spitalgasse 3 · D-79713 Bad Säckingen · Germany
 Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299
 Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446
www.vita-zahnfabrik.com · info@vita-zahnfabrik.com
 facebook.com/vita.zahnfabrik