

Completely natural. Highly functional. Prostheses using VITAPAN PLUS® and VITA LINGOFORM®

The perfect combination for anterior esthetics and posterior function



Date of issue: 09.15




VITA shade, VITA made.

VITA

A natural appearance, superior wear comfort and functionality: these are the requirements that patients want their dentures to meet. With VITAPAN PLUS anterior teeth and VITA LINGOFORM posterior teeth, you can achieve highly esthetic results in the anterior area as well as outstanding posterior

functionality — particularly in the case of mucosa-borne and **implant-supported full dentures**. See for yourself just how impressive the material properties of both tooth ranges are, and what benefits they offer when it comes to processing: Their perfect esthetics and unique VITA^lity will inspire you.



VITAPAN PLUS® – perfect anterior esthetics

Perfect for implant esthetics

- ▶ Thanks to the optimally-sized tooth-neck width

Clever multi-layer concept

- ▶ Optimally adapted opalescence and luminescence for a natural play of colors rich in nuances
- ▶ Lifelike translucency that lends the teeth a youthful appearance

Reflection of light in 3 dimensions

- ▶ Thanks to an innovative and completely natural-looking surface morphology

High abrasion resistance and material strength

- ▶ Thanks to densely cross-linked MRP composite material (**M**icrofiller **R**einforced **P**olyacrylic material)



VITA LINGOFORM® – the posterior tooth for reliable centric relation

- ▶ Combine the VITAPAN PLUS anterior teeth with VITA LINGOFORM posterior teeth for **fast and reliable setup**.
- ▶ Particularly suited for simple and precise setup in lingualized occlusion

- ▶ VITA LINGOFORM posterior teeth also ensure reliable occlusion and articulation for all setup concepts.

Ideal combination: VITAPAN PLUS® anterior & VITA LINGOFORM® posterior teeth

Indication	VITAPAN PLUS® Anterior tooth	VITA LINGOFORM® Posterior tooth
Full dentures	✓✓✓	✓✓✓
Implant prostheses	✓✓✓	✓✓✓
Temporary dentures	✓✓	✓✓✓
Model cast prostheses	✓	✓✓✓
Partial prostheses	✓	✓✓✓
Attachment Prostheses	✓	✓✓✓
Telescopic/ double crown ret. dentures	✓	✓✓✓

VITAPAN PLUS anterior teeth and VITA LINGOFORM posterior teeth are ideally suited for full dentures as well as implant-supported dentures.

- ✓✓✓ recommended
- ✓✓ very suitable
- ✓ suitable



Urban Christen, dentist / prosthodontist,
Switzerland

Perfect combination in practical applications

The prosthetic specialist Urban Christen has been using VITAPAN PLUS anterior teeth and VITA LINGOFORM posterior teeth for years in his practice. "VITAPAN PLUS offers a highly-esthetic anterior tooth that my patients also like right away.

In my opinion, VITA LINGOFORM offers the ultimate in posterior functionality: with its anatomically pre-abraded occlusal surfaces, I can achieve perfect occlusal function instantly, regardless of the intended occlusal concept."



To sum it up: VITAPAN PLUS and VITA LINGOFORM offer me the ideal combination of anterior esthetics and posterior functionality.



VITA