

VITA Easyshade® V -

Instructions for carrying out software updates with the VITA Easyshade V

Dear customer,

With this document, we would like to inform you about the new software update for your VITA Easyshade V and ask you to install the update on your instrument according to the following instructions.

If you have any questions concerning this update, please consult your regional service partner for VITA Equipment. VITA Equipment is also available at:

- Tel.: +49 (0)77 61 / 562 – 101, -105, -106 or -114
- Fax: +49 (0)77 61 / 562 – 102
- Email: instruments-service@vita-zahnfabrik.com

1. Download and installation prerequisites

- a) Internet access
- b) Bluetooth-capable Windows PC or Bluetooth dongle from the VITA Easyshade V assortment

2. General Procedure




Expected duration of the update is approx. 8 mins

1. Download software update from the Internet
2. Activate Bluetooth on the VITA Easyshade V
3. Select software update, device browsing
4. Select the desired VITA Easyshade V by clicking on it
5. Update window opens
6. Begin the update by pressing "Start"
7. Do not carry out any measuring operations etc. with the VITA Easyshade V during the update.

VITA Easyshade® V -




Instructions for carrying out software updates with the VITA Easyshade V

1. Download new software update from the Internet and save it to the computer

Step	Input by user	Function
1	Enter www.vita-zahnfabrik.com in your web browser.	
2	Select "Service" + "VITA Update Messenger"	
3	Select "To newsletters". Then select "VITA Update Messenger 01/2016" and confirm with "To newsletter".	





VITA Easyshade® V -

Instructions for carrying out software updates with the VITA Easyshade V

Step	Input by user	Function
4	Select "Download installation instructions" or "Download".	
5	Select "Save".	
6	Save the file in the directory of your choice and start VITA_EasyshadeV_Release_V505h_Pkg.exe by double-clicking. The update program will be copied into a directory of your choice.	
7	Activate Bluetooth on the VITA Easyshade V.	
8	Double-click the file on the selected drive.	

VITA Easyshade® V -

Instructions for carrying out software updates with the VITA Easyshade V

Step	Input by user	Function
9	Scan for devices.	 
10	Select the desired VITA Easyshade V by clicking on it. If only one VITA Easyshade V is available, selection is automatic.	
11	The update window opens; click on "Start" to begin. (Expected duration approx. 8 mins) During the update, the VITA Easyshade V's display remains dark (and cannot be operated).	

VITA Easyshade® V -

Instructions for carrying out software updates with the VITA Easyshade V

Step	Input by user	Function
12	Close the program once the update has finished. Restart the VITA Easyshade V.	