


VITA ENAMIC® Polishing Set clinical/technical

For simple pre- and high-gloss polishing of VITA ENAMIC hybrid ceramic



Date of issue: 09.16

VITA shade, VITA made. **VITA**

Product description

The VITA ENAMIC polishing sets were developed for reliable, efficient and material-specific surface treatment of hybrid ceramic restorations in dental practices and laboratories. The sets include various polishing instruments for pre- and

high-gloss polishing. These instruments are suitable for careful and gentle polishing of occlusal surfaces, cusps, fissures and contact points of the restoration and produce surfaces with exceptional gloss.

Product advantages

Excellent final results

High Gloss and plaque-resistant surfaces can be produced with these instruments. Precise concentricity, matched grit sizes and the individual geometries of the instruments guarantee results with superior precision.

Simple and safe handling

The instruments provide superior material removal performance. Good handling and the use without polishing paste enable simple and fast processing.

Gentle and careful processing

These instruments were developed specifically for the VITA ENAMIC to ensure gentle and careful polishing, reducing the risk of possible microcrack formation.

Product data

VITA ENAMIC Polishing Set clinical: includes a total of 6 instruments for the contra-angle; 3 instruments for pre-polishing and 3 instruments for high-gloss polishing.

VITA ENAMIC Polishing Set technical: includes a total of 8 instruments for the handpiece; 4 instruments for pre-polishing and 4 instruments for high-gloss polishing.

VITA ENAMIC® sets

Designation	Description	Content / pieces	Prod. No.
VITA ENAMIC Polishing Set clinical	2-stage polishing system with pre- and high-gloss polishers for the contra-angle	8	EENPSETCV1
VITA ENAMIC Polishing Set technical	2-stage polishing system with pre- and high-gloss polishers for the handpiece	8	EENPSETT

VITA ENAMIC® Polishing Set clinical – refill packs

Designation	Description	Content / pieces	Prod. No.
Pre-polishing (pink)	SiC instruments (silicon carbide)		
Polishing instrument, tip, VI-EB14m		3	ERWEB14M3
Polishing instrument, lens, VI-EL10m		6	ERWEL10M6
Polishing instrument, cup, VI-EC7m		6	ERWEC7M6
High-gloss polishing (grey)	Diamond instruments		
Polishing instrument, tip, VI-EB14f		3	ERWEB14F3
Polishing instrument, lens, VI-EL10f		6	ERWEL10F6
Polishing instrument, cup, VI-EC6f		6	ERWEC6F6

VITA ENAMIC® Polishing Set technical – refill packs

Designation	Description	Content / pieces	Prod. No.
Pre-polishing (pink)	SiC instruments (silicon carbide)		
Polishing instrument, tip, VI-ES5m		6	ERHES5M6
Polishing instrument, torpedo, VI-ET6m		6	ERHET6M6
Polishing instrument, wheel, VI-ER15m		6	ERHER15M6
Polishing instrument, lens, VI-EL16m		6	ERHEL16M6
High-gloss polishing (grey)	Diamond instruments		
Polishing instrument, tip, VI-ES5f		6	ERHES5F6
Polishing instrument, torpedo, VI-ET6f		6	ERHET6F6
Polishing instrument, wheel, VI-ER15f		6	ERHER15F6
Polishing instrument, lens, VI-EL16f		6	ERHEL16F6

More information about VITA ENAMIC is available at www.vita-enamic.de / www.vita-enamic.com

VITA