

> VITA ENAMIC® multiColor

> CAD/CAM blank made of hybrid ceramics



IDS
Highlights
2017

> FACTS & FEATURES

VITA ENAMIC® multiColor

- CAD/CAM blank with integrated shade gradient from cervical to incisal area

You benefit from:

- Natural play of light thanks to the blank's integrated shade gradient
- Fast processing, since the restoration can be seated as soon as polishing is complete
- Enormous load capacity, since masticatory forces are absorbed

VITA ENAMIC® multiColor is:

- The hybrid ceramic blank for reproduction of the natural shade gradient at the push of a button

Available shades and geometries:

- 1M1-HT, 1M2-HT, 2M2-HT, 3M2-HT, 4M2-HT
- **EMC-14:** 12 x 14 x 18 mm



> SYSTEM COMPONENTS

