

# VITA IMPLANT SOLUTIONS –

Perfect solutions for implant-supported restorations



VITA shade determination

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 05.17

VITA – perfect match.

**VITA**

## Product description

VITA IMPLANT SOLUTIONS (IS) are blanks made of hybrid ceramic, glass ceramic or composite, which feature an integrated interface to an adhesive/titanium base (e.g. TiBase)

and are used for CAD/CAM fabrication of implant-supported restorations. They are available in three versions: VITA CAD-Temp IS, VITA ENAMIC IS and VITA SUPRINITY IS.

# VITA IMPLANT SOLUTIONS

## Indication

For tooth-colored mesostructures (MS) and/or abutment crowns (AC)

- temporary restorations/soft tissue management:  
VITA CAD-Temp IS
- definitive superstructures:  
VITA ENAMIC IS/VITA SUPRINITY IS

## Benefits

### VITA ENAMIC IS:

- absorbing material properties thanks to dentin-like elasticity
- highly accurate results in thinning margin areas
- faster processing since superstructure can be seated immediately after CAM fabrication/polishing

### VITA SUPRINITY IS:

- high-strength, zirconia-reinforced glass ceramic for good reliability

### VITA CAD-Temp IS:

- efficient fabrication of temporaries
- simple design of the emergence profile

## Geometries

VITA IMPLANT SOLUTIONS blanks are available

- in IS-14 geometry: 18 x 14 x 12 (for mesostructures)
- in IS-16 geometry: 18 x 16 x 18 (for abutment crowns)
- with integrated S- or L-interface for each geometry

## Shades

VITA IMPLANT SOLUTIONS blanks are available

- as VITA CAD-Temp IS-16 in 3 VITA SYSTEM 3D-MASTER shades
- as VITA ENAMIC IS-14 in 5 VITA SYSTEM 3D-MASTER shades, in Translucent (T)
- as VITA ENAMIC IS-16 in 5 VITA SYSTEM 3D-MASTER shades, in High Translucent (HT)
- as VITA SUPRINITY IS-14 in 4 VITA classical/  
6 VITA SYSTEM 3D-MASTER shades, in Translucent (T)
- as VITA SUPRINITY IS-16 in 4 VITA classical/  
6 VITA SYSTEM 3D-MASTER shades, in High Translucent (HT)

## VITA IMPLANT SOLUTIONS blanks

### VITA ENAMIC® IS Translucent (T)

Shade	Designation	Geometry	Content	Prod. No.
1M1-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4EM3069605
1M1-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4EM3069725
1M2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4EM3079605
1M2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4EM3079725
2M2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4EM3129605
2M2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4EM3129725
3M2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4EM3209605
3M2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4EM3209725
4M2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4EM3289605
4M2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4EM3289725

### VITA ENAMIC® IS High Translucent (HT)

Shade	Designation	Geometry	Content	Prod. No.
1M1-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4EM4069365
1M1-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4EM4069485
1M2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4EM4079365
1M2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4EM4079485
2M2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4EM4129365
2M2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4EM4129485
3M2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4EM4209365
3M2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4EM4209485
4M2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4EM4289365
4M2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4EM4289485

**VITA SUPRINITY® IS Translucent (T)\***

Shade	Designation	Geometry	Content	Prod. No.
A1-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03369605
A1-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03369725
A2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03379605
A2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03379725
A3-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03389605
A3-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03389725
A3.5-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03399605
A3.5-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03399725
0M1-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03029605
0M1-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03029725
1M1-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03069605
1M1-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03069725
1M2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03079605
1M2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03079725
2M2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03129605
2M2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03129725
3M2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03209605
3M2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03209725
4M2-T	IS-14S	18 x 14 x 12 mm	5 pieces	EC4S03289605
4M2-T	IS-14L	18 x 14 x 12 mm	5 pieces	EC4S03289725

**VITA SUPRINITY® IS High Translucent (HT)\***

Shade	Designation	Geometry	Content	Prod. No.
A1-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04369365
A1-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04369485
A2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04379365
A2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04379485
A3-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04389365
A3-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04389485
A3.5-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04399365
A3.5-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04399485
0M1-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04029365
0M1-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04029485
1M1-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04069365
1M1-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04069485
1M2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04079365
1M2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04079485
2M2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04129365
2M2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04129485
3M2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04209365
3M2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04209485
4M2-HT	IS-16S	18 x 16 x 18 mm	5 pieces	EC4S04289365
4M2-HT	IS-16L	18 x 16 x 18 mm	5 pieces	EC4S04289485

\*) Note: Product information about VITA SUPRINITY IS is preliminary, the blanks with integrated interface will probably be available in 2018.  
VITA SUPRINITY IS will not be available in the U.S.A.

### VITA CAD-Temp® IS

Shade	Designation	Geometry	Content	Prod. No.
1M2T	IS-16S	18 x 16 x 18 mm	5 pieces	EC4CB0079365
1M2T	IS-16L	18 x 16 x 18 mm	5 pieces	EC4CB0079485
2M2T	IS-16S	18 x 16 x 18 mm	5 pieces	EC4CB0129365
2M2T	IS-16L	18 x 16 x 18 mm	5 pieces	EC4CB0129485
3M2T	IS-16S	18 x 16 x 18 mm	5 pieces	EC4CB0209365
3M2T	IS-16L	18 x 16 x 18 mm	5 pieces	EC4CB0209485

### VITA IMPLANT SOLUTIONS – SYSTEM SOLUTIONS

VITA offers VITA IMPLANT SOLUTIONS blanks with a specific holder system for the CAD/CAM system:

- CEREC/inLab (Sirona Dental Systems GmbH)

**VITA**