

VITA IMPLANT SOLUTIONS

VITA CAD-Temp® IS

VITA ENAMIC® IS

Kurzinformation

Stand: 05.17

Indikation

Für zahnfarbene Meso-Strukturen (MS) und/oder Abutment-Kronen (AC) mit integrierter Schnittstelle zu einer Klebe-/Titanbasis (z.B. Sirona TiBase).

- Provisorische Versorgungen/Weichgewebsmanagement:
VITA CAD-Temp IS
- Definitive Suprakonstruktionen:
VITA ENAMIC IS



VITA – perfect match.

VITA



Softwarevoraussetzungen

CEREC: Software ≥ 4.4

inLab: Software ≥ 15.0



Bitte Gebrauchsanweisung beachten

Mehr Informationen

- VITA IMPLANT SOLUTIONS Verarbeitungsanleitung Nr. 10150
- VITA IMPLANT SOLUTIONS Produktblatt Nr. 10154



Unter www.vita-zahnfabrik.com/VITA_ENAMIC_IS#titan_compatibility bzw. www.vita-zahnfabrik.com/VITA_CADTemp_IS#titan_compatibility finden Sie die aktuellen länderspezifischen Listen der für die VITA IMPLANT SOLUTIONS Rohlinge zugelassenen Titanbasen.



Die VITA Zahnfabrik ist nach der Medizinprodukterichtlinie zertifiziert und folgende Produkte tragen die Kennzeichnung  0124: **VITA ENAMIC® IS, VITA CAD-Temp® IS**

CEREC® und inLab® sind eingetragene Marken der Firma Sirona Dental Systems GmbH, D-Bensheim

VITA Zahnfabrik H. Rauter GmbH & Co.KG

Spitalgasse 3 | D-79704 Bad Säckingen | Germany

Tel. +49 (0) 7761/562-0 | Fax +49 (0) 7761/562-299

Hotline: Tel. +49 (0) 7761/562-222 | Fax +49 (0) 7761/562-446

www.vita-zahnfabrik.com | info@vita-zahnfabrik.com

 facebook.com/vita.zahnfabrik



VITA – perfect match.

VITA

VITA IMPLANT SOLUTIONS

VITA CAD-Temp® IS

VITA ENAMIC® IS

Summary information

Date of issue: 05.17

Indication

For tooth-colored mesostructures (MS) and/or abutment crowns (AC) with an integrated interface to an adhesive/ titanium base (such as Sirona TiBase).

- temporary restorations/soft tissue management:
VITA CAD-Temp IS
- definitive superstructures:
VITA ENAMIC IS



VITA – perfect match.

VITA



Software requirements

CEREC: software ≥ 4.4

inLab: software ≥ 15.0



See instructions for use

Further information

- VITA IMPLANT SOLUTIONS working instructions no. 10150
- VITA IMPLANT SOLUTIONS product information no. 10154



You will find the current country-specific lists of the titanium bases approved for use with the VITA IMPLANT SOLUTIONS blanks at www.vita-zahnfabrik.com/VITA_ENAMIC_IS#titan_compatibility or www.vita-zahnfabrik.com/VITA_CADTemp_IS#titan_compatibility



VITA Zahnfabrik has been certified in accordance with the Medical Device Directive and the following products bear the CE mark **CE 0124: VITA ENAMIC® IS, VITA CAD-Temp® IS**

Rx only

CEREC® and inLab® are registered trademarks of Sirona Dental Systems GmbH, Bensheim, Germany

VITA Zahnfabrik H. Rauter GmbH & Co.KG

Spitalgasse 3 | D-79704 Bad Säckingen | Germany

Tel. +49 (0) 7761/562-0 | Fax +49 (0) 7761/562-299

Technical Hotline: Tel. +49 (0) 7761/562-222 | Fax +49 (0) 7761/562-446

www.vita-zahnfabrik.com | info@vita-zahnfabrik.com

[facebook.com/vita.zahnfabrik](https://www.facebook.com/vita.zahnfabrik)

VITA – perfect match.

VITA