

Safety Data Sheet

according to WHMIS

VITA SPRAY-ON INDICATOR LIQUID

Revision date: 27.01.2017

Product code: 042-CA

Page 1 of 7

1. Identification

Product identifier

VITA SPRAY-ON INDICATOR LIQUID

Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Use as laboratory reagent

Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H. Rauter GmbH & Co. KG

Post-office box: 1338

79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

Emergency telephone number: +49-(0)761-19240

2. Hazard identification

Classification of the substance or mixture

WHMIS 2015

Hazard categories:

Flammable liquid: Flam. Liq. 2

Serious eye damage/eye irritation: Eye Irrit. 2A

Hazard Statements:

Highly flammable liquid and vapour.

Causes serious eye irritation.

Label elements

WHMIS 2015

Signal word: Danger**Pictograms:**

Hazard statements

Highly flammable liquid and vapour.

Causes serious eye irritation.

Precautionary statements

Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

Keep container tightly closed.

Wear protective gloves/protective clothing/eye protection/face protection.

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Store in a well-ventilated place. Keep cool.

Other hazards

No information available.

3. Composition/information on ingredients

Mixtures

Safety Data Sheet

according to WHMIS

VITA SPRAY-ON INDICATOR LIQUID

Revision date: 27.01.2017

Product code: 042-CA

Page 2 of 7

Hazardous components

CAS No	Chemical name	Quantity
64-17-5	ethyl alcohol	85 - < 90 %
	Spray opaque liquid (USA) Alcohol <2 % Proprietary (not hazardous)geschützte Lösemittelmischung (kein Gefahrstoff) > 98%	10 - < 15 %

4. First-aid measures

Description of first aid measures

After inhalation

In case of inhalation of decomposition products, affected person should be moved into fresh air and kept still.

After contact with skin

After contact with skin, wash immediately with plenty of water and soap.

After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

After ingestion

Rinse mouth immediately and drink plenty of water.

Indication of immediate medical attention and special treatment needed

Treat symptomatically.

5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media

Carbon dioxide (CO₂), Foam, Extinguishing powder.

Unsuitable extinguishing media

Water.

Specific hazards arising from the hazardous product

Highly flammable. Vapours can form explosive mixtures with air.

Special protective equipment and precautions for fire-fighters

In case of fire: Wear self-contained breathing apparatus.

Additional information

Use water spray jet to protect personnel and to cool endangered containers. Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Remove all sources of ignition. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

Environmental precautions

Do not allow uncontrolled discharge of product into the environment. Danger of explosion

Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Safety Data Sheet

according to WHMIS

VITA SPRAY-ON INDICATOR LIQUID

Revision date: 27.01.2017

Product code: 042-CA

Page 3 of 7

Disposal: see section 13

7. Handling and storage

Precautions for safe handling

Advice on safe handling

No special measures are necessary.

Advice on protection against fire and explosion

Keep away from sources of ignition. - No smoking. Take precautionary measures against static discharges. Vapours can form explosive mixtures with air.

Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container tightly closed. Keep in a cool, well-ventilated place. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

Advice on storage compatibility

Do not store together with: Oxidising agent. Pyrophoric or self-heating substances.

8. Exposure controls/Personal protection

Control parameters

Exposure limits (ACGIH)

CAS No	Chemical name	ppm	mg/m ³	F/ml	Category	Origin
64-17-5	Ethanol	-	-		TWA (8 h)	ACGIH-2016
		1000			STEL (15 min)	ACGIH-2016

Exposure controls



Protective and hygiene measures

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

Eye/face protection

Suitable eye protection: goggles.

Hand protection

Recommended glove articles KCL Dermatrill P NBR (Nitrile rubber) Breakthrough time (maximum wearing time) 30 min When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

Skin protection

Wear suitable protective clothing.

Respiratory protection

Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

Safety Data Sheet

according to WHMIS

VITA SPRAY-ON INDICATOR LIQUID

Revision date: 27.01.2017

Product code: 042-CA

Page 4 of 7

9. Physical and chemical properties

Information on basic physical and chemical properties

Physical state: liquid
 Colour: pink
 Odour: characteristic

Test method

pH-Value: not determined

Changes in the physical state

Melting point: not determined
 Initial boiling point and boiling range: 78 °C
 Flash point: 15 °C

Flammability

Solid: not applicable
 Gas: not applicable
 Lower explosive limits: 3,5 vol. %
 Upper explosive limits: 15 vol. %
 Ignition temperature: 425 °C

Auto-ignition temperature

Solid: not applicable
 Gas: not applicable
 Decomposition temperature: not determined

Oxidizing properties

Not oxidizing.
 Vapour pressure: <=1100 hPa
 (at 50 °C)
 Density: not determined

Solubility in other solvents

not determined
 Partition coefficient: not determined
 Vapour density: not determined
 Evaporation rate: not determined

Other information

Solid content: 0 %

10. Stability and reactivity

Conditions to avoid

Keep away from sources of heat (e.g. hot surfaces), sparks and open flames. Vapours can form explosive mixtures with air.

Incompatible materials

No information available.

Hazardous decomposition products

No known hazardous decomposition products.

Safety Data Sheet

according to WHMIS

VITA SPRAY-ON INDICATOR LIQUID

Revision date: 27.01.2017

Product code: 042-CA

Page 5 of 7

11. Toxicological information

Information on toxicological effects

Acute toxicity

Based on available data, the classification criteria are not met.

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
64-17-5	ethyl alcohol				
	dermal	LD50 7060 mg/kg			
	Spray opaque liquid (USA) Alcohol <2 % Proprietary (not hazardous)geschützte Lösemittelmischung (kein Gefahrstoff) > 98%				
	oral	LD50 20 mg/kg	Ratte		
	dermal	LD50 20,8 mg/kg	Kaninchen		
	inhalative	Data lacking			

Irritation and corrosivity

Causes serious eye irritation.

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Sensitizing effects

Based on available data, the classification criteria are not met.

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Additional information on tests

This mixture is classified as hazardous according to regulation (EC) No. 1272/2008 [CLP].

12. Ecological information

Ecotoxicity

The product is not: Ecotoxic.

Persistence and degradability

The product has not been tested.

Bioaccumulative potential

The product has not been tested.

Mobility in soil

The product has not been tested.

Other adverse effects

No information available.

Further information

Avoid release to the environment.

13. Disposal considerations

Safety Data Sheet

according to WHMIS

VITA SPRAY-ON INDICATOR LIQUID

Revision date: 27.01.2017

Product code: 042-CA

Page 6 of 7

Waste treatment methods

Advice on disposal

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

Contaminated packaging

Wash with plenty of water. Handle contaminated packages in the same way as the substance itself.

14. Transport information

Marine transport (IMDG)

UN number: UN 1993
United Nations proper shipping name: FLAMMABLE LIQUID, N.O.S. (ethyl alcohol)
Transport hazard class(es): 3
Packing group: II
 Hazard label: 3



Limited quantity: 1 L
 Excepted quantity: E2
 EmS: F-E, S-E

Air transport (ICAO-TI/IATA-DGR)

UN number: UN 1993
United Nations proper shipping name: FLAMMABLE LIQUID, N.O.S. (ethyl alcohol)
Transport hazard class(es): 3
Packing group: II
 Hazard label: 3



Limited quantity Passenger: 1 L
 Passenger LQ: Y341
 Excepted quantity: E2
 IATA-packing instructions - Passenger: 353
 IATA-max. quantity - Passenger: 5 L
 IATA-packing instructions - Cargo: 364
 IATA-max. quantity - Cargo: 60 L

Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

15. Regulatory information

Canadian regulations

16. Other information

Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product

Safety Data Sheet

according to WHMIS

VITA SPRAY-ON INDICATOR LIQUID

Revision date: 27.01.2017

Product code: 042-CA

Page 7 of 7

named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)