

Safety Data Sheet

according to 29 CFR 1910.1200(g)

VITAVM LC MODELLING LIQUID

Revision date: 16.02.2017

Product code: 151-US

Page 1 of 7

1. Identification

Product identifier

VITAVM LC MODELLING LIQUID

Recommended use of the chemical and restrictions on use

Use of the substance/mixture

Use as laboratory reagent

Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H. Rauter GmbH & Co. KG

Post-office box: 1338

79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

Emergency phone number: +49-(0)761-19240

Further Information

medical device

2. Hazard(s) identification

Classification of the chemical

29 CFR Part 1910.1200

Hazard categories:

Skin corrosion/irritation: Skin Irrit. 2

Serious eye damage/eye irritation: Eye Irrit. 2A

Respiratory or skin sensitization: Skin Sens. 1

Specific target organ toxicity single exposure: STOT SE 3

Hazard Statements:

Causes skin irritation

Causes serious eye irritation

May cause an allergic skin reaction

May cause respiratory irritation

Label elements

29 CFR Part 1910.1200

Signal word: Warning**Pictograms:**

Hazard statements

Causes skin irritation

May cause an allergic skin reaction

Causes serious eye irritation

May cause respiratory irritation

Precautionary statements

Avoid breathing dust/fume/gas/mist/vapors/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

If on skin: Wash with plenty of water.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

Safety Data Sheet

according to 29 CFR 1910.1200(g)

VITAVM LC MODELLING LIQUID

Revision date: 16.02.2017

Product code: 151-US

Page 2 of 7

Hazards not otherwise classified

No information available.

3. Composition/information on ingredients

Mixtures

Hazardous components

| CAS No | Components | Quantity |
|------------|-----------------------------------|----------|
| 109-16-0 | triethylene glycol dimethacrylate | 79.411 % |
| 90551-76-1 | Methacrylic ester | 19.803 % |
| 2867-47-2 | 2-dimethylaminoethyl methacrylate | 0.63 % |

4. First-aid measures

Description of first aid measures

After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

After contact with skin

After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

After ingestion

Rinse mouth immediately and drink plenty of water.

Most important symptoms and effects, both acute and delayed

No information available.

Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

Specific hazards arising from the chemical

Non-flammable.

Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately.

Do not allow entering drains or surface water.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

Environmental precautions

Do not allow to enter into surface water or drains.

Safety Data Sheet

according to 29 CFR 1910.1200(g)

VITAVM LC MODELLING LIQUID

Revision date: 16.02.2017

Product code: 151-US

Page 3 of 7

Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

7. Handling and storage

Precautions for safe handling

Advice on safe handling

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

Advice on protection against fire and explosion

No special fire protection measures are necessary.

Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container tightly closed.

Advice on storage compatibility

No special measures are necessary.

8. Exposure controls/personal protection

Control parameters

Exposure controls



Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

Protective and hygiene measures

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

Eye/face protection

Suitable eye protection: goggles.

Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. NBR (Nitrile rubber) Recommended glove articles KCL Camtril Velour Breakthrough time (maximum wearing time) 30 min

Skin protection

Wear suitable protective clothing.

Respiratory protection

Provide adequate ventilation as well as local exhaustion at critical locations. Technical ventilation of

Safety Data Sheet

according to 29 CFR 1910.1200(g)

VITAVM LC MODELLING LIQUID

Revision date: 16.02.2017

Product code: 151-US

Page 4 of 7

workplace

9. Physical and chemical properties

Information on basic physical and chemical properties

Physical state: liquid
 Color: colourless
 Odor: characteristic

Test method

pH-Value: not determined

Changes in the physical state

Melting point/freezing point: not determined

Initial boiling point and boiling range: ?

Flash point: ?

Flammability

Solid: not applicable

Gas: not applicable

Lower explosion limits: not determined

Upper explosion limits: not determined

Auto-ignition temperature

Solid: not applicable

Gas: not applicable

Decomposition temperature: not determined

Oxidizing properties

Not oxidizing.

Vapor pressure: <=1100 hPa
 (at 50 °C)

Density: 1,06000 g/cm³

Water solubility: No

Solubility in other solvents

not determined

Partition coefficient: not determined

Vapor density: not determined

Evaporation rate: not determined

Other information

Solid content: 0,0 %

10. Stability and reactivity

Reactivity

No hazardous reaction when handled and stored according to provisions.

Chemical stability

The product is stable under storage at normal ambient temperatures.

Possibility of hazardous reactions

No known hazardous reactions.

Safety Data Sheet

according to 29 CFR 1910.1200(g)

VITAVM LC MODELLING LIQUID

Revision date: 16.02.2017

Product code: 151-US

Page 5 of 7

Conditions to avoid

none

Incompatible materials

No information available.

Hazardous decomposition products

No known hazardous decomposition products.

11. Toxicological information

Information on toxicological effects

Acute toxicity

Based on available data, the classification criteria are not met.

| CAS No | Components | | | | |
|------------|-----------------------------------|---------------|---------|--------|--------|
| | Exposure route | Dose | Species | Source | Method |
| 90551-76-1 | Methacrylic ester | | | | |
| | dermal | LD50 mg/kg | 3000 | | |
| 2867-47-2 | 2-dimethylaminoethyl methacrylate | | | | |
| | oral | ATE mg/kg | 500 | | |
| | dermal | ATE mg/kg | 1100 | | |

Irritation and corrosivity

Causes skin irritation
Causes serious eye irritation

Sensitizing effects

May cause an allergic skin reaction (triethylene glycol dimethacrylate; 2-dimethylaminoethyl methacrylate)

Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

Specific target organ toxicity (STOT) - single exposure

May cause respiratory irritation (triethylene glycol dimethacrylate)

Specific target organ toxicity (STOT) - repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Additional information on tests

This mixture is classified as hazardous according to regulation (EC) No. 1272/2008 [CLP].

12. Ecological information

Ecotoxicity

The product is not: Ecotoxic.

Persistence and degradability

The product has not been tested.

Bioaccumulative potential

The product has not been tested.

Mobility in soil

The product has not been tested.

Safety Data Sheet

according to 29 CFR 1910.1200(g)

VITAVM LC MODELLING LIQUID

Revision date: 16.02.2017

Product code: 151-US

Page 6 of 7

Other adverse effects

No information available.

Further information

Avoid release to the environment.

13. Disposal considerations

Waste treatment methods

Advice on disposal

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

14. Transport information

Marine transport (IMDG)

| | |
|------------------------------------|--|
| <u>UN number:</u> | No dangerous good in sense of this transport regulation. |
| <u>UN proper shipping name:</u> | No dangerous good in sense of this transport regulation. |
| <u>Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| <u>Packing group:</u> | No dangerous good in sense of this transport regulation. |

Air transport (ICAO-TI/IATA-DGR)

| | |
|------------------------------------|--|
| <u>UN number:</u> | No dangerous good in sense of this transport regulation. |
| <u>UN proper shipping name:</u> | No dangerous good in sense of this transport regulation. |
| <u>Transport hazard class(es):</u> | No dangerous good in sense of this transport regulation. |
| <u>Packing group:</u> | No dangerous good in sense of this transport regulation. |

Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

Special precautions for user

No information available.

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

not applicable

15. Regulatory information

U.S. Regulations

National regulatory information

SARA Section 311/312 Hazards:

- triethylene glycol dimethacrylate (109-16-0): Immediate (acute) health hazard
- Methacrylic ester (90551-76-1): Immediate (acute) health hazard
- 2-dimethylaminoethyl methacrylate (2867-47-2): Immediate (acute) health hazard

State Regulations

Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)

This product contains no chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

16. Other information

Revision date: 16.02.2017
Revision No: 1

Safety Data Sheet

according to 29 CFR 1910.1200(g)

VITAVM LC MODELLING LIQUID

Revision date: 16.02.2017

Product code: 151-US

Page 7 of 7

Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route
(European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
IATA: International Air Transport Association
GHS: Globally Harmonized System of Classification and Labelling of Chemicals
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service
LC50: Lethal concentration, 50%
LD50: Lethal dose, 50%

Other data

The information is based on present level of our knowledge. It does not, however, give assurances of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)