

Safety Data Sheet

according to WHMIS

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 01.03.2017

Product code: 218-CA

Page 1 of 7

1. Identification

Product identifier

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Use as laboratory reagent

Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H. Rauter GmbH & Co. KG

Post-office box: 1338

79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

Emergency telephone number: +49-(0)761-19240

Further Information

medical device

2. Hazard identification

Classification of the substance or mixture

WHMIS 2015

Hazard categories:

Skin corrosion/irritation: Skin Irrit. 2

Serious eye damage/eye irritation: Eye Irrit. 2A

Respiratory or skin sensitization: Skin Sens. 1

Reproductive toxicity: Repr. 2

Hazard Statements:

Causes skin irritation.


Causes serious eye irritation.

May cause an allergic skin reaction.

Suspected of damaging fertility or the unborn child.

Label elements

WHMIS 2015

Signal word: Warning**Pictograms:**

Hazard statements

Causes skin irritation.

May cause an allergic skin reaction.

Causes serious eye irritation.

Suspected of damaging fertility or the unborn child.

Precautionary statements

Obtain special instructions before use.

Avoid breathing dust/fume/gas/mist/vapours/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Safety Data Sheet

according to WHMIS

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 01.03.2017

Product code: 218-CA

Page 2 of 7

Other hazards

No information available.

3. Composition/information on ingredients

Mixtures

Hazardous components

CAS No	Chemical name	Quantity
94108-97-1	Ditrimethylolpropane Tetraacrylate	5 - < 10 %
2867-47-2	2-dimethylaminoethyl methacrylate	5 - < 10 %
10373-78-1	Camphorquinone	5 - < 10 %
75980-60-8	diphenyl(2,4,6- trimethylbenzoyl)phosphine oxide	1 - < 5 %

4. First-aid measures

Description of first aid measures

After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

After contact with skin

After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

After ingestion

Rinse mouth immediately and drink plenty of water.

Indication of immediate medical attention and special treatment needed

Treat symptomatically.

5. Fire-fighting measures

Extinguishing media

Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

Specific hazards arising from the hazardous product

Non-flammable.

Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

6. Accidental release measures

Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

Environmental precautions

Do not allow to enter into surface water or drains.

Safety Data Sheet

according to WHMIS

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 01.03.2017

Product code: 218-CA

Page 3 of 7

Methods and material for containment and cleaning up

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

7. Handling and storage

Precautions for safe handling

Advice on safe handling

No special measures are necessary.

Advice on protection against fire and explosion

No special fire protection measures are necessary.

Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

Keep container tightly closed.

Advice on storage compatibility

No special measures are necessary.

8. Exposure controls/Personal protection

Control parameters

Exposure controls



Protective and hygiene measures

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

Eye/face protection

Suitable eye protection: goggles.

Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL DermatrILP NBR (Nitrile rubber) Breakthrough time (maximum wearing time) 480 min

Skin protection

Wear suitable protective clothing.

Respiratory protection

Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

Safety Data Sheet

according to WHMIS

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 01.03.2017

Product code: 218-CA

Page 4 of 7

9. Physical and chemical properties

Information on basic physical and chemical properties

Physical state:

Colour:

pH-Value:

not determined

Test method

Changes in the physical state

Melting point:

not determined

Initial boiling point and boiling range:

?

Flash point:

151 °C

Flammability

Solid:

not determined

Gas:

not applicable

Lower explosive limits:

not determined

Upper explosive limits:

not determined

Auto-ignition temperature

Solid:

not determined

Gas:

not applicable

Decomposition temperature:

not determined

Oxidizing properties

Not oxidizing.

Vapour pressure:

<=1100 hPa

(at 50 °C)

Density:

not determined

Water solubility:

No

Solubility in other solvents

not determined

Partition coefficient:

not determined

Vapour density:

not determined

Evaporation rate:

not determined

Other information

Solid content:

49,4 %

10. Stability and reactivity

Conditions to avoid

none

Incompatible materials

No information available.

Hazardous decomposition products

No known hazardous decomposition products.

11. Toxicological information

Information on toxicological effects

Safety Data Sheet

according to WHMIS

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 01.03.2017

Product code: 218-CA

Page 5 of 7

Acute toxicity

Based on available data, the classification criteria are not met.

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE mg/kg	500		
	dermal	ATE mg/kg	1100		

Irritation and corrosivity

Causes skin irritation.

Causes serious eye irritation.

Sensitizing effects

May cause an allergic skin reaction. (2-dimethylaminoethyl methacrylate)

Carcinogenic/mutagenic/toxic effects for reproduction

Suspected of damaging fertility or the unborn child. (diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide)

Germ cell mutagenicity: Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

STOT-single exposure

Based on available data, the classification criteria are not met.

STOT-repeated exposure

Based on available data, the classification criteria are not met.

Aspiration hazard

Based on available data, the classification criteria are not met.

Additional information on tests

This mixture is classified as hazardous according to regulation (EC) No. 1272/2008 [CLP].

12. Ecological information

Ecotoxicity

Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Persistence and degradability

The product has not been tested.

Bioaccumulative potential

The product has not been tested.

Mobility in soil

The product has not been tested.

Other adverse effects

No information available.

Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

13. Disposal considerations

Waste treatment methods

Advice on disposal

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

Safety Data Sheet

according to WHMIS

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 01.03.2017

Product code: 218-CA


Page 6 of 7

Contaminated packaging


Handle contaminated packages in the same way as the substance itself. Waste codes/waste designations according to EWC/AVV

14. Transport information

Marine transport (IMDG)

UN number:	UN 3243
United Nations proper shipping name:	SOLIDS CONTAINING TOXIC LIQUID, N.O.S. ((1-methylethylidene)bis(4,1-phenyleneoxy(2-hydroxy-3,1-propanediyl))bismethacrylate; Aliphatic urethane acrylate)
Transport hazard class(es):	6.1
Packing group:	II
Hazard label:	6.1
	
Limited quantity:	500 g
Excepted quantity:	E4
EmS:	F-A, S-A

Air transport (ICAO-TI/IATA-DGR)

UN number:	UN 3243
United Nations proper shipping name:	SOLIDS CONTAINING TOXIC LIQUID, N.O.S. ((1-methylethylidene)bis(4,1-phenyleneoxy(2-hydroxy-3,1-propanediyl))bismethacrylate; Aliphatic urethane acrylate)
Transport hazard class(es):	6.1
Packing group:	II
Hazard label:	6.1
	
Limited quantity Passenger:	1 kg
Passenger LQ:	Y644
Excepted quantity:	E4
IATA-packing instructions - Passenger:	669
IATA-max. quantity - Passenger:	25 kg
IATA-packing instructions - Cargo:	676
IATA-max. quantity - Cargo:	100 kg

Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

15. Regulatory information

Canadian regulations

16. Other information

Further Information

The information is based on present level of our knowledge. It does not, however, give assurances of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for

Safety Data Sheet

according to WHMIS

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 01.03.2017

Product code: 218-CA

Page 7 of 7

adhering to existing laws and regulations.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)