

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITAVM LC PAINT

Revision date: 23.02.2017

Product code: 162-US

Page 1 of 7

## 1. Identification

### Product identifier

VITAVM LC PAINT

### Recommended use of the chemical and restrictions on use

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H. Rauter GmbH &amp; Co. KG

Post-office box: 1338

79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

**Emergency phone number:** +49-(0)761-19240

### Further Information

medical device

## 2. Hazard(s) identification

### Classification of the chemical

#### 29 CFR Part 1910.1200

Hazard categories:

Skin corrosion/irritation: Skin Irrit. 2

Serious eye damage/eye irritation: Eye Irrit. 2A

Respiratory or skin sensitization: Skin Sens. 1

Hazard Statements:

Causes skin irritation

Causes serious eye irritation

May cause an allergic skin reaction

### Label elements

#### 29 CFR Part 1910.1200

**Signal word:** Warning**Pictograms:**

### Hazard statements

Causes skin irritation

May cause an allergic skin reaction

Causes serious eye irritation

### Precautionary statements

Avoid breathing dust/fume/gas/mist/vapors/spray.

Contaminated work clothing must not be allowed out of the workplace.

Wear protective gloves/protective clothing/eye protection/face protection.

If on skin: Wash with plenty of water.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

If eye irritation persists: Get medical advice/attention.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITAVM LC PAINT

Revision date: 23.02.2017

Product code: 162-US

Page 2 of 7

#### Hazards not otherwise classified

No information available.

### 3. Composition/information on ingredients

#### Mixtures

#### Hazardous components

CAS No	Components	Quantity
109-16-0	triethylene glycol dimethacrylate	17 %
2867-47-2	2-dimethylaminoethyl methacrylate	0.2 %

### 4. First-aid measures

#### Description of first aid measures

##### **After inhalation**

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

##### **After contact with skin**

After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary. After contact with skin, wash immediately with plenty of water and soap.

##### **After contact with eyes**

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

##### **After ingestion**

Rinse mouth immediately and drink plenty of water.

#### Most important symptoms and effects, both acute and delayed

No information available.

#### Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### 5. Fire-fighting measures

#### Extinguishing media

##### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

#### Specific hazards arising from the chemical

Non-flammable.

#### Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

#### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

### 6. Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

#### Environmental precautions

Do not allow to enter into surface water or drains.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITAVM LC PAINT

Revision date: 23.02.2017

Product code: 162-US

Page 3 of 7

#### Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## 7. Handling and storage

#### Precautions for safe handling

##### **Advice on safe handling**

No special measures are necessary.

##### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

#### Conditions for safe storage, including any incompatibilities

##### **Requirements for storage rooms and vessels**

Keep container tightly closed.

##### **Advice on storage compatibility**

No special measures are necessary.

## 8. Exposure controls/personal protection

#### Control parameters

#### Exposure controls



#### **Protective and hygiene measures**

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

#### **Eye/face protection**

Suitable eye protection: goggles.

#### **Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Camatril Velour Breakthrough time (maximum wearing time) 30 min NBR (Nitrile rubber)

#### **Skin protection**

Wear suitable protective clothing.

#### **Respiratory protection**

Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITAVM LC PAINT

Revision date: 23.02.2017

Product code: 162-US

Page 4 of 7

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state:

Color:

Odor: characteristic

#### Test method

pH-Value: not determined

### Changes in the physical state

Melting point/freezing point: not determined

Initial boiling point and boiling range: X

Flash point: X

### Flammability

Solid: not applicable

Gas: not applicable

Lower explosion limits: not determined

Upper explosion limits: not determined

### Auto-ignition temperature

Solid: not applicable

Gas: not applicable

Decomposition temperature: not determined

### Oxidizing properties

Not oxidizing.

Vapor pressure: <=1100 hPa  
(at 50 °C)

Density: not determined

Water solubility: No

### Solubility in other solvents

not determined

Partition coefficient: not determined

Vapor density: not determined

Evaporation rate: not determined

### Other information

Solid content: 1,10 %

## 10. Stability and reactivity

### Reactivity

No hazardous reaction when handled and stored according to provisions.

### Chemical stability

The product is stable under storage at normal ambient temperatures.

### Possibility of hazardous reactions

No known hazardous reactions.

### Conditions to avoid

none

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITAVM LC PAINT

Revision date: 23.02.2017

Product code: 162-US

Page 5 of 7

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

### Information on toxicological effects

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

CAS No	Components				
	Exposure route	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE mg/kg	500		
	dermal	ATE mg/kg	1100		

#### **Irritation and corrosivity**

Causes skin irritation  
Causes serious eye irritation

#### **Sensitizing effects**

May cause an allergic skin reaction (triethylene glycol dimethacrylate; 2-dimethylaminoethyl methacrylate)

#### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

#### **Specific target organ toxicity (STOT) - single exposure**

Based on available data, the classification criteria are not met.

#### **Specific target organ toxicity (STOT) - repeated exposure**

Based on available data, the classification criteria are not met.

Carcinogenicity (OSHA): No ingredient of this mixture is listed.

Carcinogenicity (IARC): No ingredient of this mixture is listed.

Carcinogenicity (NTP): No ingredient of this mixture is listed.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### **Additional information on tests**

This mixture is classified as hazardous according to regulation (EC) No. 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

#### Mobility in soil

The product has not been tested.

#### Other adverse effects

No information available.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITAVM LC PAINT

Revision date: 23.02.2017

Product code: 162-US

Page 6 of 7

#### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

### 13. Disposal considerations

#### Waste treatment methods

##### Advice on disposal

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

### 14. Transport information

#### Marine transport (IMDG)

<b>UN number:</b>	No dangerous good in sense of this transport regulation.
<b>UN proper shipping name:</b>	No dangerous good in sense of this transport regulation.
<b>Transport hazard class(es):</b>	No dangerous good in sense of this transport regulation.
<b>Packing group:</b>	No dangerous good in sense of this transport regulation.

#### Air transport (ICAO-TI/IATA-DGR)

<b>UN number:</b>	No dangerous good in sense of this transport regulation.
<b>UN proper shipping name:</b>	No dangerous good in sense of this transport regulation.
<b>Transport hazard class(es):</b>	No dangerous good in sense of this transport regulation.
<b>Packing group:</b>	No dangerous good in sense of this transport regulation.

#### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

#### Special precautions for user

No information available.

#### Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

not applicable

### 15. Regulatory information

#### U.S. Regulations

##### National regulatory information

SARA Section 311/312 Hazards:

triethylene glycol dimethacrylate (109-16-0): Immediate (acute) health hazard

2-dimethylaminoethyl methacrylate (2867-47-2): Immediate (acute) health hazard

##### State Regulations

##### Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)

This product contains no chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

### 16. Other information

Revision date: 23.02.2017

Revision No: 1

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITAVM LC PAINT

Revision date: 23.02.2017

Product code: 162-US

Page 7 of 7

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

#### Other data

The information is based on present level of our knowledge. It does not, however, give assurances of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*