

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC PRE OPAQUE

Revision date: 03.03.2017

Product code: 263-US

Page 1 of 7

## 1. Identification

### Product identifier

VITA VM LC PRE OPAQUE

### Recommended use of the chemical and restrictions on use

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H. Rauter GmbH &amp; Co. KG

Post-office box: 1338

79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

**Emergency phone number:** +49-(0)761-19240

### Further Information

medical device

## 2. Hazard(s) identification

### Classification of the chemical

#### 29 CFR Part 1910.1200

Hazard categories:

Serious eye damage/eye irritation: Eye Irrit. 2A

Respiratory or skin sensitization: Skin Sens. 1

Hazard Statements:

Causes serious eye irritation

May cause an allergic skin reaction

### Label elements

#### 29 CFR Part 1910.1200

**Signal word:** Warning**Pictograms:**

### Hazard statements

May cause an allergic skin reaction

Causes serious eye irritation

### Precautionary statements

Avoid breathing dust/fume/gas/mist/vapors/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

If on skin: Wash with plenty of water.

If skin irritation or rash occurs: Get medical advice/attention.

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

If eye irritation persists: Get medical advice/attention.

### Hazards not otherwise classified

No information available.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC PRE OPAQUE

Revision date: 03.03.2017

Product code: 263-US

Page 2 of 7

### 3. Composition/information on ingredients

#### Mixtures

#### Hazardous components

CAS No	Components	Quantity
94108-97-1	Ditrimethylolpropane Tetraacrylate	11.4 %
2867-47-2	2-dimethylaminoethyl methacrylate	0.46 %

### 4. First-aid measures

#### Description of first aid measures

##### **After inhalation**

Provide fresh air. If breathing is irregular or stopped, administer artificial respiration. Medical treatment necessary.

##### **After contact with skin**

Wash with plenty of water. Take off immediately all contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention.

##### **After contact with eyes**

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

##### **After ingestion**

Rinse mouth immediately and drink plenty of water.

#### Most important symptoms and effects, both acute and delayed

No information available.

#### Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### 5. Fire-fighting measures

#### Extinguishing media

##### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

#### Specific hazards arising from the chemical

Non-flammable.

#### Special protective equipment and precautions for fire-fighters

In case of fire: Wear self-contained breathing apparatus.

#### **Additional information**

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

### 6. Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

#### Environmental precautions

Do not allow to enter into surface water or drains.

#### Methods and material for containment and cleaning up

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC PRE OPAQUE

Revision date: 03.03.2017

Product code: 263-US

Page 3 of 7

#### Reference to other sections

Safe handling: see section 7  
Personal protection equipment: see section 8  
Disposal: see section 13

## 7. Handling and storage

#### Precautions for safe handling

##### **Advice on safe handling**

No special measures are necessary.

##### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

#### Conditions for safe storage, including any incompatibilities

##### **Requirements for storage rooms and vessels**

Keep container tightly closed.

##### **Advice on storage compatibility**

No special measures are necessary.

## 8. Exposure controls/personal protection

#### Control parameters

#### Exposure controls



#### **Protective and hygiene measures**

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

#### **Eye/face protection**

Suitable eye protection: goggles.

#### **Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatrill P Breakthrough time (maximum wearing time) 30 min NBR (Nitrile rubber)

#### **Skin protection**

Wear suitable protective clothing.

#### **Respiratory protection**

Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

## 9. Physical and chemical properties

#### Information on basic physical and chemical properties

Physical state:

Color:

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC PRE OPAQUE

Revision date: 03.03.2017

Product code: 263-US

Page 4 of 7

#### Test method

pH-Value: not determined

#### Changes in the physical state

Melting point/freezing point: not determined

Initial boiling point and boiling range: ?

Flash point: 151 °C

#### Flammability

Solid: not determined

Gas: not applicable

Lower explosion limits: not determined

Upper explosion limits: not determined

#### Auto-ignition temperature

Solid: not determined

Gas: not applicable

Decomposition temperature: not determined

#### Oxidizing properties

Not oxidizing.

Vapor pressure: <=1100 hPa  
(at 50 °C)

Density: not determined

Water solubility: No

#### Solubility in other solvents

not determined

Partition coefficient: not determined

Vapor density: not determined

Evaporation rate: not determined

#### Other information

Solid content: 49,7 %

## 10. Stability and reactivity

#### Reactivity

No hazardous reaction when handled and stored according to provisions.

#### Chemical stability

The product is stable under storage at normal ambient temperatures.

#### Possibility of hazardous reactions

No known hazardous reactions.

#### Conditions to avoid

none

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC PRE OPAQUE

Revision date: 03.03.2017

Product code: 263-US

Page 5 of 7

#### Information on toxicological effects

##### Acute toxicity

Based on available data, the classification criteria are not met.

CAS No	Components				
	Exposure route	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE mg/kg	500		
	dermal	ATE mg/kg	1100		

##### Irritation and corrosivity

Causes serious eye irritation

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

##### Sensitizing effects

May cause an allergic skin reaction (2-dimethylaminoethyl methacrylate)

##### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

##### Specific target organ toxicity (STOT) - single exposure

Based on available data, the classification criteria are not met.

##### Specific target organ toxicity (STOT) - repeated exposure

Based on available data, the classification criteria are not met.

##### Aspiration hazard

Based on available data, the classification criteria are not met.

##### Additional information on tests

This mixture is classified as hazardous according to regulation (EC) No. 1272/2008 [CLP].

## 12. Ecological information

##### Ecotoxicity

Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

##### Persistence and degradability

The product has not been tested.

##### Bioaccumulative potential

The product has not been tested.

##### Mobility in soil

The product has not been tested.

##### Other adverse effects

No information available.

##### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

## 13. Disposal considerations

##### Waste treatment methods

##### Advice on disposal

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Handle contaminated packages in the same way as the substance itself. Waste codes/waste designations

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC PRE OPAQUE

Revision date: 03.03.2017

Product code: 263-US

Page 6 of 7

according to EWC/AVV

#### 14. Transport information

##### Marine transport (IMDG)

<b>UN number:</b>	No dangerous good in sense of this transport regulation.
<b>UN proper shipping name:</b>	No dangerous good in sense of this transport regulation.
<b>Transport hazard class(es):</b>	No dangerous good in sense of this transport regulation.
<b>Packing group:</b>	No dangerous good in sense of this transport regulation.

##### Air transport (ICAO-TI/IATA-DGR)

<b>UN number:</b>	No dangerous good in sense of this transport regulation.
<b>UN proper shipping name:</b>	No dangerous good in sense of this transport regulation.
<b>Transport hazard class(es):</b>	No dangerous good in sense of this transport regulation.
<b>Packing group:</b>	No dangerous good in sense of this transport regulation.

##### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

##### Special precautions for user

No information available.

##### Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

not applicable

#### 15. Regulatory information

##### U.S. Regulations

###### National regulatory information

SARA Section 311/312 Hazards:

Ditrimethylolpropane Tetraacrylate (94108-97-1): Immediate (acute) health hazard  
2-dimethylaminoethyl methacrylate (2867-47-2): Immediate (acute) health hazard

###### State Regulations

###### Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)

This product contains no chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

#### 16. Other information

Revision date: 03.03.2017

Revision No: 1

##### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%

##### Other data

The information is based on present level of our knowledge. It does not, however, give assurances of product

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC PRE OPAQUE

Revision date: 03.03.2017

Product code: 263-US

Page 7 of 7

properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*