# **VITA VIONIC® SOLUTIONS**

For digital prosthetics with open CAD/CAM systems



Date of issue: 12.17

VITA – perfect match.



#### **Product description**

VITA VIONIC SOLUTIONS enables the efficient fabrication of CAD/CAM wax setups. First, a time-saving, digital prosthetic setup is generated using CAD software. Second, a precise wax base or a fully formed wax setup made from VITA VIONIC WAX is produced from the denture design


via CAM technology. Once a successful try-in is completed, the denture is conventionally produced in the final stage. With VITA VIONIC, laboratories benefit from higher levels of precision while saving time.



With what?

- VITA VIONIC WAX in pink for a precise wax base
- VITA VIONIC WAX in white for a full-format wax setup

#### **BENEFITS**




**Efficient wax setup** 

Precise wax setup

Comfortable wax try-in

- Thanks to digital design and CAM production
- Thanks to precisely milled molds based on the digital construction
- Thanks to fabrication of complete CAD/CAM denture out of white wax

## **CAD DESIGN**



What?

• CAD software with VITA-specific, coordinated design options

What for?


#### Simple, patient-specific reconstruction of the function $\dots$

- through the integrated buccal and lingual occlusal concepts
- through the integrated skeletal classes (I + II)
- through the integrated occlusal positions (normal/cross bite)


What for?

• for the digital, patient-specific denture setup at the touch of a button

#### **WORKFLOW**



# **BASIC INFORMATION**



#### **Characteristics**

#### VITA VIONIC WAX

• are shape-stable, temperature-resistant synthetic-wax blanks

#### **Available shades**

- Pink
- White





## **Available geometries**

• Ø 98.4 mm/H 30 mm







## **SYSTEM COMPATIBILITY**




#### **CAD** technology partner for VITA Digital Prosthetics

• 3Shape AS


#### VITA VIONIC WAX is available for the following systems\*:

- CORiTEC line (imes-icore GmbH)
- DG SHAPE DWX line (DG SHAPE)



<sup>\*)</sup> Validated VITA system partner for the processing of the selected VITA VIONIC components.

# More information about VITA VIONIC SOLUTIONS is available at www.vita-zahnfabrik.com/prosthetics



Please note: Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product. Furthermore, our liability for the accuracy of this information is independent of the legal basis and, in as far as legally permissible, shall always be limited to the value as invoiced of the goods supplied, excluding value-added tax. In particular, as far as legally permissible, we do not assume any liability for loss of earnings, indirect damages, ensuing damages or for third-party claims against the purchaser. Claims for damages based on fault liability (culpa in contrahendo, breach of contract, unlawful acts, etc.) can only be made in the case of intent or gross negligence. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information: 12.17

After the publication of these information for use any previous versions become obsolete. The current version can be found at www.vita-zahnfabrik.com

VITA Zahnfabrik has been certified in accordance to the Medical Device Directive and the following products bear the CE mark C  $\, {\bf C}_{\,}\,$  0124 :

VITA VIONIC®WAX



VITA Zahnfabrik H. Rauter GmbH & Co.KG
Spitalgasse 3 · D-79713 Bad Säckingen · Germany
Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299
Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446
www.vita-zahnfabrik.com · info@vita-zahnfabrik.com
facebook.com/vita.zahnfabrik