

> VITA YZ® SOLUTIONS

> VITA YZ® XT – Extra Translucent ZrO₂

VITA YZ® ST – Super Translucent ZrO₂

VITA YZ® HT – High Translucent ZrO₂

VITA YZ® T – Translucent ZrO₂



> FACTS & FEATURES

VITA YZ® XT is

- an extra translucent zirconia and is available in non-colored, as well as monochromatic and multicolored versions.

You benefit from

- esthetic monolithic and partially veneered restorations with exceptional vibrancy, thanks to extremely high translucency.
- efficient fabrication, thanks to tooth-colored, multichromatic blanks with an integrated tooth shade gradient.

VITA YZ® XT is

- the solution for three-unit anterior bridges with natural shade effects.

Available shades:

- VITA YZ® XT^{White} non-colored, white
- VITA YZ® XT^{Color} precolored in 16 VITA classical A1–D4 shades
- VITA YZ® XT^{Multicolor} precolored in 7 VITA classical A1–D4 shades*
- VITA YZ® XT SHADE LIQUID in 16 VITA classical A1–D4 shades

Geometry sizes:

- VITA YZ® XT disc geometries: 98.4 mm in heights: 14 mm, 18 mm, 20 mm, 22 mm (exclusively for Multicolor), 25 mm

> SYSTEM COMPONENTS



VITA YZ® XT SHADE LIQUID
Coloring

VITA AKZENT® Plus
Characterization

* VITA YZ® XT^{Multicolor} versions are available beginning in summer 2018 in heights 14 mm, 18 mm and 22 mm, as well as in VITA classical A1–D4 shades A1, A2, A3, A3.5, B2, C2 and D2.

> VITA YZ® SOLUTIONS

VITA YZ® XT – Extra Translucent ZrO₂

> **VITA YZ® ST** – Super Translucent ZrO₂

VITA YZ® HT – High Translucent ZrO₂

VITA YZ® T – Translucent ZrO₂



> FACTS & FEATURES

VITA YZ® ST is

- a super translucent zirconia and is available in non-colored, as well as monochromatic and multicolored versions.

You benefit from

- efficient fabrication, thanks to tooth-colored blanks.
- high-strength, multi-unit bridge restorations with incredibly lifelike appearance.
- monolithic bridge restorations that preserve natural tooth structure and partially veneered bridge restorations.

VITA YZ® ST is

- the solution for high-strength bridge restorations in the posterior area.

Available shades:

- **VITA YZ® ST^{White}** non-colored, white
- **VITA YZ® ST^{Color}** precolored in 16 VITA classical A1–D4 shades
- **VITA YZ® ST^{Multicolor}** precolored in 7 VITA classical A1–D4 shades*
- **VITA YZ® ST SHADE LIQUID** in 16 VITA classical A1–D4 colors

Geometry sizes:

- **VITA YZ® ST disc geometries:**
98.4 mm in heights: 14 mm, 18 mm, 20 mm, 22 mm (exclusively for Multicolor), 25 mm

> SYSTEM COMPONENTS



VITA YZ® ST SHADE LIQUID
Coloring

VITA AKZENT® Plus
Characterization

* **VITA YZ® ST^{Multicolor}** versions are available beginning in summer 2018 in heights 14 mm, 18 mm and 22 mm, as well as in VITA classical A1–D4 shades A1, A2, A3, A3.5, B2, C2 and D2.

> VITA YZ® SOLUTIONS

VITA YZ® XT – Extra Translucent ZrO₂

VITA YZ® ST – Super Translucent ZrO₂

> **VITA YZ® HT** – High Translucent ZrO₂

VITA YZ® T – Translucent ZrO₂



> FACTS & FEATURES

VITA YZ® HT is

- a highly translucent zirconia with enormous strength and is available in non-colored and monochromatic versions.

You benefit from

- highly esthetic restorations, thanks to perfectly matched veneering ceramics.
- efficient shade reproduction, thanks to perfectly matched components (e.g., liquids/stains/veneering ceramics).

VITA YZ® HT is

- for fully veneered, four-unit substructures.

Available shades:

- **VITA YZ® HT^{White}** non-colored, white
- **VITA YZ® HT^{Color}** precolored in 3 VITA classical A1–D4 shades
- **VITA YZ® HT SHADE LIQUID** in 16 VITA classical A1–D4 shades

Geometry sizes:

- **VITA YZ® HT disc geometries***:
98.4 mm in the following heights: 12 mm, 14 mm, 16 mm, 18 mm, 20 mm, 25 mm
- **VITA YZ® HT block geometries***:
20/19, 40/19, 55/19, 65/25

> SYSTEM COMPONENTS



* The disc geometries 12 mm and 16 mm, as well as all block geometries, are only available in white.

> VITA YZ® SOLUTIONS

VITA YZ® XT – Extra Translucent ZrO₂

VITA YZ® ST – Super Translucent ZrO₂

VITA YZ® HT – High Translucent ZrO₂

> VITA YZ® T – Translucent ZrO₂



> FACTS & FEATURES

VITA YZ® T is

- a precision-fit, high-strength and translucent zirconia and is available in non-colored and monochromatic versions.

You benefit from

- accurate results, thanks to high edge stability and precision-fit sintering.
- high reliability, thanks to enormous strength and 15 years of clinical experience.
- efficient shade reproduction, thanks to perfectly matched veneering ceramic.

VITA YZ® T is

- the solution for high-strength, fully veneered bridge substructures.


Available shades:

- **VITA YZ® T^{White}** non-colored, white
- **VITA YZ® T^{Color}** precolored, LL1/light, LL2/medium, LL3/intense
- **VITA YZ® T COLORINIG LIQUID** light/pale, medium, intense, neutral

Geometry sizes:

- **VITA YZ® T disc geometries*:**
98.4 mm in the following heights: 12 mm, 14 mm, 16 mm, 18 mm, 20 mm, 25 mm
- **VITA YZ® T block geometries:**
20/19, 40/19, 55, 65/25

> SYSTEM COMPONENTS



VITA YZ® T COLORING LIQUID
Coloring

VITAVM®9
Individualization

VITA AKZENT® Plus
Characterization

* Sizes 12 mm and 16 mm are only available in white.