

# VITA LUMIN<sup>®</sup> VACUUM

VITA Teeth



VITA shade determination

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 08.18

VITA – perfect match.

**VITA**

Formenkarte · Keramikzähne


Mould Chart · Ceramic denture teeth

Carte de formes · Dents en céramique


Carta de formas · Dientes de cerámica

Tabella delle forme · Denti di ceramica


VITA LUMIN® VACUUM Upper anteriors




schaufelförmig  
ovoid  
en forme de pelle  
en forma de pala  
a forma di pala






dreieckig  
tapering  
triangulaire  
triangular  
triangolare






Lower anteriors




rechteckig  
 rectangular  
 rectangulaire  
 rectangular  
 rettangolare



quadratisch  
 square  
 carré  
 cuadrado  
 quadrato



# Upper and lower posteriors



## CUSPIFORM






## SYNOFORM



## VITA LUMIN® VACUUM Combination Chart

LUMIN® VACUUM Kombinationstabelle  
 LUMIN® VACUUM Combination Chart  
 Tableau de Combinaison LUMIN® VACUUM  
 Tabla de Combinación LUMIN® VACUUM  
 Tabella di Combinazione LUMIN® VACUUM

No.	 mm	No.	 mm	mm	
				No.	 mm
<b>180</b>	44.7 50.7	<b>639</b>	33.6	<b>41C</b> 29.5 31.8	<b>11S</b> 29.3 31.2
<b>182</b>	43.0 48.6	<b>539</b>	32.3	<b>40C</b> 28.3 29.5	<b>10S</b> 28.2 29.7
<b>185</b>	46.1 52.2	<b>730</b>	35.6	<b>44C</b> 32.3 34.3	<b>13S</b> 30.3 32.7
<b>187</b>	44.6 50.2	<b>760</b>	33.6	<b>41C</b> 29.5 31.8	<b>11S</b> 29.3 31.2
<b>190</b>	42.6 48.1	<b>539</b>	32.3	<b>40C</b> 28.3 29.5	<b>10S</b> 28.2 29.7
<b>513</b>	40.1 45.8	<b>539</b>	32.3	<b>40C</b> 28.3 29.5	<b>10S</b> 28.2 29.7
<b>565</b>	45.3 51.4	<b>730</b>	35.6	<b>44C</b> 32.3 34.3	<b>13S</b> 30.3 32.7
<b>612</b>	43.9 49.9	<b>539</b>	32.3	<b>41C</b> 29.5 31.8	<b>11S</b> 29.3 31.2
<b>622</b>	44.1 50.1	<b>639</b>	33.6	<b>41C</b> 29.5 31.8	<b>11S</b> 29.3 31.2
<b>627</b>	42.3 47.8	<b>760</b>	33.6	<b>40C</b> 28.3 29.5	<b>10S</b> 28.2 29.7
<b>632</b>	42.5 48.3	<b>730</b>	35.6	<b>41C</b> 29.5 31.8	<b>11S</b> 29.3 31.2
<b>722</b>	43.7 49.8	<b>639</b>	33.6	<b>41C</b> 29.5 31.8	<b>11S</b> 29.3 31.2
<b>732</b>	44.8 50.5	<b>760</b>	33.6	<b>41C</b> 29.5 31.8	<b>11S</b> 29.3 31.2
<b>734</b>	43.3 49.4	<b>760</b>	33.6	<b>41C</b> 29.5 31.8	<b>11S</b> 29.3 31.2
<b>737</b>	45.7 51.6	<b>730</b>	35.6	<b>44C</b> 32.3 34.3	<b>13S</b> 30.3 32.7
<b>742</b>	46.2 52.2	<b>740</b>	35.1	<b>44C</b> 32.3 34.3	<b>13S</b> 30.3 32.7
<b>743</b>	46.7 52.9	<b>740</b>	35.1	<b>44C</b> 32.3 34.3	<b>13S</b> 30.3 32.7
<b>773</b>	47.1 53.2	<b>639</b>	33.6	<b>44C</b> 32.3 34.3	<b>13S</b> 30.3 32.7
<b>843</b>	47.6 53.6	<b>639</b>	33.6	<b>44C</b> 32.3 34.3	<b>13S</b> 30.3 32.7
<b>854</b>	48.0 54.6	<b>55</b>	36.1	<b>44C</b> 32.3 34.3	<b>13S</b> 30.3 32.7

**VITA LUMIN VACUUM Anteriores und Posteriores sind in den VITA classical Farben A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3 erhältlich.**

Für weitere Farben fallen Lieferzeiten an.

**VITA LUMIN VACUUM Anteriores and Posteriores are available in the VITA classical shades A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3.**

Further colours are subject to delivery times.

**VITA LUMIN VACUUM Anteriores et Posteriores sont disponibles dans les teintes VITA classical A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3.**

Pour toutes les couleurs supplémentaires s'applique un délai de livraison.


**I denti VITA LUMIN VACUUM Anteriores e Posteriores sono disponibili in colorazione VITA classical A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3.**

Colori supplementari hanno un tempo di consegna.

**La gama VITA LUMIN VACUUM Anteriores y Posteriores está disponible en los colores VITA classical A1, A2, A3, A3.5, B2, B3, C1, C2, D2, D3.**

Otros colores tendrán plazo de entrega.

With the unique VITA SYSTEM 3D-MASTER, all natural tooth shades can be systematically determined and perfectly reproduced.



**Please note:** Our products must be used in accordance with the instructions for use. We accept no liability for any damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials and equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage. The VITA Modulbox is not necessarily a component of the product. Date of issue of this information: 08.18

After the publication of this information for use any previous versions become obsolete. The current version can be found at [www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com)

VITA Zahnfabrik has been certified in accordance to the Medical Device Directive and the following products bear the CE mark **CE** 0124 :

VITA LUMIN VACUUM Teeth fulfil the requirements of ISO 4824 : 1993.

# VITA

VITA Zahnfabrik H. Rauter GmbH & Co.KG  
Spitalgasse 3 · D-79713 Bad Säckingen · Germany  
Tel. +49(0)7761-562-0 · Fax +49(0)7761-562-299  
Hotline: Tel. +49(0)7761-562-222 · Fax +49(0)7761-562-446  
[www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) · [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com)  
[facebook.com/vita.zahnfabrik](https://www.facebook.com/vita.zahnfabrik)