

# VITA ADIVA® IA-CEM

Ultra opaque luting solution for implant restorations



VITA shade determination

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue: 2022-10

VITA – perfect match.

**VITA**

## Product description

- What?**
- VITA ADIVA IA-CEM is an ultra opaque, strongly-masking, dual-curing luting composite for implant prosthetics.
- What for?**
- It is used for definitive, adhesive bonding of (abutment) crowns and mesostructures on titanium adhesive bases, individual titanium and zirconia abutments as well as crowns on discolored tooth stumps.
- With what?**
- There are material-specific, matched bonding agents for the luting composite to ensure durable bonding of the restoration and the abutment.

## PRODUCT BENEFITS

### Reliable masking

- Reliable masking of titanium abutments, thanks to ultra opaque, tooth-colored luting composite.

### Wide range of applications

- The only ultra opaque luting composite suitable for intra- and extraoral use.

### Safe adhesive bond

- Durable bond, thanks to high adhesion values to titanium and (hybrid) ceramic restorative materials.

### Efficient luting

- Material consumption is approx. 30 % lower than for conventional mixing tips, thanks to the automix syringe with a short mixing tip\*.

### Time-efficient treatment

- The luting material ensures time-efficient treatment, since self curing is completed extraorally after 7 minutes and intraorally after 4 minutes.

## PRODUCT USE

- VITA ADIVA IA-CEM offers universal and reliable use for bonding of all restorations made of feldspar, glass and hybrid ceramics, zirconia and composite.



### VITA ADIVA® IA-CEM SET Prod. No. FAIACEMSET

Designation	Quantity	Available as	Prod. No.***
VITA ADIVA IA-CEM, shade: ultra opaque	5 ml	Syringe + 15 T-Mixers 1:1, brown	FAIACEMU05
VITA ADIVA CERA-ETCH	6 ml	Bottle	FACE6
VITA ADIVA C-PRIME	5 ml	Bottle	FACP5
VITA ADIVA ZR-PRIME	5 ml	Bottle	FAZRP5
VITA ADIVA IA-CEM ACCESSORIES: VITA ADIVA T-Mixer, VITA ADIVA Microbrush, VITA ADIVA mixing trays			

### VITA ADIVA® IA-CEM ACCESSORIES

Designation	Units	Prod. No.***
VITA ADIVA T-Mixer 1:1	15 pieces	FATM11
VITA ADIVA Microbrush	25 pieces	FAMB25
VITA ADIVA mixing trays	25 pieces	FAAT25

\*) Source: T-Mixer Themenblatt, Sulzer Mixpac AG, Haag, Switzerland

\*\*) Conditioning the titanium abutment with Monobond® Plus. Monobond® Plus is a registered trademark of Ivoclar Vivadent AG, Schaan, Liechtenstein.

\*\*\*) All VITA ADIVA components can be (re)ordered individually using the Prod. No. specified.

**VITA**