

## Safety Data Sheet

according to WHMIS

### VITA VM LC Chroma Plus EFFECT LINER

Revision date: 12.08.2019

Product code: 217

Page 1 of 7

## 1. Identification

### Product identifier

VITA VM LC Chroma Plus EFFECT LINER

### Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H.Rauter GmbH &amp; Co.KG

Post-office box: 1338  
79704 Bad Säckingen

Telephone: +49(0)7761-562-0

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

Telefax: +49(0)7761-562-299

### Emergency telephone number:

+49-(0)761-19240

### Further Information

medical device

## 2. Hazard identification

### Classification of the substance or mixture

#### WHMIS 2015

Respiratory or skin sensitization: Skin Sens. 1

### Label elements

#### WHMIS 2015

Signal word: Warning

Pictograms:



#### Hazard statements

May cause an allergic skin reaction.

#### Precautionary statements

Avoid breathing dust/fume/gas/mist/vapours/spray.

Contaminated work clothing should not be allowed out of the workplace.

Wear protective gloves/protective clothing/eye protection/face protection.

### Other hazards

No information available.

## 3. Composition/information on ingredients

### Mixtures

## Safety Data Sheet

according to WHMIS

### VITA VM LC Chroma Plus EFFECT LINER

Revision date: 12.08.2019

Product code: 217

Page 2 of 7

#### Hazardous components

CAS No	Chemical name	Quantity
72869-86-4	7,7,9-Trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)	15 - < 40% (*)
109-16-0	2,2'-ethylenedioxydiethyl dimethacrylate	3 - < 7% (*)
2867-47-2	2-dimethylaminoethyl methacrylate	0.1 - < 1% (*)

(\*) The actual concentration is withheld as a trade secret.

#### 4. First-aid measures

##### Description of first aid measures

###### **After inhalation**

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

###### **After contact with skin**

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

###### **After contact with eyes**

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

###### **After ingestion**

Rinse mouth immediately and drink plenty of water.

##### Most important symptoms and effects, whether acute or delayed

No information available.

##### Indication of immediate medical attention and special treatment needed

Treat symptomatically.

#### 5. Fire-fighting measures

##### Extinguishing media

###### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

##### Specific hazards arising from the hazardous product

Non-flammable.

##### Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

##### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately.

Do not allow entering drains or surface water.

#### 6. Accidental release measures

##### Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

##### Environmental precautions

Do not allow to enter into surface water or drains.

##### Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

## Safety Data Sheet

according to WHMIS

### VITA VM LC Chroma Plus EFFECT LINER

Revision date: 12.08.2019

Product code: 217

Page 3 of 7

#### Reference to other sections

- Safe handling: see section 7
- Personal protection equipment: see section 8
- Disposal: see section 13

## 7. Handling and storage

### Precautions for safe handling

#### **Advice on safe handling**

No special measures are necessary.

#### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

### Conditions for safe storage, including any incompatibilities

#### **Requirements for storage rooms and vessels**

Keep container tightly closed.

#### **Hints on joint storage**

No special measures are necessary.

## 8. Exposure controls/Personal protection

### Control parameters

### Exposure controls



#### **Protective and hygiene measures**

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

#### **Eye/face protection**

Wear eye protection/face protection.

#### **Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatril P Breakthrough time (maximum wearing time) 10 min NBR (Nitrile rubber)

#### **Skin protection**

Use of protective clothing.

#### **Respiratory protection**

Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state:	Liquid
Colour:	
Odour:	characteristic
pH-Value:	not determined

## Safety Data Sheet

according to WHMIS

### VITA VM LC Chroma Plus EFFECT LINER

Revision date: 12.08.2019

Product code: 217

Page 4 of 7

#### Changes in the physical state

Melting point:	not determined
Initial boiling point and boiling range:	139 °C
Flash point:	251 °C

#### Flammability

Solid:	not applicable
Gas:	not applicable

#### Explosive properties

The product is not: Explosive.

Lower explosive limits:	not determined
Upper explosive limits:	not determined

#### Auto-ignition temperature

Solid:	not applicable
Gas:	not applicable

Decomposition temperature:	not determined
----------------------------	----------------

#### Oxidizing properties

Not oxidising.

Vapour pressure: (at 50 °C)	<=1100 hPa
--------------------------------	------------

Density:	not determined
----------	----------------

Water solubility:	No
-------------------	----

#### Solubility in other solvents

not determined

Partition coefficient:	not determined
------------------------	----------------

Vapour density:	not determined
-----------------	----------------

Evaporation rate:	not determined
-------------------	----------------

#### Other information

Solid content:	70,0 %
----------------	--------

### 10. Stability and reactivity

#### Reactivity

No hazardous reaction when handled and stored according to provisions.

#### Chemical stability

The product is stable under storage at normal ambient temperatures.

#### Possibility of hazardous reactions

No known hazardous reactions.

#### Conditions to avoid

none

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

### 11. Toxicological information

#### Information on toxicological effects

## Safety Data Sheet

according to WHMIS

### VITA VM LC Chroma Plus EFFECT LINER

Revision date: 12.08.2019

Product code: 217

Page 5 of 7

#### Acute toxicity

Based on available data, the classification criteria are not met.

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE 500 mg/kg			
	dermal	ATE 1100 mg/kg			

#### Irritation and corrosivity

Based on available data, the classification criteria are not met.

#### Sensitizing effects

May cause an allergic skin reaction.

(7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers); 2,2'-ethylenedioxydiethyl dimethacrylate; 2-dimethylaminoethyl methacrylate)

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### STOT-single exposure

Based on available data, the classification criteria are not met.

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### Additional information on tests

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

Harmful to aquatic life with long lasting effects.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

#### Mobility in soil

The product has not been tested.

#### Other adverse effects

No information available.

#### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

## 13. Disposal considerations

#### Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Handle contaminated packages in the same way as the substance itself. Waste codes/waste designations according to EWC/AVV

## Safety Data Sheet

according to WHMIS

### VITA VM LC Chroma Plus EFFECT LINER

Revision date: 12.08.2019

Product code: 217

Page 6 of 7

#### 14. Transport information

##### Marine transport (IMDG)

<u>UN number:</u>	No dangerous good in sense of this transport regulation.
<u>United Nations proper shipping name:</u>	No dangerous good in sense of this transport regulation.
<u>Transport hazard class(es):</u>	No dangerous good in sense of this transport regulation.
<u>Packing group:</u>	No dangerous good in sense of this transport regulation.

##### Air transport (ICAO-TI/IATA-DGR)

<u>UN number:</u>	No dangerous good in sense of this transport regulation.
<u>United Nations proper shipping name:</u>	No dangerous good in sense of this transport regulation.
<u>Transport hazard class(es):</u>	No dangerous good in sense of this transport regulation.
<u>Packing group:</u>	No dangerous good in sense of this transport regulation.

##### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

#### 15. Regulatory information

##### Canadian regulations

#### 16. Other information

##### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%  
CLP: Classification, labelling and Packaging  
REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways

## Safety Data Sheet

according to WHMIS

### VITA VM LC Chroma Plus EFFECT LINER

Revision date: 12.08.2019

Product code: 217

Page 7 of 7

(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

VOC: Volatile Organic Compounds

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*