

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 12.08.2019

Product code: 218

Page 1 of 8

## 1. Identification

### Product identifier

VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

### Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name:	VITA Zahnfabrik H.Rauter GmbH & Co.KG	
Street:	Spitalgasse 3	
Place:	D-79713 Bad Säckingen	
Post-office box:	1338	
	D-79704 Bad Säckingen	
Telephone:	+49(0)7761-562-0	Telefax: +49(0)7761-562-299
e-mail:	info@vita-zahnfabrik.com	
Contact person:	regulatory affairs	
e-mail:	info@vita-zahnfabrik.com	
Internet:	www.vita-zahnfabrik.com	
Responsible Department:	Regulatory Affairs	

**Emergency telephone number:** +49-(0)761-19240

### Further Information

medical device

## 2. Hazard identification

### Classification of the substance or mixture

#### UN GHS (ST/SG/AC.10/11/Rev.8)

Hazard categories:

Skin corrosion/irritation: Skin Irrit. 2

Serious eye damage/eye irritation: Eye Irrit. 2

Respiratory or skin sensitization: Skin Sens. 1

Reproductive toxicity: Repr. 2

Hazardous to the aquatic environment: Aquatic Chronic 3

Hazard Statements:

Causes skin irritation.

Causes serious eye irritation.

May cause an allergic skin reaction.

Suspected of damaging fertility or the unborn child.

Harmful to aquatic life with long lasting effects.

### Label elements

#### UN GHS (ST/SG/AC.10/11/Rev.8)

#### Hazard components for labelling

7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)

2-dimethylaminoethyl methacrylate

Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

**Signal word:** Warning

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)


### VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 12.08.2019

Product code: 218

Page 2 of 8

#### Pictograms:



#### Hazard statements

H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H361	Suspected of damaging fertility or the unborn child.
H412	Harmful to aquatic life with long lasting effects.

#### Precautionary statements

P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P273	Avoid release to the environment.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P302+P352	IF ON SKIN: Wash with plenty of water.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

#### Other hazards

No information available.

### 3. Composition/information on ingredients

#### Mixtures

##### Hazardous components

CAS No	Chemical name	Quantity
	Classification (UN GHS (ST/SG/AC.10/11/Rev.8))	
72869-86-4	7,7,9-Trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)	15 - < 20 %
	Skin Sens. 1B, Aquatic Chronic 3; H317 H412	
94108-97-1	Ditrimethylolpropane Tetraacrylate	5 - < 10 %
	Eye Irrit. 2, Aquatic Chronic 2; H319 H411	
2867-47-2	2-dimethylaminoethyl methacrylate	5 - < 10 %
	Acute Tox. 4, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1; H312 H302 H315 H319 H317	
10373-78-1	Camphorquinone	1 - < 5 %
	Skin Irrit. 2, Eye Irrit. 2, STOT SE 3; H315 H319 H335	
75980-60-8	Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	1 - < 5 %
	Repr. 2, Skin Sens. 1, Aquatic Chronic 2; H361f H317 H411	

### 4. First-aid measures

#### Description of first aid measures

##### After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

##### After contact with skin

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

##### After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 12.08.2019

Product code: 218

Page 3 of 8

ophthalmologist immediately.

#### **After ingestion**

Rinse mouth immediately and drink plenty of water.

#### **Most important symptoms and effects, whether acute or delayed**

No information available.

#### **Indication of immediate medical attention and special treatment needed**

Treat symptomatically.

## 5. Fire-fighting measures

### **Extinguishing media**

#### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

### **Specific hazards arising from the hazardous product**

Non-flammable.

### **Special protective equipment and precautions for fire-fighters**

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

### **Additional information**

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

## 6. Accidental release measures

### **Personal precautions, protective equipment and emergency procedures**

#### **General advice**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

### **Environmental precautions**

Do not allow to enter into surface water or drains.

### **Methods and material for containment and cleaning up**

#### **Other information**

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

### **Reference to other sections**

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## 7. Handling and storage

### **Precautions for safe handling**

#### **Advice on safe handling**

No special measures are necessary.

#### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

#### **Advice on general occupational hygiene**

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

### **Conditions for safe storage, including any incompatibilities**

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 12.08.2019

Product code: 218

Page 4 of 8

#### Requirements for storage rooms and vessels

Keep container tightly closed.

#### Hints on joint storage

No special measures are necessary.

### 8. Exposure controls/Personal protection

#### Control parameters

#### Exposure controls



#### Individual protection measures, such as personal protective equipment

##### Eye/face protection

Suitable eye protection: goggles.

##### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL DermatrilP NBR (Nitrile rubber) Breakthrough time: 480 min

##### Skin protection

Use of protective clothing.

##### Respiratory protection

Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

### 9. Physical and chemical properties

#### Information on basic physical and chemical properties

Physical state:

Colour:

##### Changes in the physical state

Melting point/freezing point: not determined

Boiling point or initial boiling point and boiling range: ?

Flash point: 151 °C

##### Flammability

Solid/liquid: not determined

Gas: not applicable

##### Explosive properties

The product is not: Explosive.

Lower explosive limits: not determined

Upper explosive limits: not determined

##### Self-ignition temperature

Solid: not determined

Gas: not applicable

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 12.08.2019

Product code: 218

Page 5 of 8

Decomposition temperature:	not determined
pH-Value:	not determined
Water solubility:	No
<b>Solubility in other solvents</b>	
not determined	
Partition coefficient n-octanol/water:	not determined
Vapour pressure: (at 50 °C)	<=1100 hPa
Density:	not determined
Relative vapour density:	not determined

#### Other information

##### Information with regard to physical hazard classes

Oxidizing properties  
Not oxidising.

##### Other safety characteristics

Solid content:	49,4 %
Evaporation rate:	not determined

## 10. Stability and reactivity

#### Reactivity

No hazardous reaction when handled and stored according to provisions.

#### Chemical stability

The product is stable under storage at normal ambient temperatures.

#### Possibility of hazardous reactions

No known hazardous reactions.

#### Conditions to avoid

none

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

#### Information on toxicological effects

##### Acute toxicity

Based on available data, the classification criteria are not met.

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE 500 mg/kg			
	dermal	ATE 1100 mg/kg			

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 12.08.2019

Product code: 218

Page 6 of 8

#### **Irritation and corrosivity**

Causes skin irritation.  
Causes serious eye irritation.

#### **Sensitizing effects**

May cause an allergic skin reaction.  
(7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers);  
2-dimethylaminoethyl methacrylate; Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide)

#### **Carcinogenic/mutagenic/toxic effects for reproduction**

Suspected of damaging fertility or the unborn child. (Diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide)  
Germ cell mutagenicity: Based on available data, the classification criteria are not met.  
Carcinogenicity: Based on available data, the classification criteria are not met.

#### **STOT-single exposure**

Based on available data, the classification criteria are not met.

#### **STOT-repeated exposure**

Based on available data, the classification criteria are not met.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### **Additional information on tests**

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### **Ecotoxicity**

Harmful to aquatic life with long lasting effects.

#### **Persistence and degradability**

The product has not been tested.

#### **Bioaccumulative potential**

The product has not been tested.

#### **Mobility in soil**

The product has not been tested.

#### **Endocrine disrupting properties**

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### **Other adverse effects**

No information available.

#### **Further information**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

## 13. Disposal considerations

#### **Waste treatment methods**

##### **Disposal recommendations**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### **Contaminated packaging**

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself. Waste codes/waste designations according to EWC/AVV

## 14. Transport information

#### **Marine transport (IMDG)**

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 12.08.2019

Product code: 218

Page 7 of 8

**UN number or ID number:** UN 3243  
**United Nations proper shipping name:** SOLIDS CONTAINING TOXIC LIQUID, N.O.S.  
**Transport hazard class(es):** 6.1  
**Packing group:** II  
 Hazard label: 6.1



Special Provisions: 217, 274  
 Limited quantity: 500 g  
 Excepted quantity: E4  
 EmS: F-A, S-A

#### Air transport (ICAO-TI/IATA-DGR)

**UN number or ID number:** UN 3243  
**United Nations proper shipping name:** SOLIDS CONTAINING TOXIC LIQUID, N.O.S.  
**Transport hazard class(es):** 6.1  
**Packing group:** II  
 Hazard label: 6.1



Special Provisions: A50  
 Limited quantity Passenger: 1 kg  
 Passenger LQ: Y644  
 Excepted quantity: E4  
 IATA-packing instructions - Passenger: 669  
 IATA-max. quantity - Passenger: 25 kg  
 IATA-packing instructions - Cargo: 676  
 IATA-max. quantity - Cargo: 100 kg

#### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

#### Special precautions for user

No information available.

#### Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

not applicable

### 15. Regulatory information

#### National regulatory information

Employment restrictions: Observe employment restrictions for young people.  
 Skin resorption/Sensitization: Causes allergic hypersensitivity reactions.

### 16. Other information

#### Changes

This data sheet contains changes from the previous version in section(s): 1,2,4,6,7,8,9,11,12,13,14,15,16.

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITA VM LC OPAQUE PASTE und VM LC GINGIVA OPAQUE PASTE

Revision date: 12.08.2019

Product code: 218

Page 8 of 8

(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LL50: Lethal loading, 50%

EL50: Effect loading, 50%

EC50: Effective Concentration 50%

ErC50: Effective Concentration 50%, growth rate

NOEC: No Observed Effect Concentration

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative

RID: Regulations concerning the international carriage of dangerous goods by rail

ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

VOC: Volatile Organic Compounds

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*