

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC WINDOW

Revision date: 13.08.2019

Product code: 264

Page 1 of 7

## 1. Identification

### Product identifier

VITA VM LC WINDOW

### Recommended use of the chemical and restrictions on use

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H.Rauter GmbH &amp; Co.KG

Post-office box: 1338  
79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

**Emergency phone number:** +49-(0)761-19240

### Further Information

medical device

## 2. Hazard(s) identification

### Classification of the chemical

#### 29 CFR Part 1910.1200

Serious eye damage/eye irritation: Eye Irrit. 2A

Respiratory or skin sensitization: Skin Sens. 1

### Label elements

#### 29 CFR Part 1910.1200

Signal word: Warning

Pictograms:



### Hazard statements

May cause an allergic skin reaction

Causes serious eye irritation

### Precautionary statements

Avoid breathing dust/fume/gas/mist/vapors/spray.

Wear protective gloves/protective clothing/eye protection/face protection.

### Hazards not otherwise classified

No information available.

## 3. Composition/information on ingredients

### Mixtures

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC WINDOW

Revision date: 13.08.2019

Product code: 264

Page 2 of 7

#### Hazardous components

CAS No	Components	Quantity
72869-86-4	7,7,9-Trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)	32.87 %
94108-97-1	Ditrimethylolpropane Tetraacrylate	11.5 %
2867-47-2	2-dimethylaminoethyl methacrylate	0.1188 %

#### 4. First-aid measures

##### Description of first aid measures

###### **After inhalation**

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

###### **After contact with skin**

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary. After contact with skin, wash immediately with plenty of water and soap.

###### **After contact with eyes**

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

###### **After ingestion**

Rinse mouth immediately and drink plenty of water.

##### Most important symptoms and effects, both acute and delayed

No information available.

##### Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

#### 5. Fire-fighting measures

##### Extinguishing media

###### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

##### Specific hazards arising from the chemical

Non-flammable.

##### Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

##### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

#### 6. Accidental release measures

##### Personal precautions, protective equipment and emergency procedures

Provide adequate ventilation. Do not breathe gas/fume/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

##### Environmental precautions

Do not allow to enter into surface water or drains.

##### Methods and material for containment and cleaning up

Take up mechanically. Treat the recovered material as prescribed in the section on waste disposal.

##### Reference to other sections

Safe handling: see section 7

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC WINDOW

Revision date: 13.08.2019

Product code: 264

Page 3 of 7

Personal protection equipment (PPE): see section 8

Disposal: see section 13

## 7. Handling and storage

### Precautions for safe handling

#### Advice on safe handling

No special measures are necessary.

#### Advice on protection against fire and explosion

No special fire protection measures are necessary.

### Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Keep container tightly closed.

#### Hints on joint storage

No special measures are necessary.

## 8. Exposure controls/personal protection

### Control parameters

### Exposure controls



#### Protective and hygiene measures

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.

Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

#### Eye/face protection

Suitable eye protection: goggles.

#### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles Dermatril P Breakthrough time (maximum wearing time) 30 min NBR (Nitrile rubber)

#### Skin protection

Use of protective clothing.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state:

Color:

Odor: characteristic

pH-Value: not determined

#### Changes in the physical state

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC WINDOW

Revision date: 13.08.2019

Product code: 264

Page 4 of 7

Melting point/freezing point: not determined  
 Initial boiling point and boiling range: 251 °C  
 Flash point: > 250 °C

#### Flammability

Solid: not determined  
 Gas: not applicable

#### Explosive properties

The product is not: Explosive.

Lower explosion limits: not determined  
 Upper explosion limits: not determined

#### Auto-ignition temperature

Solid: not determined  
 Gas: not applicable

Decomposition temperature: not determined

#### Oxidizing properties

Not oxidising.

Vapor pressure: <=1100 hPa  
 (at 50 °C)

Density: not determined

Water solubility: No

#### Solubility in other solvents

not determined

Partition coefficient: not determined

Vapor density: not determined

Evaporation rate: not determined

#### Other information

Solid content: 98,5 %

## 10. Stability and reactivity

#### Reactivity

No hazardous reaction when handled and stored according to provisions.

#### Chemical stability

The product is stable under storage at normal ambient temperatures.

#### Possibility of hazardous reactions

No known hazardous reactions.

#### Conditions to avoid

none/none

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

#### Information on toxicological effects

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC WINDOW

Revision date: 13.08.2019

Product code: 264

Page 5 of 7

#### Acute toxicity

Based on available data, the classification criteria are not met.

CAS No	Components				
	Exposure route	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE 500 mg/kg			
	dermal	ATE 1100 mg/kg			

#### Irritation and corrosivity

Causes serious eye irritation

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

#### Sensitizing effects

May cause an allergic skin reaction

(7,7,9-Trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers); 2-dimethylaminoethyl methacrylate)

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### Specific target organ toxicity (STOT) - single exposure

Based on available data, the classification criteria are not met.

#### Specific target organ toxicity (STOT) - repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### Additional information on tests

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

Harmful to aquatic life with long lasting effects.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

#### Mobility in soil

The product has not been tested.

#### Other adverse effects

No information available.

#### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

## 13. Disposal considerations

#### Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Handle contaminated packages in the same way as the substance itself. Waste codes/waste designations

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC WINDOW

Revision date: 13.08.2019

Product code: 264

Page 6 of 7

according to EWC/AVV

#### 14. Transport information

##### Marine transport (IMDG)

**UN number:** No dangerous good in sense of this transport regulation.  
**UN proper shipping name:** No dangerous good in sense of this transport regulation.  
**Transport hazard class(es):** No dangerous good in sense of this transport regulation.  
**Packing group:** No dangerous good in sense of this transport regulation.

##### Air transport (ICAO-TI/IATA-DGR)

**UN number:** No dangerous good in sense of this transport regulation.  
**UN proper shipping name:** No dangerous good in sense of this transport regulation.  
**Transport hazard class(es):** No dangerous good in sense of this transport regulation.  
**Packing group:** No dangerous good in sense of this transport regulation.

##### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: no

##### Special precautions for user

No information available.

##### Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

not applicable

#### 15. Regulatory information

##### U.S. Regulations

###### National regulatory information

SARA Section 311/312 Hazards:

7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)  
 (72869-86-4): Immediate (acute) health hazard  
 Ditrimethylolpropane Tetraacrylate (94108-97-1): Immediate (acute) health hazard  
 2-dimethylaminoethyl methacrylate (2867-47-2): Immediate (acute) health hazard

##### State Regulations

###### Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)

This product can not expose you to chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

#### 16. Other information

Revision date: 13.08.2019

Revision No: 2

##### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
 (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
 IMDG: International Maritime Code for Dangerous Goods  
 IATA: International Air Transport Association  
 GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
 EINECS: European Inventory of Existing Commercial Chemical Substances  
 ELINCS: European List of Notified Chemical Substances  
 CAS: Chemical Abstracts Service  
 LC50: Lethal concentration, 50%  
 LD50: Lethal dose, 50%  
 CLP: Classification, labelling and Packaging

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA VM LC WINDOW

Revision date: 13.08.2019

Product code: 264

Page 7 of 7

REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
EmS: Emergency Schedules  
MFAG: Medical First Aid Guide  
ICAO: International Civil Aviation Organization  
MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
IBC: Intermediate Bulk Container  
VOC: Volatile Organic Compounds  
SVHC: Substance of Very High Concern  
For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Other data

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*