

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITAVM LC PAINT

Revision date: 06.08.2019

Product code: 162

Page 1 of 8

## 1. Identification

### Product identifier

VITAVM LC PAINT

### Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H.Rauter GmbH & Co.KG  
Street: Spitalgasse 3  
Place: D-79713 Bad Säckingen  
Post-office box: 1338  
D-79704 Bad Säckingen  
Telephone: +49(0)7761-562-0  
e-mail: info@vita-zahnfabrik.com  
Contact person: regulatory affairs  
e-mail: info@vita-zahnfabrik.com  
Internet: www.vita-zahnfabrik.com  
Responsible Department: Regulatory Affairs

Telefax: +49(0)7761-562-299

**Emergency telephone number:** +49-(0)761-19240

### Further Information

medical device

## 2. Hazard identification

### Classification of the substance or mixture

#### UN GHS (ST/SG/AC.10/11/Rev.8)

Hazard categories:

Respiratory or skin sensitization: Skin Sens. 1

Hazardous to the aquatic environment: Aquatic Chronic 3

Hazard Statements:

May cause an allergic skin reaction.

Harmful to aquatic life with long lasting effects.

### Label elements

#### UN GHS (ST/SG/AC.10/11/Rev.8)

#### Hazard components for labelling

2,2'-ethylenedioxydiethyl dimethacrylate

7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)

**Signal word:** Warning

#### Pictograms:



#### Hazard statements

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

#### Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITAVM LC PAINT

Revision date: 06.08.2019

Product code: 162

Page 2 of 8

P273 Avoid release to the environment.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection.  
 P302+P352 IF ON SKIN: Wash with plenty of water.

#### Other hazards

No information available.

### 3. Composition/information on ingredients

#### Mixtures

##### Hazardous components

CAS No	Chemical name	Quantity
	Classification (UN GHS (ST/SG/AC.10/11/Rev.8))	
72869-86-4	7,7,9-Trimethyl-4,13-dioxo-3,14-dioxo-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers)	40 - < 45 %
	Skin Sens. 1B, Aquatic Chronic 3; H317 H412	
109-16-0	2,2'-ethylenedioxydiethyl dimethacrylate	15 - < 20 %
	Skin Sens. 1B; H317	
2867-47-2	2-dimethylaminoethyl methacrylate	< 1 %
	Acute Tox. 4, Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Skin Sens. 1; H312 H302 H315 H319 H317	

### 4. First-aid measures

#### Description of first aid measures

##### **After inhalation**

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

##### **After contact with skin**

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary. After contact with skin, wash immediately with plenty of water and soap.

##### **After contact with eyes**

Rinse immediately carefully and thoroughly with eye-bath or water. In case of eye irritation consult an ophthalmologist.

##### **After ingestion**

Rinse mouth immediately and drink plenty of water.

#### Most important symptoms and effects, whether acute or delayed

No information available.

#### Indication of immediate medical attention and special treatment needed

Treat symptomatically.

### 5. Fire-fighting measures

#### Extinguishing media

##### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

#### Specific hazards arising from the hazardous product

Non-flammable.

#### Special protective equipment and precautions for fire-fighters

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITAVM LC PAINT

Revision date: 06.08.2019

Product code: 162

Page 3 of 8

#### Additional information

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

#### 6. Accidental release measures

##### Personal precautions, protective equipment and emergency procedures

###### General advice

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

###### Environmental precautions

Do not allow to enter into surface water or drains.

###### Methods and material for containment and cleaning up

###### Other information

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

###### Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

#### 7. Handling and storage

##### Precautions for safe handling

###### Advice on safe handling

No special measures are necessary.

###### Advice on protection against fire and explosion

No special fire protection measures are necessary.

###### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat, drink, smoke, sniff.

##### Conditions for safe storage, including any incompatibilities

###### Requirements for storage rooms and vessels

Keep container tightly closed.

###### Hints on joint storage

No special measures are necessary.

#### 8. Exposure controls/Personal protection

##### Control parameters

##### Exposure controls



##### Individual protection measures, such as personal protective equipment

###### Eye/face protection

Wear eye/face protection.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITAVM LC PAINT

Revision date: 06.08.2019

Product code: 162

Page 4 of 8

#### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Camatril Velour Breakthrough time: 30 min NBR (Nitrile rubber)

#### Skin protection

Use of protective clothing.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state:

Colour:

Odour: characteristic

#### Changes in the physical state

Melting point/freezing point: not determined

Boiling point or initial boiling point and boiling range: X

Flash point: X

#### Flammability

Solid/liquid: not applicable

Gas: not applicable

#### Explosive properties

The product is not: Explosive.

Lower explosive limits: not determined

Upper explosive limits: not determined

#### Self-ignition temperature

Solid: not applicable

Gas: not applicable

Decomposition temperature: not determined

pH-Value: not determined

Water solubility: No

#### Solubility in other solvents

not determined

Partition coefficient n-octanol/water: not determined

Vapour pressure: <=1100 hPa  
(at 50 °C)

Density: not determined

Relative vapour density: not determined

### Other information

#### Information with regard to physical hazard classes

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITAVM LC PAINT

Revision date: 06.08.2019

Product code: 162

Page 5 of 8

Oxidizing properties

Not oxidising.

#### Other safety characteristics

Solid content:

1,10 %

Evaporation rate:

not determined

## 10. Stability and reactivity

### Reactivity

No hazardous reaction when handled and stored according to provisions.

### Chemical stability

The product is stable under storage at normal ambient temperatures.

### Possibility of hazardous reactions

No known hazardous reactions.

### Conditions to avoid

none

### Incompatible materials

No information available.

### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

### Information on toxicological effects

#### Acute toxicity

Based on available data, the classification criteria are not met.

CAS No	Chemical name				
	Route of exposure	Dose	Species	Source	Method
2867-47-2	2-dimethylaminoethyl methacrylate				
	oral	ATE 500 mg/kg			
	dermal	ATE 1100 mg/kg			

#### Irritation and corrosivity

Based on available data, the classification criteria are not met.

#### Sensitizing effects

May cause an allergic skin reaction.

(7,7,9-Trimethyl-4,13-dioxo-3,14-dioxa-5,12-diaza-hexadecan-1,16-diol-dimethacrylat (mixture of isomers); 2,2'-ethylenedioxydiethyl dimethacrylate; 2-dimethylaminoethyl methacrylate)

#### Carcinogenic/mutagenic/toxic effects for reproduction

Based on available data, the classification criteria are not met.

#### STOT-single exposure

Based on available data, the classification criteria are not met.

#### STOT-repeated exposure

Based on available data, the classification criteria are not met.

#### Aspiration hazard

Based on available data, the classification criteria are not met.

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITAVM LC PAINT

Revision date: 06.08.2019

Product code: 162

Page 6 of 8

#### Additional information on tests

The mixture is classified as hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

Harmful to aquatic life with long lasting effects.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

#### Mobility in soil

The product has not been tested.

#### Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### Other adverse effects

No information available.

#### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

## 13. Disposal considerations

#### Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Handle contaminated packages in the same way as the substance itself.

## 14. Transport information

#### Marine transport (IMDG)

##### UN number or ID number:

No dangerous good in sense of this transport regulation.

##### United Nations proper shipping name:

No dangerous good in sense of this transport regulation.

##### Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### Packing group:

No dangerous good in sense of this transport regulation.

#### Air transport (ICAO-TI/IATA-DGR)

##### UN number or ID number:

No dangerous good in sense of this transport regulation.

##### United Nations proper shipping name:

No dangerous good in sense of this transport regulation.

##### Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### Packing group:

No dangerous good in sense of this transport regulation.

#### Environmental hazards

ENVIRONMENTALLY HAZARDOUS: No

#### Special precautions for user

No information available.

#### Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

not applicable

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITAVM LC PAINT

Revision date: 06.08.2019

Product code: 162

Page 7 of 8

## 15. Regulatory information

### National regulatory information

Employment restrictions: Observe employment restrictions for young people.  
Skin resorption/Sensitization: Causes allergic hypersensitivity reactions.

## 16. Other information

### Changes

This data sheet contains changes from the previous version in section(s): 1,2,4,6,7,8,9,11,12,13,14,15,16.

### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%  
CLP: Classification, labelling and Packaging  
REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
EmS: Emergency Schedules  
MFAG: Medical First Aid Guide  
ICAO: International Civil Aviation Organization  
MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
IBC: Intermediate Bulk Container  
VOC: Volatile Organic Compounds  
SVHC: Substance of Very High Concern  
For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

### Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for

## Safety Data Sheet

according to UN GHS (ST/SG/AC.10/11/Rev.8)

### VITAVM LC PAINT

Revision date: 06.08.2019

Product code: 162

Page 8 of 8

adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*