

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA AKZENT Plus PASTE FLUID

Revision date: 13.08.2019

Product code: 277

Page 1 of 6

## 1. Identification

### Product identifier

VITA AKZENT Plus PASTE FLUID

### Recommended use of the chemical and restrictions on use

#### Use of the substance/mixture

Use as laboratory reagent

### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H. Rauter GmbH & Co. KG

Post-office box: 1338  
79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

**Emergency phone number:** +49-(0)7761-562-0

### Further Information

medical device

## 2. Hazard(s) identification

### Classification of the chemical

#### 29 CFR Part 1910.1200

This mixture is not classified as hazardous in accordance with Regulation 29 CFR 1910.1200(d).

### Label elements

### Hazards not otherwise classified

No information available.

## 3. Composition/information on ingredients

### Mixtures

#### Chemical characterization

Substance, organic

#### Hazardous components

CAS No	Components	Quantity
55965-84-9	reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1)	0.0009 %

## 4. First-aid measures

### Description of first aid measures

#### After inhalation

Provide fresh air. If breathing is irregular or stopped, administer artificial respiration. Medical treatment necessary.

#### After contact with skin

After contact with skin, wash immediately with plenty of water and soap. Take off immediately all contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention.

#### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

#### After ingestion

Rinse mouth immediately and drink plenty of water.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA AKZENT Plus PASTE FLUID

Revision date: 13.08.2019

Product code: 277

Page 2 of 6

#### Most important symptoms and effects, both acute and delayed

No information available.

#### Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### 5. Fire-fighting measures

#### Extinguishing media

##### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

#### Specific hazards arising from the chemical

Non-flammable.

#### Special protective equipment and precautions for fire-fighters

In case of fire: Wear self-contained breathing apparatus.

#### **Additional information**

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

### 6. Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

Use personal protection equipment.

#### Environmental precautions

Do not allow to enter into surface water or drains.

#### Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment (PPE): see section 8

Disposal: see section 13

### 7. Handling and storage

#### Precautions for safe handling

##### **Advice on safe handling**

No special measures are necessary.

##### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

#### Conditions for safe storage, including any incompatibilities

##### **Requirements for storage rooms and vessels**

Keep container tightly closed.

##### **Hints on joint storage**

No special measures are necessary.

### 8. Exposure controls/personal protection

#### Control parameters

#### **Exposure limits**

CAS No.	Substance	ppm	mg/m <sup>3</sup>	f/cc	Category	Origin
56-81-5	Glycerin (mist) Respirable fraction	-	5		TWA (8 h)	PEL

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA AKZENT Plus PASTE FLUID

Revision date: 13.08.2019

Product code: 277

Page 3 of 6

#### Exposure controls



##### **Protective and hygiene measures**

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat, drink, smoke, sniff.

##### **Eye/face protection**

Wear eye protection/face protection.

##### **Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Suitable material: NBR (Nitrile rubber)

##### **Skin protection**

Use of protective clothing.

##### **Respiratory protection**

Technical ventilation of workplace Open windows to ensure natural ventilation.

## 9. Physical and chemical properties

#### Information on basic physical and chemical properties

Physical state:	Liquid
Color:	colorless
Odor:	characteristic
pH-Value:	not determined

##### **Changes in the physical state**

Melting point/freezing point:	not determined
Initial boiling point and boiling range:	207 °C
Flash point:	109 °C

##### **Flammability**

Solid:	not applicable
Gas:	not applicable

##### **Explosive properties**

The product is not: Explosive.

Lower explosion limits:	not determined
Upper explosion limits:	not determined

##### **Auto-ignition temperature**

Solid:	not applicable
Gas:	not applicable
Decomposition temperature:	not determined

##### **Oxidizing properties**

Not oxidising.

Vapor pressure:	not determined
Density:	1,06200 g/cm <sup>3</sup>

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA AKZENT Plus PASTE FLUID

Revision date: 13.08.2019

Product code: 277

Page 4 of 6

#### Solubility in other solvents

not determined

Partition coefficient:

not determined

Vapor density:

not determined

Evaporation rate:

not determined

#### Other information

Solid content:

0,0 %

## 10. Stability and reactivity

#### Reactivity

No hazardous reaction when handled and stored according to provisions.

#### Chemical stability

The product is stable under storage at normal ambient temperatures.

#### Possibility of hazardous reactions

No known hazardous reactions.

#### Conditions to avoid

none/none

#### Incompatible materials

No information available.

#### Hazardous decomposition products

No known hazardous decomposition products.

## 11. Toxicological information

#### Information on toxicological effects

##### Acute toxicity

CAS No	Components				
	Exposure route	Dose	Species	Source	Method
55965-84-9	reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1)				
	oral	ATE 100 mg/kg			
	dermal	ATE 50 mg/kg			
	inhalation vapour	ATE 0,5 mg/l			
	inhalation aerosol	ATE 0,05 mg/l			

#### Additional information on tests

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

#### Ecotoxicity

The product is not: Ecotoxic.

#### Persistence and degradability

The product has not been tested.

#### Bioaccumulative potential

The product has not been tested.

#### Mobility in soil

The product has not been tested.

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA AKZENT Plus PASTE FLUID

Revision date: 13.08.2019

Product code: 277

Page 5 of 6

#### Other adverse effects

No information available.

#### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

### 13. Disposal considerations

#### Waste treatment methods

##### **Disposal recommendations**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### **Contaminated packaging**

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

### 14. Transport information

#### **Marine transport (IMDG)**

<b><u>UN number:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>UN proper shipping name:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>Transport hazard class(es):</u></b>	No dangerous good in sense of this transport regulation.
<b><u>Packing group:</u></b>	No dangerous good in sense of this transport regulation.

#### **Air transport (ICAO-TI/IATA-DGR)**

<b><u>UN number:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>UN proper shipping name:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>Transport hazard class(es):</u></b>	No dangerous good in sense of this transport regulation.
<b><u>Packing group:</u></b>	No dangerous good in sense of this transport regulation.

#### **Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: no

#### **Special precautions for user**

No information available.

#### **Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code**

not applicable

### 15. Regulatory information

#### **U.S. Regulations**

##### **National regulatory information**

SARA Section 311/312 Hazards:  
 reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (55965-84-9): Immediate (acute) health hazard

##### **State Regulations**

##### **Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)**

This product can not expose you to chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

### 16. Other information

Revision date: 13.08.2019  
 Revision No: 2

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA AKZENT Plus PASTE FLUID

Revision date: 13.08.2019

Product code: 277

Page 6 of 6

#### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%  
CLP: Classification, labelling and Packaging  
REACH: Registration, Evaluation and Authorization of Chemicals  
GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
UN: United Nations  
DNEL: Derived No Effect Level  
DMEL: Derived Minimal Effect Level  
PNEC: Predicted No Effect Concentration  
ATE: Acute toxicity estimate  
LL50: Lethal loading, 50%  
EL50: Effect loading, 50%  
EC50: Effective Concentration 50%  
ErC50: Effective Concentration 50%, growth rate  
NOEC: No Observed Effect Concentration  
BCF: Bio-concentration factor  
PBT: persistent, bioaccumulative, toxic  
vPvB: very persistent, very bioaccumulative  
RID: Regulations concerning the international carriage of dangerous goods by rail  
ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
EmS: Emergency Schedules  
MFAG: Medical First Aid Guide  
ICAO: International Civil Aviation Organization  
MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
IBC: Intermediate Bulk Container  
VOC: Volatile Organic Compounds  
SVHC: Substance of Very High Concern  
For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Other data

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*