

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA LOW FUSING MODELLING LIQUID

Revision date: 05.08.2019

Product code: 160

Page 1 of 6

#### 1. Identification

##### Product identifier

VITA LOW FUSING MODELLING LIQUID

##### Recommended use of the chemical and restrictions on use

###### Use of the substance/mixture

Use as laboratory reagent

##### Details of the supplier of the safety data sheet

Company name: VITA Zahnfabrik H. Rauter GmbH &amp; Co. KG

Post-office box: 1338  
79704 Bad Säckingen

Telephone: +49(0)7761-562-0

Telefax: +49(0)7761-562-299

e-mail: info@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

**Emergency phone number:** +49-(0)7761-562-0

##### Further Information

medical device

#### 2. Hazard(s) identification

##### Classification of the chemical

###### 29 CFR Part 1910.1200

This mixture is not classified as hazardous in accordance with Regulation 29 CFR 1910.1200(d).

##### Label elements

##### Hazards not otherwise classified

No information available.

#### 3. Composition/information on ingredients

##### Mixtures

###### Chemical characterization

Mixtures Product/Substance is inorganic. Substance, organic

#### 4. First-aid measures

##### Description of first aid measures

###### After inhalation

Provide fresh air.

###### After contact with skin

Wash with plenty of water. Take off contaminated clothing and wash it before reuse.

###### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water.

###### After ingestion

Rinse mouth immediately and drink plenty of water.

##### Most important symptoms and effects, both acute and delayed

No information available.

##### Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

#### 5. Fire-fighting measures

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA LOW FUSING MODELLING LIQUID

Revision date: 05.08.2019

Product code: 160

Page 2 of 6

#### Extinguishing media

##### **Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

#### Specific hazards arising from the chemical

Non-flammable.

#### Special protective equipment and precautions for fire-fighters

In case of fire: Wear self-contained breathing apparatus.

### 6. Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

Use personal protection equipment.

#### Environmental precautions

No special environmental measures are necessary. Clean contaminated articles and floor according to the environmental legislation.

#### Methods and material for containment and cleaning up

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

#### Reference to other sections

Safe handling: see section 7

Personal protection equipment (PPE): see section 8

Disposal: see section 13

### 7. Handling and storage

#### Precautions for safe handling

##### **Advice on safe handling**

No special measures are necessary.

##### **Advice on protection against fire and explosion**

No special fire protection measures are necessary.

#### Conditions for safe storage, including any incompatibilities

##### **Requirements for storage rooms and vessels**

Keep container tightly closed.

##### **Hints on joint storage**

No special measures are necessary.

### 8. Exposure controls/personal protection

#### Control parameters

#### Exposure controls

##### **Protective and hygiene measures**

Take off contaminated clothing. Wash hands before breaks and after work. When using do not eat, drink, smoke, sniff.

##### **Eye/face protection**

Wear eye protection/face protection.

##### **Hand protection**

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves. Recommended glove articles KCL Dermatrill P NBR (Nitrile rubber)

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA LOW FUSING MODELLING LIQUID

Revision date: 05.08.2019

Product code: 160

Page 3 of 6

#### Skin protection

Use of protective clothing.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection. Provide adequate ventilation as well as local exhaust at critical locations. Technical ventilation of workplace

## 9. Physical and chemical properties

### Information on basic physical and chemical properties

Physical state:

Color: colorless

Odor: odorless

pH-Value: 5,5

#### Changes in the physical state

Melting point/freezing point: not determined

Initial boiling point and boiling range: 100 °C

Flash point: ?

#### Flammability

Solid: not applicable

Gas: not applicable

#### Explosive properties

The product is not: Explosive.

Lower explosion limits: not determined

Upper explosion limits: not determined

#### Auto-ignition temperature

Solid: not applicable

Gas: not applicable

Decomposition temperature: not determined

#### Oxidizing properties

Not oxidising.

Vapor pressure: <=1100 hPa  
(at 50 °C)

Density: 1,05 g/cm<sup>3</sup>

Water solubility: No

#### Solubility in other solvents

not determined

Partition coefficient: not determined

Vapor density: not determined

Evaporation rate: not determined

### Other information

Solid content: 0

## 10. Stability and reactivity

### Reactivity

No hazardous reaction when handled and stored according to provisions.

### Chemical stability

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA LOW FUSING MODELLING LIQUID

Revision date: 05.08.2019

Product code: 160

Page 4 of 6

The product is stable under storage at normal ambient temperatures.

#### **Possibility of hazardous reactions**

No known hazardous reactions.

#### **Conditions to avoid**

none/none

#### **Incompatible materials**

No information available.

#### **Hazardous decomposition products**

No known hazardous decomposition products.

## 11. Toxicological information

### **Information on toxicological effects**

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

#### **Irritation and corrosivity**

Based on available data, the classification criteria are not met.

#### **Sensitizing effects**

Based on available data, the classification criteria are not met.

#### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

#### **Specific target organ toxicity (STOT) - single exposure**

Based on available data, the classification criteria are not met.

#### **Specific target organ toxicity (STOT) - repeated exposure**

Based on available data, the classification criteria are not met.

Carcinogenicity (OSHA): No ingredient of this mixture is listed.

Carcinogenicity (IARC): No ingredient of this mixture is listed.

Carcinogenicity (NTP): No ingredient of this mixture is listed.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### **Additional information on tests**

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

## 12. Ecological information

### **Ecotoxicity**

The product is not: Ecotoxic.

### **Persistence and degradability**

The product has not been tested.

### **Bioaccumulative potential**

The product has not been tested.

### **Mobility in soil**

The product has not been tested.

### **Other adverse effects**

No information available.

### **Further information**

Avoid release to the environment.

## 13. Disposal considerations

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA LOW FUSING MODELLING LIQUID

Revision date: 05.08.2019

Product code: 160

Page 5 of 6

#### Waste treatment methods

##### **Disposal recommendations**

Dispose of waste according to applicable legislation.

##### **Contaminated packaging**

Wash with plenty of water. Completely emptied packages can be recycled.

#### 14. Transport information

##### **Marine transport (IMDG)**

**UN number:** No dangerous good in sense of this transport regulation.

**UN proper shipping name:** No dangerous good in sense of this transport regulation.

**Transport hazard class(es):** No dangerous good in sense of this transport regulation.

**Packing group:** No dangerous good in sense of this transport regulation.

##### **Air transport (ICAO-TI/IATA-DGR)**

**UN number:** No dangerous good in sense of this transport regulation.

**UN proper shipping name:** No dangerous good in sense of this transport regulation.

**Transport hazard class(es):** No dangerous good in sense of this transport regulation.

**Packing group:** No dangerous good in sense of this transport regulation.

##### **Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: no

##### **Special precautions for user**

No information available.

##### **Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code**

not applicable

#### 15. Regulatory information

##### **U.S. Regulations**

##### **State Regulations**

##### **Safe Drinking Water and Toxic Enforcement Act of 1986 (Proposition 65, State of California)**

This product can not expose you to chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

#### 16. Other information

Revision date: 05.08.2019

Revision No: 2

##### **Abbreviations and acronyms**

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

CLP: Classification, labelling and Packaging

REACH: Registration, Evaluation and Authorization of Chemicals

## Safety Data Sheet

according to 29 CFR 1910.1200(g)

### VITA LOW FUSING MODELLING LIQUID

Revision date: 05.08.2019

Product code: 160

Page 6 of 6

GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals

UN: United Nations

DNEL: Derived No Effect Level

DMEL: Derived Minimal Effect Level

PNEC: Predicted No Effect Concentration

ATE: Acute toxicity estimate

LL50: Lethal loading, 50%

EL50: Effect loading, 50%

EC50: Effective Concentration 50%

ErC50: Effective Concentration 50%, growth rate

NOEC: No Observed Effect Concentration

BCF: Bio-concentration factor

PBT: persistent, bioaccumulative, toxic

vPvB: very persistent, very bioaccumulative

RID: Regulations concerning the international carriage of dangerous goods by rail

ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
(Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)

EmS: Emergency Schedules

MFAG: Medical First Aid Guide

ICAO: International Civil Aviation Organization

MARPOL: International Convention for the Prevention of Marine Pollution from Ships

IBC: Intermediate Bulk Container

SVHC: Substance of Very High Concern

For abbreviations and acronyms, see table at <http://abbrev.esdscom.eu>

#### Other data

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*